


# cardiofaxVET

Veterinary Electrocardiograph

ECG-1950K



*Fighting Disease with Electronics*

 **NIHON KOHDEN**

## Product Overview



cardiofaxVET, ECG-1950K, is a highly reliable 6-lead electrocardiograph for dogs and cats. It provides not only 6 channel resting recording but also 1 or 3 channel rhythm recording on 110 mm paper from a built-in thermal array printer.

## Features



### Software

cardiofaxVET includes Nihon Kohden's comprehensive ECAPS12V interpretation program.

This easy-to-use program provides analysis with over 60 clinical findings. Depending on the body position, size and type of animal (dog or cat), this sophisticated program provides accurate analysis of acquired signals.

### Color Display (320 x 240 QVGA)

The vivid, color LCD screen shows all 6 traces clearly. This facilitates quicker diagnosis of abnormal heart rhythms.

### Compact and Lightweight

The slim and lightweight 2.0 kg body guarantees stress-free mobility in veterinary clinics and hospitals. With its compact size, cardiofaxVET can be used anywhere.

### Data Storage

Up to 40 ECG files can be stored in internal memory. Approximately 3,000 ECG files can be stored in a 64MB SD card.

### Data Entry

You can easily enter the name of the dog or cat, ID, and other information with the alphanumeric keys on the panel.



## Options



- Animal-use electrodes
- Patient cable (IEC/AHA compliant)
- Power cord
- Battery
- Cart
- Patient cable hanger
- 110 mm recording paper
- SD memory card

## Recording Examples



- 3 channel resting + 1 channel rhythm recording and interpretation



- 1 channel rhythm recording

This brochure may be revised or replaced by Nihon Kohden at any time without notice.



NIHON KOHDEN CORPORATION  
1-31-4 Nishiochiai, Shinjuku-ku, Tokyo 161-8560, Japan  
Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100  
[www.nihonkohden.com](http://www.nihonkohden.com)