

Fighting Disease with Electronics



ACCESSORIES

NEUROLOGY

2 0 1 4

Introduction

Dear Customer,

For more than 60 years NIHON KOHDEN has been one of Japan's leading players in medical technology development and manufacturing and has subsidiaries and representative offices in Europe, the US, as well as in the Middle and Far East.

"Fighting Disease with Electronics" – our corporate philosophy – reflects the high standards we have set ourselves in combating disease by developing medical systems, contributing to higher quality of life for patients.

Our large product range of accessories and consumables is known for its outstanding quality, ease of use and reliability. Please note, that only the use of original NIHON KOHDEN approved accessories and consumables guarantee optimal examination results and make an important contribution towards ensuring that your device will continue to function impeccably and reliably even in years to come, while preserving the device's features and functional capabilities.

As your partner in neurophysiological diagnostics, NIHON KOHDEN provides only tested and approved product solutions from classic electrodes for neuromuscular disorders and brain activity detection to specific products for capnography or amplitude-integrated EEG (aEEG). Our commitment to ensure continually highest quality allows you the best care of your patients and enhances the performances of your electromyographs or electroencephalographs.

NIHON KOHDEN invites you to discover the quality and features of our complete accessories and consumables portfolio in this catalogue. Call today for more information or order directly!

Best regards,

NIHON KOHDEN EUROPE GmbH



Kristin Baar

Product Manager Accessories & Consumables

Accessories & Consumables Neurology

	User Guide	6
	Highest Safety and Quality Standards	7
	Quick Guide Electromyographs, Electroencephalographs, Polysomnography Boxes	8
I.	Accessories & Consumables for Neurophysiological Diagnostic Examinations	9
1.	Needle Electrodes for Electromyography (EMG) Studies	10
1.1.	Disposable Concentric EMG Needle Electrodes	11
1.1.1.	Premium Line	11
1.1.2.	Classic Line	12
1.2.	Reusable Concentric EMG Needle Electrodes	13
1.3.	Reusable Single Fibre & Macro EMG Needle Electrodes	14
1.4.	Reusable Bipolar Concentric Needle Electrodes	15
1.5.	Monopolar Needle Electrodes	16
1.5.1.	Disposable Monopolar & Subdermal Needle Electrodes	16
1.5.2.	Reusable Monopolar & Subdermal Needle Electrodes	17
2.	Electrodes for Evoked Potential (EP) Studies	18
2.1.	EEG/EP Electrodes	19
2.2.	EEG Electrodes	21
3.	Electrodes for Nerve Conduction Studies (NCS)	22
3.1.	Stimulation Electrodes	26
3.2.	Finger Electrodes	28
3.3.	Surface Electrodes	29
3.3.1.	Disposable Surface Electrodes	29
3.3.2.	Reusable Surface Electrodes	31
4.	EMG Ground Electrodes	32
4.1.	Disposable EMG Ground Electrodes	33
4.2.	Reusable EMG Ground Electrodes	33
5.	Additional EMG Accessories	35
5.1.	Acoustic Stimulation	36
5.2.	Visual Stimulation	38
5.3.	Temperature Measurement	38
5.4.	Pastes & Creams	39
5.5.	Recording Paper	40
5.6.	Cables	41
5.7.	Adapter	43

II.	Accessories & Consumables for Electroencephalography (EEG)	44
6.	Accessory Kit for EEG	45
7.	Bridge Electrodes	47
8.	EEG Surface Electrodes	50
8.1.	EEG Disk Electrodes	51
8.2.	EEG Cup Electrodes	52
9.	Subdermal Needle Electrodes	54
10.	Earlobe Electrodes	56
11.	ECG Measurement	58
11.1.	Clip Electrodes	59
11.2.	ECG Electrodes	61
12.	EEG Head Caps	62
12.1.	NIHON KOHDEN Caps	63
12.2.	ElectroCaps	64
12.2.1	ElectroCap Systems	64
12.2.2.	ElectroCaps for Adults	67
12.2.3.	ElectroCaps for Children/Adults	68
12.2.4.	ElectroCaps for Infants	70
12.2.5.	ElectroCap Accessories	71
12.3.	BraiNet Caps	75
13.	Additional EEG Accessories	77
13.1.	Pastes & Creams	78
13.2.	Bags	79
13.3.	Recording Paper	80
13.4.	Recording Ink	81
13.5.	Recording Pens	82
III.	Accessories & Consumables for amplitude-integrated Electroencephalography (aEEG)	83
14.	aEEG Measurement	84
14.1.	aEEG Starter Kits	86
14.2.	aEEG Accessories	89

IV.	Accessories & Consumables for Polysomnography (PSG)	92
<hr/>		
15.	Sleep Study Accessory Kits	93
<hr/>		
16.	Sleep Study Accessories	99
16.1.	EEG Electrodes	100
16.2.	Pastes & Creams	101
16.3.	ECG Electrodes	102
16.4.	Snore Sensors	103
16.5.	Body Position Sensors	104
16.6.	Respiratory Effort Sensors	105
16.6.1.	Respiratory Piezo Crystal Effort Sensors	105
16.6.2.	Respiratory Inductive Plethysmography (RIP) Belts	111
16.7.	Flow Sensors	113
16.7.1.	Reusable Thermistors	113
16.7.2.	Reusable Thermocouples	114
16.7.3.	Airflow Pressure Sensors	116
16.7.4.	Pressure Cannulas	118
16.7.5.	Pressure Transducer	120
16.8.	Limb Movement Sensors	122
16.9.	NIHON KOHDEN SpO ₂ Measurement	124
16.9.1.	NIHON KOHDEN SpO ₂ Connection Cables	124
16.9.2.	NIHON KOHDEN SpO ₂ Probes	125
16.9.2.1.	Reusable SpO ₂ Probes	126
16.9.2.2.	Disposable SpO ₂ Probes	128
16.9.2.2.1.	Disposable SpO ₂ Probes 24 pcs/box – cable length: 0.8 m	128
16.9.2.2.2.	Disposable SpO ₂ Probes 24 pcs/box – cable length: 1.6 m	129
16.9.2.2.3.	Disposable SpO ₂ Probes 5 pcs/box – cable length: 0.8 m	129
16.9.2.2.4.	Disposable SpO ₂ Probes 5 pcs/box – cable length: 1.6 m	131
16.9.2.2.5.	Disposable Multi-Site SpO ₂ Probes 5 pcs/box – cable length: 2.0 m	132
16.10.	NONIN SpO ₂ Probes	133
16.11.	ETCO ₂ Measurement (Capnography)	134
16.11.1.	Intubated Expired and Inspired ETCO ₂ cap-ONE – „Opening a New Era in Capnography“	136
16.11.2.	Intubated and Non-intubated Expired ETCO ₂ cap-ONE – „Opening a New Era in Capnography“	137
16.12.	Sensor Tester	139
<hr/>		
	Accessory Order Form	140
<hr/>		
	Index as per Supply Code	141
<hr/>		
	Index as per Model Code	144
<hr/>		

User Guide

Care of Environment

NIHON KOHDEN is making every possible effort to ensure that the actions of the company and its employees contribute to preserving the environment. In order to conserve paper resources, we decided a more environment friendly approach which allows us to update the catalogue more frequently in the future. Therefore the new accessories and consumables catalogue was designed for view on your computer. It allows you to zoom pictures in order to see more details of the product. If you need indeed certain chapters in print version, you can easily select and print them.

Find the Right Accessories in the Catalogue

You do not know exactly which article you need for your examination?

Following this user guide you will find the content of this catalogue which is organized into parameter categories. In a particular category you will see a list of accessories and consumables with corresponding supply codes, model codes, descriptions and pictures. If you are not sure if one of the listed products will fit to your requirements, please contact us or one of our regional representatives.

You already know the article that you need for your examination?

In the index of the catalogue, from page 141 to page 145 the articles are listed alphabetically as per supply code; from page 144 to page 145, the products are listed alphabetically as per model code. They will help you to find directly the according product page.

If you are using the PDF version, simply use the "search function" of this by entering the description, supply or model code. Or just click on the related chapter on the left side of the document.

Order Information

To order NIHON KOHDEN accessories and consumables or to get more information, please refer to the indicated regional NIHON KOHDEN representatives below.

For your convenience, we have a printable order form at the end of this catalogue so that you can fax in your order directly or mail/e-mail it to us. Please make sure to fill out all requested information as well as the required supply code and according quantity. Please check the general terms and conditions including delivery time, payment and general warranty with your responsible regional NIHON KOHDEN representatives:

For Germany, Austria, Switzerland and Benelux

NIHON KOHDEN EUROPE GmbH, Vertrieb Deutschland
Albert-Einstein-Ring 9, 14532 Kleinmachnow, Deutschland
Fax: +49 (0) 33203/573-19, E-Mail: bestellung@nke.de

For Spain and Portugal

NIHON KOHDEN IBERICA S.L.
C/Ulises, 75A, 28043 Madrid, Spain
Fax: +34 (0) 91 300 46 76, E-Mail: info@nkib.es

For United Kingdom

NIHON KOHDEN UNITED KINGDOM Ltd.
Tolworth Tower – 19th Floor, Ewell Road, Surbiton,
Surrey KT6 7EL, United Kingdom
Fax: +44 (0) 20 83 90 46 75, E-Mail: mail@nihonkohden.co.uk

For France

NIHON KOHDEN FRANCE SARL
8, rue François Delage, 94230 Cachan, France
Fax: +33 (0) 1 49 08 93 32, E-Mail: info@nkfrance.fr

For Italy

NIHON KOHDEN ITALIA S.r.l.
Via Fratelli Bronzetti, 28, 24124 Bergamo, Italy
Fax: +39 (0) 35-232546, E-Mail: info@nkitaly.com

For all other European Countries, Russia, Turkey and Israel

NIHON KOHDEN EUROPE GmbH
Raiffeisenstraße 10, 61191 Rosbach v. d. H., Germany
Fax: +49 (0) 6003/827-599, E-Mail: info@nke.de

Highest Safety and Quality Standards

Comprehensive quality assurance systems ensure that each of the more than 4,000 NIHON KOHDEN staff members makes his/her contribution towards development, manufacture, marketing and customer services such that maximum product quality is guaranteed.

NIHON KOHDEN is certified according to ISO 9001, ISO 13485 and ISO 14001.

NIHON KOHDEN designs and manufactures medical equipment in conformity with the standards of IEC 60601-1 (Medical Electrical Equipment: General Requirements for Safety).

NIHON KOHDEN products have earned a reputation for quality and long life. One reason is that all export products and major components are manufactured in our own ISO 9001 and IEC 13485 certified factories.

In 1984, NIHON KOHDEN implemented a comprehensive companywide quality system, IQC – Integrated Quality Control. Under this system, quality begins at the earliest stage of every product's life. Our IQC applies to the entire manufacturing process, from the planning stage through product development, design, parts procurement, production, inspection, storage, sales, distribution and after-sales support. IQC also involves nearly every employee.

Our rigid test procedures are structured to exceed the provisions of the GMP, Good Manufacturing Practice, codes enforced by the US Food and Drug Administration. NIHON KOHDEN also has 130 in-house quality control inspection staff to ensure that standards are implemented and maintained in every department of the company.

IQC extends beyond production and into sales and support at each subsidiary. This staff is responsible for ongoing contact with customers to determine their needs and making sure that these changing requirements are reflected in new products. Service engineers are responsible for communicating any problem that may arise in the use of NIHON KOHDEN products, and assuring the customers of NIHON KOHDEN's commitment to their satisfaction.

In 1995, NIHON KOHDEN received ISO 9001 certification, the international standard for quality systems. In addition, NIHON KOHDEN's four manufacturing facilities in Japan are the first factories in the medical equipment industry to be certified by the CSA, Canadian Standards Association, under their category certification system. Most products sent to North America carry the CSA label.

Quality improvement is and always has been an ongoing process for NIHON KOHDEN.

Convince yourself of our expertise and the high quality of our products and services! Order directly!

Quick Guide

Electromyographs, Electroencephalographs, Polysomnography Boxes

ELECTROMYOGRAPHS



MEE-1000



MEB-9400



MEB-9200



MEB-2300



MEB-2200



MEB-9100



Combo



MEB-7102



MEB-5500



MEB-5200



MEB-4200

ELECTROENCEPHALOGRAPHS



EEG-1200



WEE-1000



EEG-9200



EEG-9100



EEG-1100



Combo



EEG-8300



EEG-7400




EEG-7300



EEG-5500



EEG-5600



EEG-4400

POLYSOMNOGRAPHY BOXES



PSG-1100



JE-921A



JE-912AK



NF-SCREENER3

I. Accessories & Consumables for Neurophysiological Diagnostic Examinations

Labour-intensive neurophysiological diagnostic examinations place particular demands on the ergonomics of the diagnostic measurement systems used. The development of the first electromyographs by NIHON KOHDEN in 1956 laid the foundations for a comprehensive product range. Reliable EMG measurements depend on consistent impedances (electrode-skin junction). The quality of the electrode is therefore extremely important to obtain a noise-free high signal quality. Inappropriate working electrodes cause additional manual adjustments by the user which increases the time of measurement. Since the beginning NIHON KOHDEN electrodes are known for their outstanding quality which allows a quick and reliable recording. Thanks to our extensive experience, we are able to offer you ideal product solutions for hospitals and medical practices.

1. Needle Electrodes for Electromyography (EMG) Studies

Electromyography (EMG) studies are a technique to verify physical condition of the muscles and the nerves that manage the muscles. For this purpose either surface (or skin) electrodes or minimal invasive electrodes (needle electrodes or wires) can be used. Whereas the first types of electrodes are mainly used for general muscle activity observation, the second types are usually employed to identify electrical activities of a nerve root.

To record the electrical behaviour of the muscles a very thin needle electrode is percutaneously inserted through the skin into the muscle. EMG studies are helpful to detect muscle disorders during rest or activity.

1.1. Disposable Concentric EMG Needle Electrodes

1.2. Reusable Concentric EMG Needle Electrodes

1.3. Reusable Single Fibre & Macro EMG Needle Electrodes

1.4. Reusable Bipolar Concentric Needle Electrodes

1.5. Monopolar Needle Electrodes

1.1. Disposable Concentric EMG Needle Electrodes

Facts

STRAIGHTFORWARD, COMFORTABLE, EFFECTIVE

NIHON KOHDEN disposable concentric EMG needle electrodes are characterized by an easy handling. The ergonomic colour-coded hub design allows an easy identification of the needle size. For your different kinds of applications, NIHON KOHDEN provides you a complete portfolio of suitable needle lengths and diameters. A high patient comfort is assured by an excellent sharpness of the needle tip for a soft skin penetration and a minimized patients' pain. Additionally the tip geometry makes less penetration force necessary and the orientation-free connection to the connection cable simplifies the usage. High signal quality is achieved by low impedances. The stable handling and good grip due to the ribbed hub design guarantee confident and reliable recordings.

1.1.1. Premium Line

Supply Code	Model	Description
NKD-TR4537P		Premium Line Needle Trial Kit, incl. 6 NKD-NM4537P & 1 NKD-BM113A For current products: MEB-9400/2300 Packing unit: 1 kit
NKD-NM4525P		Disposable Concentric EMG Needle Electrode, Premium Line, colour: violet, size: 0.45 x 25 mm, requires NKD-BM113A For current products: MEB-9400/2300 Packing unit: 25 pcs D
NKD-NM4537P		Disposable Concentric EMG Needle Electrode, Premium Line, colour: orange, size: 0.45 x 37 mm, requires NKD-BM113A For current products: MEB-9400/2300 Packing unit: 25 pcs D
NKD-NM4550P		Disposable Concentric EMG Needle Electrode, Premium Line, colour: pink, size: 0.45 x 50 mm, requires NKD-BM113A For current products: MEB-9400/2300 Packing unit: 25 pcs D
NKD-BM113A		Connection Cable for Disposable EMG Needles NKD-NMxxxx, 5-PIN DIN connector, cable length: 1.5 m For current products: MEB-9400/2300 Packing unit: 1 pcs R



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

1.1.2.

Classic Line

Supply Code	Model	Description
NKD-TR4537		Classic Line Needle Trial Kit, incl. 10 NKD-NM4537 & 1 NKD-BM113A For current products: MEB-9400/2300 Packing unit: 1 kit
NKD-NM3525		Disposable Concentric EMG Needle Electrode, Classic Line, colour: dark grey, size: 0.35 x 25 mm, requires NKD-BM113A For current products: MEB-9400/2300 Packing unit: 25 pcs D
NKD-NM4525		Disposable Concentric EMG Needle Electrode, Classic Line, colour: black, size: 0.45 x 25 mm, requires NKD-BM113A For current products: MEB-9400/2300 Packing unit: 25 pcs D
NKD-NM4537		Disposable Concentric EMG Needle Electrode, Classic Line, colour: blue, size: 0.45 x 37 mm, requires NKD-BM113A For current products: MEB-9400/2300 Packing unit: 25 pcs D
NKD-NM4545		Disposable Concentric EMG Needle Electrode, Classic Line, colour: yellow, size: 0.45 x 45 mm, requires NKD-BM113A For current products: MEB-9400/2300 Packing unit: 25 pcs D
NKD-NM4540		Disposable Concentric EMG Needle Electrode, Classic Line, colour: red, size: 0.45 x 40 mm, requires NKD-BM113A For current products: MEB-9400/2300 Packing unit: 25 pcs D
NKD-NM4550		Disposable Concentric EMG Needle Electrode, Classic Line, colour: grey, size: 0.45 x 50 mm, requires NKD-BM113A For current products: MEB-9400/2300 Packing unit: 25 pcs D
NKD-NM5565		Disposable Concentric EMG Needle Electrode, Classic Line, colour: beige, size: 0.55 x 65 mm, requires NKD-BM113A For current products: MEB-9400/2300 Packing unit: 25 pcs D
NKD-BM113A		Connection Cable for Disposable EMG Needles NKD-NMxxxx, 5-PIN DIN connector, cable length: 1.5 m For current products: MEB-9400/2300 Packing unit: 1 pcs R



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

1.2. Reusable Concentric EMG Needle Electrodes

Facts

Before use, autoclave the needle electrode with the tip of the needle covered with protective tube at 121 °C and 1 bar overpressure for 30 minutes. If the needle tip is deformed or the teflon coating on the needle is damaged, stop using the electrode and replace it with a new one.

Supply Code	Model	Description
H630	NM-121T	Reusable Concentric EMG Needle Electrode, platinum, size: 0.45 x 20 mm, requires K611/K623 For current products: MEB-9400/2300 Packing unit: 4 pcs 
H631	NM-151T	Reusable Concentric EMG Needle Electrode, platinum, size: 0.45 x 50 mm, requires K611/K623 For current products: MEB-9400/2300 Packing unit: 4 pcs 
H650	NM-131T	Reusable Concentric EMG Needle Electrode, platinum, size: 0.45 x 30 mm, requires K611/K623 For current products: MEB-9400/2300 Packing unit: 4 pcs 
H651	NM-320T	Reusable Concentric EMG Needle Electrode, stainless steel, size: 0.32 x 20 mm, requires K611/K623 For current products: MEB-9400/2300 Packing unit: 4 pcs 
H652	NM-330T	Reusable Concentric EMG Needle Electrode, stainless steel, size: 0.32 x 30 mm, requires K611/K623 For current products: MEB-9400/2300 Packing unit: 4 pcs 
K611	BM-121S	Connection Cable for H630/631/650/651/652/660, 6-PIN round connector, cable length: 1.0 m For current products: MEB-9400/2300 Packing unit: 1 pcs 
K623	BM-115S	Connection Cable for H630/631/650/651/652/660, 6-PIN round connector, cable length: 1.8 m For current products: MEB-9400/2300 Packing unit: 1 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!



 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

1.3. Reusable Single Fibre & Macro EMG Needle Electrodes

Facts

Before use, autoclave the needle electrode with the tip of the needle covered with protective tube at 121 °C and 1 bar overpressure for 30 minutes. If the needle tip is deformed or the teflon coating on the needle is damaged, stop using the electrode and replace it with a new one.

Supply Code	Model	Description
H650A	NM-130S	Reusable Concentric Macro EMG Electrode, platinum, size: 0.45 x 30 mm, teflon isolation: 15 mm, requires K619 For current products: MEB-9400/2300 Packing unit: 1 pcs 
H660	NM-640T	Reusable Single Fibre Electrode, platinum/wolfram, size: 0.45 x 30 mm, internal dia: 0.03 mm, requires K611/K623 For current products: MEB-9400/2300 Packing unit: 1 pcs 
H661	NM-640S	Reusable Single Fibre Macro EMG Electrode, platinum/wolfram, size: 0.45 x 30 mm, teflon isolation: 15 mm, internal dia: 0.03 mm, requires K619 For current products: MEB-9400/2300 Packing unit: 1 pcs 
K619	BM-840S	Connection Cable for H650A/661, 6-PIN DIN round connector, cable length: 1.0 m For current products: MEB-9400/2300 Packing unit: 1 pair 
K611	BM-121S	Connection Cable for H630/631/650/651/652/660, 6-PIN DIN round connector, cable length: 1.0 m For current products: MEB-9400/2300 Packing unit: 1 pcs 
K623	BM-115S	Connection Cable for H630/631/650/651/652/660, 6-PIN DIN round connector, cable length: 1.8 m For current products: MEB-9400/2300 Packing unit: 1 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!


 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

1.4. Reusable Bipolar Concentric Needle Electrodes

Facts

Before use, autoclave the needle electrode with the tip of the needle covered with protective tube at 121 °C and 1 bar over-pressure for 30 minutes. If the needle tip is deformed or the teflon coating on the needle is damaged, stop using the electrode and replace it with a new one.

Supply Code	Model	Description
H632	NM-220T	Reusable Bipolar Concentric Needle Electrode, size: 0.40 x 20 mm, requires K612/K624 For current products: MEB-9400/2300 Packing unit: 2 pcs 
H633	NM-250T	Reusable Bipolar Concentric Needle Electrode, size: 0.40 x 50 mm, requires K612/K624 For current products: MEB-9400/2300 Packing unit: 2 pcs 
K612	BM-211S	Connection Cable for H632/633, 6-PIN DIN round connector, cable length: 1.0 m For current products: MEB-9400/2300 Packing unit: 1 pcs 
K624	BM-215S	Connection Cable for H632/633, 6-PIN DIN round connector, cable length: 1.8 m For current products: MEB-9400/2300 Packing unit: 1 pcs 



1.5. Monopolar Needle Electrodes

1.5.1. Disposable Monopolar & Subdermal Needle Electrodes

Supply Code	Model	Description
NKD-EPD1		Disposable Monopolar Needle Electrode for EP, straight type, steel, size: 0.35 x 12 mm, 1 mm connector, cable length: 8.0 cm, sterilised, requires NKD-EPDK For current products: MEB-9400/2300; MEE-1000 Packing unit: 25 pcs 



NKD-EPD2		Disposable Monopolar Needle Electrode for EP, bent type, steel, size: 0.35 x 12 mm, 1 mm connector, cable length: 8.0 cm, sterilised, requires NKD-EPDK For current products: MEB-9400/2300; MEE-1000 Packing unit: 25 pcs 
----------	--	--



NKD-EPDK		Connection Cable for NKD-EPD1/NKD-EPD2, 1.5 mm touch proof connector to 1 mm connector, cable length: 1.5 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 10 pcs 
----------	--	--



NE-110B		Disposable Subdermal Needle Electrode Prewired, single type, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 1.0 m For current products: EEG-1200/9100; MEB-9400/2300; MEE-1000 Packing unit: 25 pcs 
---------	--	---



NE-115B		Disposable Subdermal Needle Electrode Prewired, single type, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 1.5 m For current products: EEG-1200/9100; MEB-9400/2300; MEE-1000 Packing unit: 25 pcs 
---------	--	---

NE-120B		Disposable Subdermal Needle Electrode Prewired, single type, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 2.0 m For current products: EEG-1200/9100; MEB-9400/2300; MEE-1000 Packing unit: 25 pcs 
---------	--	---

NE-125B		Disposable Subdermal Needle Electrode Prewired, single type, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 2.5 m For current products: EEG-1200/9100; MEB-9400/2300; MEE-1000 Packing unit: 25 pcs 
---------	--	---

NE-215B		Disposable Subdermal Needle Electrode Prewired, twisted pair, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 1.5 m For current products: EEG-1200/9100; MEB-9400/2300; MEE-1000 Packing unit: 20 pcs 
---------	--	--



NE-220B		Disposable Subdermal Needle Electrode Prewired, twisted pair, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 2.0 m For current products: EEG-1200/9100; MEB-9400/2300; MEE-1000 Packing unit: 20 pcs 
---------	--	--

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

1.5.2.

Reusable Monopolar & Subdermal Needle Electrodes

Supply Code **Model** **Description**
H537A **NE-224S** Reusable Subdermal EEG Needle Electrode, stainless steel, 1.5 mm touch proof connector, size: 0.22 x 22.5 mm, cable length: 1.5 m
For current products: EEG-1200/9100; MEB-9400/2300; BSM-3000/6000/9100
Packing unit: 20 pcs **R**



H665 **NM-710T** Reusable Monopolar Needle Electrode for EP/EMG, stainless steel, size: 0.25 x 13 mm, requires K616
For current products: MEB-9400/2300
Packing unit: 5 pcs **R**



H666 **NM-715S** Reusable Monopolar Needle Electrode for EP/EMG, stainless steel, teflon-coated isolation, size: 0.3 x 15 mm, requires K620
For current products: MEB-9400/2300
Packing unit: 5 pcs **R**




H667 **NM-730S** Reusable Monopolar Needle Electrode for EP/EMG, stainless steel, teflon-coated isolation, size: 0.3 x 30 mm, requires K620
For current products: MEB-9400/2300
Packing unit: 5 pcs **R**



H668 **NM-830S** Reusable Monopolar Needle Electrode for EP/EMG, stainless steel, teflon-coated isolation, size: 0.4 x 30 mm, requires K620
For current products: MEB-9400/2300
Packing unit: 5 pcs **R**



H669 **NM-850S** Reusable Monopolar Needle Electrode for EP/EMG, stainless steel, teflon-coated isolation, size: 0.4 x 50 mm, requires K620
For current products: MEB-9400/2300
Packing unit: 5 pcs **R**



K616 **BM-710S** Connection Cable for H665, shielded, 2 mm connector, cable length: 1.0 m
For old product: MEB-4200/5200/5300/5500/7102
Packing unit: 5 pcs **R**



K620 **BM-800S** Connection Cable for H666/H667/H668/H669, shielded, 6-PIN DIN round connector, cable length: 1.0 m
For current products: MEB-9400/2300
Packing unit: 1 pcs **R**



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

2. Electrodes for Evoked Potential (EP) Studies

Evoked potential (EP) studies are tests to measure electrical activity along the nerve pathways in response to visual, acoustical or electrical stimulations. Those stimulations create meticulous electrical signals that pass through the nerve tract and the spinal cord to certain regions of the brain where they are detected by electrodes. For visual evoked response (VEP) tests the electrodes are placed on the scalp. For a brainstem auditory evoked response (ABR) test electrodes are applied on the head or the earlobe whereas they are attached to the wrist, the back of the knee or other places for a somatosensory evoked response (SEP) test.

EP studies are maintained to identify sensory nerve or optic nerve disorders. In addition to this, hearing tests, the indication of brain stem tumours and multiple sclerosis as well as the observation of brain activity among coma patients and the confirmation of brain death, can be realized by those tests. NIHON KOHDEN is focusing on highest signal quality and is providing premium accessories and consumables.

2.1. EEG/EP Electrodes

2.2. EEG Electrodes

2.1. EEG/EP Electrodes

Supply Code	Model	Description
H528	NE-136B	Reusable Evoked EEG Disk Electrodes, Ag/AgCl, cable length: 0.3 m, for BM-230B For current products: MEB-2300 Packing unit: 3 pcs R
H852A	NE-132B	Reusable Evoked EEG Disk Electrodes Kit, Ag, shielded cables, keyhole connector, electrode cable length: 1.5 m Note: Kit contains 5 x disk electrodes & 1 x jumper cable K617C. For current products: MEB-9400/2300 Packing unit: 1 kit R
H854	NE-142B	Reusable Evoked EEG Disk electrodes Kit, Ag, shielded cables, keyhole connector, for ABR, electrode cable length: 0.5 m Note: Kit contains 5 x disk electrodes & 1 x jumper cable K617C. For current products: MEB-9400/2300 Packing unit: 1 kit R
H851	NE-121J	Reusable Evoked EEG Disk Electrodes Kit, Ag, shielded cables, 2 mm connector, electrode cable length: 1.5 m Note: Kit contains 5 x disk electrodes & 1 x jumper cable K617. For old products: MEB-5200/5300/7102 Packing unit: 1 kit R
H851A	NE-131B	Reusable Evoked EEG Disk Electrodes Kit, Ag, shielded cables, 1.5 mm touch proof connector, electrode cable length: 1.5 m Note: Kit contains 5 x disk electrodes & 1 x jumper cable K617B. For old products: MEB-2200/9100/9200 Packing unit: 1 kit R
H852	NE-121B	Reusable Evoked EEG Disk Electrodes Kit, Ag, shielded cables, 2 mm connector, electrode cable length: 1.5 m Note: Kit contains 5 x disk electrodes & 1 x jumper cable K617A. For old products: MEB-4200/5500 Packing unit: 1 kit R




D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
K626	BM-910B	EP Connection Cable for H854, shielded, 9-PIN DIN round connector, cable length: 3.0 m For current products: MEB-9400 Packing unit: 1 pcs 



K617	2114-064171	EMG Jumper Cable for H851, 2 mm dia 1-PIN connector, unshielded, cable length: 0.15 m For old products: MEB-5200/5300/7102 Packing unit: 5 pcs 
------	-------------	--



K617A	NE-420B	EMG Jumper Cable for H852, 2 mm dia 2-PIN connector, unshielded, cable length: 0.15 m For old products: MEB-4200/5500 Packing unit: 5 pcs 
-------	---------	---



K617B	NE-430B	EMG Jumper Cable for H851A, 1-PIN DIN round connector, unshielded, cable length: 0.15 m For current products: MEB-2200/9100/9200 Packing unit: 5 pcs 
-------	---------	---



K617C	NE-432B	EP Jumper Cable for H852A/H854, 2-PIN DIN keyhole connector, unshielded, cable length: 0.15 m For current products: MEB-9400/2300 Packing unit: 5 pcs 
-------	---------	---




2.2. EEG Electrodes

Facts

Attaching the electrodes

1. To reduce the skin impedance, clean the electrode sites with a cotton pad moistened with alcohol or rub them with skinPure skin preparation gel (F020). Thoroughly dry the skin with dry gauze.
2. Apply an optimal amount of Elefix EEG conduction paste (F509 or F510) to the cleaned sites and rub the paste around the sites in a circle not larger than 10 mm in diameter.
3. Apply Elefix EEG conduction paste about 1 mm thick to the disk (cup) of the electrode, and attach it on the pasted skin sites.
4. Cover the electrode with a piece of cotton, and press the edge of the cotton. Fasten the electrode on the site with surgical tape if necessary.
5. Connect the electrode leads to the corresponding input jacks on the electrode junction box.



Supply Code	Model	Description
H503A	NE-113A	EEG Surface Electrode for EOG & EMG, Ag, size: 11 mm dia, with brown plastic covering, 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100; MEE-1000; MEB-2300 Packing unit: 26 pcs 



H503B	NE-114A	EEG Surface Electrode for EOG & EMG, Ag, size: 11 mm dia, 1.5 mm touch proof connector, cable length: 0.8 m For current products: EEG-1200/9100; MEE-1000; MEB-2300 Packing unit: 26 pcs 
-------	---------	--

Supply Code	Model	Description
H526	NE-134A	EEG Electrode (Collodion), Ag, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100; BSM-3000/6000/9100; PSG-1100; MEB-2300; MEE-1000 Packing unit: 12 pcs 



H527	NE-136A	EEG Electrode (Collodion), Ag, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 0.7 m For current products: EEG-1200/9100; BSM-3000/6000/9100; PSG-1100; MEB-2300; MEE-1000 Packing unit: 12 pcs 
------	---------	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

3. Electrodes for Nerve Conduction Studies (NCS)


Nerve Conduction Studies (NCS) measure the conduction rapidity of an electrical impulse through a nerve. Usually the test is done with two surface electrodes. One electrode is directly attached over the identified nerve stimulates the nerve with a soft electrical impulse. The other electrode, placed over the muscle controlled by this nerve, records the time the muscle needs to contract as well as the signal quality sent between both electrodes in response to this electrical stimulation. The response speed is called the conduction velocity.

Contrary to EMG studies which detect the proper functionality of muscles on nerve stimulations (muscle disorders), NCS examine disorders of the peripheral nervous system. Both measurements combined, assist to sense the existence, place and degree of disorder that injure the nerves and muscles.

3.1. Stimulation Electrodes

3.2. Finger Electrodes

3.3. Surface Electrodes



- **Electrode position**
 - Active electrode (+): belly of abductor pollicis brevis.
 - Reference electrode (-): joint of thumb finger.
- **Stimulation point**
 - Palm stimulation point: center of palm, set cathode to proximal.
 - Distal stimulation point: between the tendons of the palmaris longus and flexor carpi ulnaris. Keep the distance equally at every study.
 - Proximal stimulation point: just lateral to the brachial artery.
- **F-wave**

Electrode position is the same as MCS. Cathode of stimulus is on proximal side.


Median nerve MCS normative data (Kimura 1983)

Stimulation point	Amplitude (mV) Mean±SD (lower limit)	Latency (ms) Mean±SD (upper limit)	Conduction velocity (m/s) Mean±SD (lower limit)
Palm	6.9±3.2 (3.5)	1.86±0.28 (2.4)	
Distal	7.0±3.0 (3.5)	3.49±0.34 (4.2)	48.8±5.3 (38)
Proximal	7.0±2.7 (3.5)	7.39±0.69 (8.8)	57.7±4.9 (48)
Axilla	7.2±2.9 (3.5)	9.81±0.89 (11.6)	63.5±6.2 (51)

Note: Amplitude is the potential between the negative peak of the waveform and the baseline.

Median nerve F-wave normative data (Kimura 1983)

Stimulation point	Latency (ms) Mean±SD (upper limit)	Conduction velocity (m/s) Mean±SD (lower limit)
Distal	26.6±2.2 (31)	65.3±4.7 (56)




- **Electrode position**
 - Active electrode (+): PIP joint of index finger.
 - Reference electrode (-): DIP point of index finger.
- **Stimulation point**
 - Palm stimulation point: center of palm, set cathode to distal.
 - Distal stimulation point: between the tendons of the palmaris longus and flexor carpi ulnaris. Keep the distance equally at every study.
 - Proximal stimulation point: just lateral to the brachial artery.

Median nerve SCS normative data (Kimura 1983)

Stimulation point	Amplitude (mV) Mean±SD (lower limit)	Latency (ms) Mean±SD (upper limit)	Conduction velocity (m/s) Mean±SD (lower limit)
Palm	39.0±16.8 (20)	1.37±0.24 (1.9)	58.8±5.8 (47)
Distal	38.5±15.6 (19)	2.84±0.34 (3.5)	56.2±5.8 (44)
Proximal	32.0±15.5 (16)	6.46±0.71 (7.9)	61.9±4.2 (53)

Note: Amplitude is the potential between the negative peak of the waveform and the baseline.



- **Electrode position**
 - Active electrode (+): Belly of ADM (Abductor Digiti Minimi).
 - Reference electrode (-): MP joint of the 5th finger.
- **Stimulation point**
 - Distal stimulation point: just over or lateral to the flexor carpi ulnaris. Keep the distance equally at every study.
 - Below elbow stimulation point: 5 - 6 cm distal to the groove for ulnar nerve.
 - Above elbow stimulation point: 5 - 6 cm proximal to the groove.
- **F-wave**

Electrode position is the same as MCS. Cathode of stimulus is on proximal side.


Ulnar nerve MCS normative data (Kimura 1983)

Stimulation point	Amplitude (mV) Mean±SD (lower limit)	Latency (ms) Mean±SD (upper limit)	Conduction velocity (m/s) Mean±SD (lower limit)
Distal	5.7±2.0 (2.8)	2.59±0.39 (3.4)	58.7±5.1 (49)
Below elbow	5.5±2.0 (2.7)	6.10±0.69 (7.5)	61.0±5.5 (50)
Above elbow	5.5±1.9 (2.7)	8.04±0.76 (9.6)	66.5±6.3 (54)
Axilla	5.6±2.1 (2.7)	9.90±0.91 (11.7)	

Note: Amplitude is the potential between the negative peak of the waveform and the baseline.

Ulnar nerve F-wave normative data (Kimura 1983)

Stimulation point	Latency (ms) Mean±SD (upper limit)	Conduction velocity (m/s) Mean±SD (lower limit)
Distal	27.6±2.2 (32)	65.3±4.8 (55)




- **Electrode position**
 - Active electrode (+): PIP joint of 5th finger.
 - Reference electrode (-): DIP joint of 5th finger.
- **Stimulation point**
 - Distal stimulation point: just over or lateral to the flexor carpi ulnaris. Keep the distance equally at every study.
 - Below elbow stimulation point: 5 - 6 cm distal to the groove for ulnar nerve.
 - Above elbow stimulation point: 5 - 6 cm proximal to the groove for ulnar nerve.

Ulnar nerve SCS normative data (Kimura 1983)


Stimulation point	Amplitude (mV) Mean±SD (lower limit)	Latency (ms) Mean±SD (upper limit)	Conduction velocity (m/s) Mean±SD (lower limit)
Distal	35.0±14.7 (18)	2.54±0.29 (3.1)	54.8±5.3 (44)
Below elbow	28.8±12.2 (15)	5.67±0.59 (6.9)	64.7±5.4 (53)
Above elbow	28.3±11.8 (14)	7.46±0.64 (8.7)	66.7±6.4 (54)

Note: Amplitude is the potential between the negative peak of the waveform and the baseline.




- **Electrode position**
 - Active electrode (+): 2 cm proximal from interdigit index finger and middle finger.
 - Reference electrode (-): PIP joint of index finger.
- **Stimulation point**
 - Median and Ulnar nerve, equal distance to the active electrode

Lumbrical muscle interosseous muscle motor conduction study is a very sensitive method for CTS (carpal tunnel syndrome). An active electrode placed as shown in the figure both picks up compound muscle action potential from 2nd lumbrical muscle innervated by median nerve and interosseous muscle innervated by ulnar nerve. Stimulation is applied an equal distance on median and ulnar nerve to the active electrode. Compare the latency of each CMAPS, then difference between them is higher than 0.5 ms suggests CTS.



- **Electrode position**
 - Active electrode (+): PIP joint of ring finger.
 - Reference electrode (-): DIP joint of ring finger.
- **Stimulation point**
 - Median and Ulnar nerve, equal distance to the active electrode.

Ring finger method is another sensitive method for CTS. A ring electrode is placed on the ring finger as shown in the figure. The radial side of ring finger is innervated by median nerve, ulnar side is by ulnar nerve. As for lumbrical muscle use the interosseous muscle method, and stimulate the median nerve and ulnar nerve at equal distance to the active electrode. Compare each SNAP latency and measure the difference. If it is higher than 0.5 ms, it suggests CTS.



- **Electrode position**
 - Active electrode (+): Belly of the abductor hallucis muscle (AH).
 - Reference electrode (-): Joint of large toe (hallux).
- **Stimulation point**
 - Articulatio talo crutails stimulation point: posterior boarder of malleolus medialis. Keep the distance between the cathode stimulation point and active electrode to compare with the distal latency.
 - Distal: 10 cm proximal from the active electrode.
 - Articulatio genuis stimulation point: The central portion of the popliteal fossa.
- **F-wave**

The electrode position is the same as that of MCS. Cathode of stimulus is on proximal side.


Tibial nerve MCS normative data (Kimura 1983)

Stimulation point	Amplitude (mV) Mean±SD (lower limit)	Latency (ms) Mean±SD (uppepr limit)	Conduction velocity (m/s) Mean±SD (lower limit)
Distal stimulation point	5.8±1.9 (2.9)	3.96±1.00 (6.0)	48.5±3.6 (41)
Proximal stimulation point	5.1±2.2 (2.5)	12.05±1.53 (15.1)	

Note: Amplitude is the potential between the negative peak of the waveform and the baseline.

Tibial nerve F-wave normative data (Kimura 1983)

Stimulation point	Latency (ms) Mean±SD (uppepr limit)	Conduction velocity (m/s) Mean±SD (lower limit)
Articulatio talo crutails	47.7±5.0 (58)	52.6±4.3 (44)



- **Electrode position**
 - Active electrode (+): Belly of the extensor digitorum brevis muscle (EDB).
 - Reference electrode (-): Joint of little toe (digitus minimus).
- **Stimulation point**
 - Articulatio talo crutails stimulation point: lateral to the tibialis anterior tendon. Keep the distance between the cathode stimulation point and active electrode to compare with the distal latency.
 - Stimulation point below the fibula head: 2 cm distal to the fibula head.
 - Stimulation point above the fibula head: 10 cm proximal to the fibula head just medial to the end of the biceps femoris tendon.
- **F-wave**

The electrode position is the same as that of MCS. Cathode of stimulus is on proximal side.


Peroneal nerve MCS normative data (Kimura 1983)

Stimulation point	Amplitude (mV) Mean±SD (lower limit)	Latency (ms) Mean±SD (uppepr limit)	Conduction velocity (m/s) Mean±SD (lower limit)
Articulatio talo crutails	5.1±2.3 (2.5)	3.77±0.86 (5.5)	48.3±3.9 (40)
Below the fibula head	5.1±2.0 (2.5)	10.79±1.06 (12.9)	52.0±6.2 (40)
Above the fibula head	5.1±1.9 (2.5)	12.51±1.17 (14.9)	

Note: Amplitude is the potential between the negative peak of the waveform and the baseline.

Peroneal nerve F-wave normative data (Kimura 1983)

Stimulation point	Latency (ms) Mean±SD (uppepr limit)	Conduction velocity (m/s) Mean±SD (lower limit)
Articulatio talo crutails	48.4±4.0 (56)	49.8±3.6 (43)



- **Electrode position**
 - Active electrode (+): Posterior to the lateral malleolus one third of the distance to the lateral malleolus along the line drawn between the heel and the malleolus.
 - Reference electrode (-): A disk electrode is placed 3 cm distal to the active electrode.
- **Stimulation point**
 - 14 cm proximal to the active electrode and 1 to 2 cm lateral to the Achilles tendon.

Sural nerve SCS normative data (Kimura 1983)


Stimulation point	Age	Amplitude (mV) Mean±SD	Latency (ms) Mean±SD	Conduction velocity (m/s) Mean±SD
Wrist joint	10-40	20.9±8.0	2.7±0.3	52.5±5.6
	41-84	17.2±6.7	2.8±0.3	51.1±5.9

Note: Amplitude is the potential between the negative peak of the waveform and the baseline.


Correlation between F-wave minimum latency and height

Ulnar nerve Correlation between F-wave minimum latency and height (Distal stimulation at upper extremity)


F-wave minimum latency has a strong correlation with height, and you have to consider the patient's height when referring to average value.



Tibial nerve Correlation between F-wave minimum latency and height (Distal stimulation at lower extremity)



Reference: Nobrega J. A. M., Pinheiro D. S., Manzano G. M. Jun Kimura, Various aspects of F-wave values in healthy population, Clinical Neurophysiology 115 (2004) 2336-2342.



**NIHON KOHDEN ME COURSE
The New Neuro Series**

**Basic Nerve Conduction
Studies with Neuropack**


For a deeper inside in nerve conduction recording you can order the NIHON KOHDEN book about Basic Nerve Con-duction Studies free of charge (supply code: SP58-006 '09.02.NV.B). Please contact us or one of our representatives.

3.1. Stimulation Electrodes

Facts

Attaching the electrodes

1. Connect the electrode connector to the stimulation output connector of the device.
2. For stable measurement, rub the sites with a cotton pad moistened with alcohol or rub them with skinPure skin preparation gel (F020).
3. Soak the felt tips in water or physiological salt solution.
4. Push the two soaked felt tips into each pad holder.
5. Press the electrode against the skin site and fasten it with a provided band. The side with the black mark is the negative electrode.



Cleaning

1. If you leave the electrode after use without cleaning, it will rust and the performance of the electrode will deteriorate. Clean the electrode immediately after use.
2. Remove the felt tips from the electrode and wash the electrode with lukewarm water.
3. Do not wash with hot water or chlorine bleach!
4. Thoroughly wipe off moisture and keep the electrode dry and clean.

Supply Code	Model	Description
H635	NM-410S	Reusable Surface Stimulation Electrode, child, electrode: felt, electrode difference: 10 mm, size: 4 mm dia, incl. 10 felt tips type H647 & 2 Velcro bands for child foot/hand, 5-PIN DIN round connector, cable length: 2.0 m Foot band: width: 7 mm, length: 58 cm Hand band: width: 7 mm, length: 28 cm Note: YZ-0172 is required for MEB-9200/9400. For current products: MEB-9400/2300 Packing unit: 1 kit R



H636	NM-420S	Reusable Surface Stimulation Electrode, adult, electrode: felt, electrode difference: 23 mm, size: 8 mm dia, incl. 10 felt tips type H646 & 2 Velcro bands for adult foot/hand, 5-PIN DIN round connector, cable length: 2.0 m Foot band: width: 1.5 cm, length: 78 cm Hand band: width: 1.5 cm, length: 38 cm Note: YZ-0172 is required for MEB-9200/9400. For current products: MEB-9400/2300 Packing unit: 1 kit R
------	---------	---



H646	6143-009402	Felt Tips for H636, large For current products: MEB-9400/2300 Packing unit: 100 pcs
------	-------------	---



H647	6143-009411	Felt Tips for H635, small For current products: MEB-9400/2300 Packing unit: 100 pcs
------	-------------	---



2114-080955		Reusable Metal Tips for Stimulation Electrode H636, for long-term stimulation in OR, alternative to H646 For current products: MEB-9400/2300 Packing unit: 1 pcs R
-------------	--	---



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Facts

Attaching the electrodes

1. Connect the electrode connector to the stimulation output connector of the device.
2. For stable measurement, rub the sites with a cotton pad moistened with alcohol or rub them with skinPure skin preparation gel (F020).
3. Apply a small amount of Elefix EEG conduction paste (F509 or F510) to the tip of the electrode.
4. Press the electrodes against the skin sites and fasten it with a provided band. The side with the black mark is the negative electrode.

Cleaning

1. If you leave the electrodes after use without cleaning, the paste will harden and the performance of the electrode will deteriorate.
2. Wipe the electrodes with a clean cloth immediately after use.

Supply Code	Model	Description
H639	NM-422B	Reusable Surface Stimulation Electrode, electrode: stainless steel, electrode difference: 23 mm, cable length: 2.0 m, Velcro bands for adult foot/hand, band width: 1.5 cm, band length: 38 cm For current products: MEB-9400/2300 Packing unit: 1 pcs 



Supply Code	Model	Description
NKD-STL		Reusable Surface Stimulation Electrode, electrode: felt, electrode difference: 20 mm incl. 10 x disposable felt tips with size: 7 x 20 mm & 1 x elastic fixing band of 42 cm with 5-PIN DIN round connector, cable length: 1.5 m For current products: MEB-9400/2300 Packing unit: 1 kit 



NKD-STL-BAND		Reusable Elastic Fixing Band for NKD-STL, 5-PIN DIN round connector, length: 42 cm, width: 2 cm For current products: MEB-9400/2300 Packing unit: 1 pcs 
--------------	--	---



NKD-STLF		Felt Tips for NKD-STL, size: 7 x 20 mm For current products: MEB-9400/2300 Packing unit: 10 pcs
----------	--	---



3.2. Finger Electrodes

Facts

Finger electrodes are used for stimulation and recording of sensory nerves in finger or toe to detect e.g. the carpal tunnel syndrome (CTS).

1. When used as a stimulation electrode:

Connect the electrode connectors to the stimulation output connectors of the device. For H653 (NM-450S) the adapter YZ-0172 is necessary.

When used as a recording electrode:

Connect the electrode connectors to the electrode lead jacks on the electrode junction box or the extension cable.

- For stable measurement, rub the attachment sites with a cotton pad moistened with alcohol or skinPure skin preparation gel (F020).
- Apply an optimal amount of Elefix EEG conduction paste (F509 or F510) to the finger and put the finger through the electrodes. The black part is the negative electrode. Attach the black part at the base of the finger.

Cleaning

- If you leave the electrode after use without cleaning, the paste will harden and the performance of the electrode will deteriorate.
- Wipe the electrode with a clean cloth immediately after use.

Supply Code	Model	Description
H653	NM-450S	<p>Reusable Finger Electrode, stainless steel, 5-PIN DIN round connector, max. loop length: 25 mm, cable length: 2.0 m Note: Suitable as recording and stimulation electrode! For current products: MEB-9400/2300 Packing unit: 1 pair </p>
H659	NM-451B	<p>Reusable Finger Electrode for BM-230B, stainless steel, 2x 1.5 mm touch proof connector, max. loop length: 140 mm, cable length: 0.3 m Note: Suitable as recording and stimulation electrode! For current products: MEB-2300 Packing unit: 1 pair </p>
YZ-0172		<p>Adapter for Stimulation/Finger Electrodes, 5-PIN DIN connector (female) to stimulation output with 8-PIN DIN connector (male) For current products: MEB-9400/2300 Packing unit: 1 pcs </p>




3.3. Surface Electrodes

Surface electrodes are used to observe the general surface muscle activation instead of single muscle fibres. A proper skin preparation is essential for good signals and no artefacts. The electrodes' position for the defined muscle should be determined with caution as the location is important to obtain accurate signals. If the electrodes are located on a nearby muscle than the targeted one, the EMG signal could show up a mixed action potential picture. The direction of the electrodes should be set with the direction of the muscle fibres and the reference electrode needs to be placed on electrically neutral position.

3.3.1. Disposable Surface Electrodes

Supply Code	Model	Description
H690	NM-317Y3	Disposable NCS Electrodes, 2 recording electrodes (1 x black, 1 x red, size of electrodes: 14 x 25 mm) & 1 ground electrode (green, size of electrode: Phi=28 mm), cable length: 0.2 m, 1.5 mm touch proof connector, cable length: 0.2 m, requires K625A/K625B/BM-230B For current products: MEB-9400/2300 Packing unit: 10 x 3 pcs 



H691	NM-319Y	Disposable NCS Electrodes, 4 recording electrodes (2 x black, 2 x red, size of electrodes: 14 x 25 mm), cable length: 0.2 m, 1.5 mm touch proof connector, requires K625A/K625B/BM-230B For current products: MEB-9400/2300 Packing unit: 10 x 4 pcs 
------	---------	--




H692	NM-316Y	Disposable NCS Electrodes, 4 recording electrodes (2 x black, 2 x red, size of electrodes: 14 x 25 mm) & 1 x ground electrode (green, size of electrode: Phi=28 mm), cable length: 0.2 m, 1.5 mm touch proof connector, requires K625A/K625B/BM-230B For current products: MEB-9400/2300 Packing unit: 10 x 5 pcs 
------	---------	---



Facts

Attaching the electrodes

- To reduce the skin impedance, clean the electrode sites with a cotton pad moistened with alcohol or rub them with skinPure skin preparation gel (F020). Thoroughly dry the skin with dry gauze.
- Thoroughly dry the skin with a clean cotton pad. Residual sweat, moisture or skinPure will reduce the adhesive properties of the electrode.
- Write down the attachment sites on the tag of the electrode and attach the tag to the lead near the connector to distinguish the electrode.
- Remove the electrode from the backing sheet. Avoid touching the adhesive surface.
- Press one edge of the electrode down on the proposed measuring site and then continue to apply the rest of the electrode in the same manner.
- Press the electrode firmly with a finger for secure contact.
- Connect the electrode lead to the connection cable.



Supply Code	Model	Description
H694A	NM-314YS	Disposable NCS Electrodes, 4 MEP/SEP electrodes, size: 27 x 61 mm, black/grey, 1.5 mm touch proof connector, cable length: 0.2 m, requires K625A/K625B/K629A Note: Intended use is recording! Connector: black: - pole, beige: + pole For current products: MEB-9400/2300; MEE-1000 Packing unit: 10 x 4 pcs 



H694B	NM-314YL	Disposable NCS Electrodes, 4 MEP/SEP electrodes, size: 27 x 61 mm, black/grey, 1.5 mm touch proof connector, cable length: 1.8 m, requires K625A/K625B/K629A Note: Intended use is recording! Connector: black: - pole, beige: + pole For current products: MEB-9400/2300; MEE-1000 Packing unit: 10 x 4 pcs 
-------	----------	---

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description	
NKD-NCSD1		Ambu Disposable NCS & EP Electrodes, Ag/AgCl, size: 15 x 20 mm, 1.5 mm touch proof connector, cable length: 0.1 m For current products: MEB-9400/2300 Packing unit: 12 pcs 	
NKD-NCSD2		Ambu Disposable NCS & EP Electrodes, Ag/AgCl, size: 15 x 20 mm, 1.5 mm touch proof connector, cable length: 0.5 m For current products: MEB-9400/2300 Packing unit: 12 pcs 	
K625A	BM-001B	NCS Connection Cable for NCS Electrodes H69x, electrode side: 3 x 1.5 mm touch proof connector, device side: 6-PIN DIN round connector, cable length: 1.5 m For current products: MEB-9400/2300 Packing unit: 1 pcs 	
K625B	BM-002B	NCS Connection Cable for NCS Electrodes H69x, electrode side: 2 x 1.5 mm touch proof connector, device side: 6-PIN DIN round connector, cable length: 1.5 m For current products: MEB-9400/2300 Packing unit: 1 pcs 	
K625C	BM-003B	NCS Connection Cable for NCS Electrodes H69x, electrode side: 3 x 1.5 mm touch proof connector, device side: 6-PIN DIN round connector, cable length: 3.5 m For current products: MEB-9400/2300 Packing unit: 1 pcs 	
K625D	BM-004B	NCS Connection Cable for NCS Electrodes H69x, electrode side: 2 x 1.5 mm touch proof connector, device side: 6-PIN DIN round connector, cable length: 3.5 m For current products: MEB-9400/2300 Packing unit: 1 pcs 	
K629A	BM-990B	Stimulation Extension Cable for NCS Electrodes H69x, electrode side: 2 x 1.5 mm touch proof connector, device side: 5-PIN DIN round connector, cable length: 3.5 m For current products: MEB-9400/2300 Packing unit: 1 pcs 	
K636	BM-005B	Extension Cable for 2 Bipolar Electrodes, 2 signal parallel cable for epidural catheter, 2 x 1.5 mm touch proof connector male type and 2 x 1.5 mm touch proof connector female type, cable length: 3.5 m For current products: MEE-1000; MEB-2300 Packing unit: 1 pcs 	

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!


 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

3.3.2. Reusable Surface Electrodes

Facts

Attaching the electrodes

1. Connect the electrode leads to the corresponding input jacks.
2. To reduce impedance, clean the electrode sites with a cotton pad moistened with alcohol or rub them with skinPure skin preparation gel (F020).
3. Apply an optimal amount of cardioCream (F010), electrolyte cream, to the cleaned sites and rub the paste around the site in a circle not larger than 10 mm in diameter.
4. Apply an optimal amount of Elefix EEG conduction paste (F509 or F510) to the disk (cup) of the electrode and attach it on the pasted skin sites as shown below. The black part is the negative electrode.
5. Cover the electrode with a piece of cotton and gently press the edge of the cotton. Fasten the electrode on the sites with surgical tape if necessary.



Supply Code	Model	Description
H634	NM-312S	Reusable Surface Electrodes, Ag/AgCl, size: 1 cm dia, 6-PIN DIN round connector, cable length: 1.43 m For current products: MEB-9400/2300 Packing unit: 1 pair R



H638	NM-315S	Reusable Surface Electrodes, Ag/AgCl, size: 1 cm dia, 6-PIN DIN round connector, cable length: 2.3 m For current products: MEB-9400/2300 Packing unit: 1 pair R
------	---------	--

Supply Code	Model	Description
K637	BM-312B	Extension Cable for Surface Electrode H634, shielded, electrode side: 2 x alligator clip, device side: 6-PIN DIN connector, cable length: 3.5 m For current products: MEB-9400/2300 Packing unit: 1 pcs R

K638A	BM-900B	Adapter Cable for H634/K611, electrode side: 6-PIN DIN connector, device side: 2 x 1.5 mm touch proof connector, cable length: 0.3 m For current products: MEB-9400/2300 Packing unit: 1 pcs R
-------	---------	---

K638B	BM-901B	Adapter Cable for H634/K611, electrode side: 6-PIN DIN connector, device side: 2 x 1.5 mm touch proof connector, cable length: 0.3 m For current products: MEE-1000 Packing unit: 1 pcs R
-------	---------	--

D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

4. EMG Ground Electrodes

The ground electrode is a reference electrode to the recording electrode. Usually this electrode is placed away from the recording electrode on a neutral position. During electrical stimulation however, the ground electrode is located between the stimulating and recording electrode to reduce artefacts. For acoustic or visual stimulation the grounding is usually done on the forehead. It is important to have the most possible skin-electrode contact. Therefore a large electrode is preferred.

4.1. Disposable EMG Ground Electrodes

4.2. Reusable EMG Ground Electrodes

4.1. Disposable EMG Ground Electrodes

Facts

Attaching the electrodes

1. To reduce the skin impedance, clean the electrode sites with a cotton pad moistened with alcohol or rub them with skinPure skin preparation gel (F020). Thoroughly dry the skin with dry gauze.
2. Thoroughly dry the skin with a clean cotton pad. Residual sweat, moisture or skinPure will reduce the adhesive properties of the electrode.
3. Write down the attachment sites on the tag of the electrode and attach the tag to the lead near the connector to distinguish the electrode.
4. Remove the electrode from the backing sheet. Avoid touching the adhesive surface.
5. Press one edge of the electrode down on the proposed measuring site and then continue to apply the rest of the electrode in the same manner.
6. Press the electrode firmly with a finger for secure contact.
7. Connect the electrode lead to the connection cable.

Supply Code	Model	Description
H693	NM-310Y	Large Disposable NCS Ground Electrode, size: 42 x 84 mm, grey, 1.5 mm touch proof connector, cable length: 1.8 m, requires K625A/K625B/K629A Note: Intended use is grounding! For current products: MEB-9400/2300; MEE-1000 Packing unit: 10 pcs 



4.2. Reusable EMG Ground Electrodes

Supply Code	Model	Description
H658	NM-550B	Reusable Metal Plate Ground Electrode, stainless steel, size: 30 mm dia, 1.5 mm touch proof connector, cable length: 1.5 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs 



H658A	NM-551B	Reusable Metal Plate Ground Electrode, stainless steel, size: 30 mm dia, 1.5 mm touch proof connector, cable length: 0.3 m, for BM-230B For current products: MEB-2300 Packing unit: 1 pcs 
-------	---------	--



H662	NM-501B	Reusable Ground Electrode, adult, wrist/ankle, keyhole connector, incl. 1 x cable 1.5 m & 1 x ground Velcro strap 45 cm For current products: MEB-9400 /2300 Packing unit: 1 pcs 
------	---------	--



H662A	NM-506B	Reusable Ground Electrode, adult, wrist/ankle, keyhole connector, incl. 1 x cable 30 cm & 1 x ground Velcro strap 45 cm, for BM-230B For current products: MEB-2300 Packing unit: 1 pcs 
-------	---------	---



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
H663	NM-502B	Reusable Ground Electrode, adult, thigh, keyhole connector, incl. 1 x cable 1.5 m & 1 x ground Velcro strap 75 cm For current products: MEB-9400/2300 Packing unit: 1 pcs 
H664	NM-503B	Reusable Ground Electrode, child, wrist, keyhole connector, incl. 1 x cable 1.5 m & 1 x ground Velcro strap 31 cm For current products: MEB-9400 /2300 Packing unit: 1 pcs 
H664A	NM-505B	Reusable Ground Electrode, child, wrist, keyhole connector, incl. 1 x cable 30 cm & 1 x ground Velcro strap 31 cm, for BM-230B For current products: MEB-2300 Packing unit: 1 pcs 
H657	NM-540B	Reusable Metal Plate Ground Electrode, stainless steel, size: 30 mm dia, 2 mm connector, cable length: 1.5 m For old products: MEB-4200/5500/9100 Packing unit: 1 pcs 
H654	NM-511S	Reusable Ground Electrode, child, wrist, size: 4 mm dia connector, incl. 1 x cable 1.5 m & 1 x ground Velcro strap 31 cm For old products: MEB-4200/5200/5300/7102 Packing unit: 1 pcs 
H655	NM-522S	Reusable Ground Electrode, adult, arm/leg, size: 4 mm dia connector, incl. 1 x cable 1.5 m & 1 x ground Velcro strap 45 cm For old products: MEB-4200/5200/5300/7102 Packing unit: 1 pcs 
H656	NM-531S	Reusable Ground Electrode, adult, thigh, size: 4 mm dia connector, incl. 1 x cable 1.5 m & 1 x ground Velcro strap 75 cm For old products: MEB-4200/5200/5300/7102 Packing unit: 1 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

5. Additional EMG Accessories

5.1. Acoustic Stimulation

5.2. Visual Stimulation

5.3. Temperature Measurement

5.4. Pastes & Creams

5.5. Recording Paper

5.6. Cables

5.7. Adapter


5.1. Acoustic Stimulation

Supply Code	Model	Description
Y849A	DR-531B-14	Headphone, adult, 5-PIN XLR connector, cable length: 3.0 m, 2-channel, band width: 160 mm, pad size: 117 x 92 mm For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs R
Y849B	DR-531B-15	Headphone, child, 5-PIN XLR connector, cable length: 3.0 m, 2-channel, band width: 110 mm, pad size: 70 x 70 mm For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs R
Y848A	DR-531B-10B	Headphone, adult, 4-PIN XLR connector, cable length: 3.0 m, 2-channel, band width: 160 mm, pad size: 117 x 92 mm For old products: MEB-2200/5500/7202 Packing unit: 1 pcs R
Y849C	RP-86004	Ear Pads for Headphone Y848A/Y849A, adult, size: 11.7 x 9.2 cm For current products: MEB-9400/2300; MEE-1000 For old products: MEB-2200 Packing unit: 2 pcs R
Y849D		Ear Pads for Headphone Y848A/Y849A, child, size: 8.0 cm dia For current products: MEB-9400/2300; MEE-1000 For old products: MEB-2200 Packing unit: 2 pcs R
HEAD-BAND/C		Headband for Headphone Y848A/Y849A, child, size: Phi=118mm For current products: MEB-9400/2300; MEE-1000 For old products: MEB-2200 Packing unit: 1 pcs R
HEAD BAND		Headband for Headphone Y848A/Y849A, adult, size: Phi=160 mm For current products: MEB-9400/2300; MEE-1000 For old products: MEB-2200 Packing unit: 1 pcs R
YE-103J		Earphones, 5-PIN DIN round connector, cable length: 3.0 m, supplied with 2 ear adapters in 5 different sizes each For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs R



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description	
YE-102J		Earphones, child/adult, 4-PIN XLR connector, cable length: 3.0 m, supplied with 2 ear adapters in 5 different sizes each For old products: MEB-2200/4104/4200/5200/5300/5500/7103 Packing unit: 1 pcs R	
K635A	BM-701B	Extension Cable for YE-103J, cable length: 3.5 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs R	
K635B	BM-702B	Extension Cable for LS-102J, cable length: 3.5 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs R	
K101	RC-86014-H5P	Extension Cable for Headphone Y849A/B, 5-PIN connector, cable length: 3 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs R	
K101B	RC-86014-R4C	Extension Cable for Headphone Y848A, 4-PIN connector, cable length: 3 m For old products: MEB-2200 Packing unit: 1 pcs R	
Y855A	YZ-002A3	Ear Pieces, large, for YE-102J/103J, requires Y855E For current products: MEB-2300; MEE-1000 Packing unit: 20 pcs D	
Y855B	YZ-002A4	Ear Pieces, medium, for YE-102J/103J, requires Y855E For current products: MEB-2300; MEE-1000 Packing unit: 20 pcs D	
Y855C	YZ-002A5	Ear Pieces, small, for YE-102J/103J, requires Y855E For current products: MEB-2300; MEE-1000 Packing unit: 20 pcs D	

D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
Y855D	YZ-002A6	Ear Pieces, extra small, for YE-102J/103J, requires Y855E For current products: MEB-2300; MEE-1000 Packing unit: 20 pcs 



Y855E	YZ-002A7	Sound Conduction Tubes, for YE-102J/103J, necessary for ear pieces Y855x For current products: MEB-2300; MEE-1000 Packing unit: 10 pcs 
-------	----------	--



5.2. Visual Stimulation

Supply Code	Model	Description
LS-102J		LED Goggle, 4-PIN DIN connector, cable length: 3.0 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs 



LS-101J		LED Goggle, 7-PIN DIN connector, cable length: 3.0 m For old products: MEB-2200/4200/5500 Packing unit: 1 pcs 
---------	--	---



5.3. Temperature Measurement

Supply Code	Model	Description
NKD-P242B	YSI-409B	Body Surface Temperature Probe, 6.3 mm jack connector, disk size: 9.5 mm dia, cable length: 3.5 m Note: Use NKD-TEMPHEARTS4009/4012 to fix this probe! For current products: all BSM; PVM-2703; TEC-8300; MEB-9200/2300 Packing unit: 1 pcs 



NKD-TEMPHEARTS4009		Tempheart Adhesive Pads, to fix temperature probes, size: 1.5 cm dia For current products: all BSM; PVM-2703; TEC-8300; MEB-9200/2300 Packing unit: 100 pcs 
--------------------	--	---



NKD-TEMPHEARTS4012		Mini-Tempheart Adhesive Pads, to fix temperature probes, size: 1.13 cm dia For current products: all BSM; PVM-2703; TEC-8300; MEB-9200/2300 Packing unit: 100 pcs 
--------------------	--	---



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

5.4. Pastes & Creams

Supply Code	Model	Description
F010	Z-101BC	cardioCream ECG Electrode Paste, 100 g tube Packing unit: 2 pcs 
<i>cardioCream</i>		



F020	YZ-0019	skinPure Skin Preparation Gel, reduces the electrode skin impedance, 135 g tube Note: Can be used with all surface applied electrodes. Packing unit: 2 pcs 
<i>skinPure</i>		



Facts

The NIHON KOHDEN electrode paste Elefix integrates outstanding sticking-qualities with an excellent electrical conductivity. Additionally Elefix has a creamy quality so that the electrodes can be attached firmly to the patient without additional fixation, on the other hand Elefix does not harden so that the electrodes can be removed without any problems even after longer examinations. The contact impedance between skin and electrode is greatly reduced. In addition, the signal transmission remains stable even under unfavourable external conditions. NIHON KOHDEN Elefix can be used with all surface applied electrodes.

Supply Code	Model	Description
F509	Z-181JE	Elefix EEG Conduction Paste, 180 g tube Packing unit: 1 pcs 

Elefix



F510	Z-401CE	Elefix EEG Conduction Paste, 400 g cup Packing unit: 1 pcs 
------	---------	---

Elefix



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

5.5. Recording Paper

Supply Code	Model	Description
A720B	FQW210-10B	Thermal Paper, Z-fold, 210 mm x 100 m (page 300 mm), 333 pages, without grid, with queue mark/perforation For old products: MEB-4200; MEM-4104/4200 Packing unit: 5 pcs 



A725	FQW110-3-150	Recording Paper, Z-fold, 110 mm x 30 m (page 150 mm), 200 pages, without grid/queue mark/perforation For old products: MEB-4300/5304/5500/7102/7202; MEM-7202 Packing unit: 10 pcs 
------	--------------	--

A726	FQW110-2-150	Thermal Paper, Z-fold, 110 mm x 20 m (page 150 mm), 133 pages, without grid, with queue mark/perforation For old products: MEB-4300/5500 Packing unit: 10 pcs 
------	--------------	---

NKD-110-150	FQW110-3-150	Thermal Paper, Z-fold, 110 mm x 30 m (page 150 mm), 300 pages, without grid, with queue mark/perforation For old products: MEB-5300/7102; MEM-7102/7112 Packing unit: 1 pcs 
-------------	--------------	---


NKD-210-10/295		Thermal Paper, Z-fold, 210 mm x 100 m (page 295 mm), 339 pages, without grid, with queue mark/perforation, equivalent to A052 & NKD-210-10A For current products: ECG-1550/9320 For old products: MEB-4200; MEM-4100/4200 Packing unit: 10 pcs 
----------------	--	--

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

5.6. Cables

Supply Code	Model	Description
K617	2114-064171	EMG Jumper Cable for H851, 2 mm dia 1-PIN connector, unshielded, cable length: 0.15 m For old products: MEB-4200/5200/5300/7102 Packing unit: 5 pcs 
K617A	NE-420B	EMG Jumper Cable for H852, 2 mm dia 2-PIN connector, unshielded, cable length: 0.15 m For old products: MEB-4200/5200/5300/7102 Packing unit: 5 pcs 
K617B	NE-430B	EMG Jumper Cable for H851A, 1-PIN DIN round connector, unshielded, cable length: 0.15 m For current products: MEB-9400/2300 Packing unit: 5 pcs 
K617C	NE-432B	EP Jumper Cable for H852A/H854, 2-PIN DIN keyhole connector, unshielded, cable length: 0.15 m For current products: MEB-9400/2300 Packing unit: 5 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
K297A	EX212RE20	Extension Cable for Electrodes, shielded, 1.5 mm touch proof connector, colour: red, cable length: 2.0 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs 
K297B	EX212YE20	Extension Cable for Electrodes, shielded, 1.5 mm touch proof connector, colour: yellow, cable length: 2.0 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs 
K297C	EX212GR20	Extension Cable for Electrodes, shielded, 1.5 mm touch proof connector, colour: green, cable length: 2.0 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs 
K297D	EX212WH20	Extension Cable for Electrodes, shielded, 1.5 mm touch proof connector, colour: white, cable length: 2.0 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs 
K297E	EX212BK20	Extension Cable for Electrodes, shielded, 1.5 mm touch proof connector, colour: black, cable length: 2.0 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs 
K297F	EX212BW20	Extension Cable for Electrodes, shielded, 1.5 mm touch proof connector, colour: brown, cable length: 2.0 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs 
K297G	EX212PR20	Extension Cable for Electrodes, shielded, 1.5 mm touch proof connector, colour: purple, cable length: 2.0 m For current products: MEB-9400/2300; MEE-1000 Packing unit: 1 pcs 
K637	BM-312B	Extension Cable for Surface Electrode H634, shielded, electrode side: 2 x alligator clip, device side: 6-PIN DIN connector, cable length: 3.5 m For current products: MEB-9400/2300 Packing unit: 1 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
K621A	BM-503B	Extension Cable for ECG Electrodes, unshielded, electrode side: alligator clip, device side: 1-PIN DIN connector, colours: red/green/yellow, cable length: 3.0 m For current products: MEE-1000 Packing unit: 3 pcs 

K622A	BM-504B	Extension Cable for ECG Electrodes, shielded, electrode side: alligator clip, device side: 2-PIN DIN connector (keyhole), colours: red/green/yellow, cable length: 3.0 m For current products: MEB-9400/2300 Packing unit: 3 pcs 
-------	---------	--



K634	BM-121B	Extension Cable for MS-120B, unshielded, multi-contact connector, cable length: 5.0 m For current products: MEE-1000 Packing unit: 1 pcs 
------	---------	--



K626	BM-910B	EP Connection Cable for H854, shielded, 9-PIN DIN round connector, cable length: 3.0 m For current products: MEB-9400 Packing unit: 1 pcs 
------	---------	--



K627	BM-911B	ARB Extension Cable, shielded, 9-PIN DIN round connector, cable length: 3.0 m For current products: MEB-9400 Packing unit: 1 pcs 
------	---------	--



5.7. Adapter

Supply Code	Model	Description
YZ-0172		Adapter for Stimulation/Finger Electrodes, 5-PIN DIN connector (female) to stimulation output with 8-PIN DIN connector (male) Note: YZ-0172 is required for MEB-9200/9400. For current products: MEB-9400/2300 Packing unit: 1 pcs 
Y085	6144-007171	Adapter with 2 mm Socket and 1.5 mm touch proof connector, to connect electrode cables with 2 mm connector to input-boxes with 1.5 mm touch proof type sockets For current products: EEG-1200/9100; MEB-9400/2300; MEE-1000 Packing unit: 5 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!





II. Accessories & Consumables for Electroencephalography (EEG)

Electroencephalographic diagnostics in hospitals and medical practices place high demands on personnel and technology. From daily routines to highly specialized examinations in long-term monitoring, polysomnography or the intensive care unit, EEG systems are required which satisfy the needs of the individual.

An Electroencephalogram (EEG) is a brain wave test that detects abnormalities in the electrical activity of the brain. This helps to diagnose: seizure disorders, head injuries, sleep disorders, unconsciousness reasons and further neurological problems. Thanks to our extensive experience, we are able to offer you the ideal solution for hospitals and medical practices. Our high quality approved EEG accessories and consumables help you to obtain best signals and to optimize your EEG recording.

Determination of the electrodes position (10/20 system)

1. Divide the longitudinal line of the head into halves and attach an electrode at the Cz point (Fig. 1).
2. Divide the distance between Cz and the nasion in proportions as shown in Fig. 1 and attach an electrode at Fz and Pz. Do not attach an electrode at Fpz and Oz.
3. Divide the transverse line of the head into proportions of 10%, 20%, 20%, 20%, 20% and 10% as shown in Fig. 2 and attach electrodes at T3, T4, C3 and C4.
4. Divide the peripheral line passing over Fpz, T4, Oz and T3 into proportions shown in Fig. 3 and attach electrodes at Fp2, F8, T6, O2, Fp1, F7, T5 and O1.
5. Take the line passing over F7, Fz and F8 and attach an electrode at the mid-point (F3) between F7 and Fz as shown in Fig. 4. Attach another electrode at the mid-point (F4) between F8 and Fz.
6. Take the line passing over T5, Pz and T6 and attach an electrode at the mid-point (P3) between T5 and Pz. Attach another electrode at the mid-point (P4) between T6 and Pz.



6. Accessory Kit for EEG

Supply Code **Model**
NKD-BE-413AG

Description

EEG Accessory Kit – touch proof on electrode side

Kit contains: 1 universal EEG Cap (H564A) incl. 1 chin strap (H551A), 1 package EEG bridge electrodes (H542B) incl. related connection cables (K518A), 1 package of felt pads (NKD-FILZ), 1 package of O-rings (H659), 3 earlobe electrodes (H543B), 1 package of 4 Fastclip electrodes (H068A), 2 Fastclip connection cables (K512A), 2 tubes of cardio-Cream ECG electrode paste (F010) and 2 tubes of skinPure skin preparation gel (F020)

Packing unit: 1 kit



H564A



H551A



H542B



K518A



NKD-FILZ



H659



H543B



H068A



K512A



F010



F020

🗑 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

🔄 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

7. Bridge Electrodes

Facts

The high quality, sintered EEG bridge electrode is adjustable in height through simply rotating. The easy handling and cleaning makes this electrode a convenient daily work electrode.

Supply Code	Model	Description
H542B		Sintered EEG Bridge Electrode, Ag/AgCl, colour: beige, flexible, with cross hole and bar-shaped stand, 2 mm connector For current products: EEG-1200/9100 Packing unit: 23 pcs 



H542C		EEG Bridge Electrode, Ag, colour: white, flexible, round stand, 2 mm connector For current products: EEG-1200/9100 Packing unit: 23 pcs 
-------	--	---



NKD-FILZ		Electrode Felt Pad Kit for H542B, size: 25 mm dia, fixation with O-Ring Kit H569 For current products: EEG-1200/9100 Packing unit: 100 pcs 
----------	--	--



H561		EEG Electrode Pad Kit for H542C, incl. 100 felt tips, 100 felt pads & 50 O-Rings For current products: EEG-1200/9100 Packing unit: 1 kit 
------	--	--



NKD-O-RING		O-Ring Kit for Sintered EEG Bridge Electrode H542B, size: 8.5 mm dia internal diameter For current products: EEG-1200/9100 Packing unit: 50 pcs 
------------	--	---




H569		O-Ring Kit for Sintered EEG Bridge Electrode H542B, size: 1.3 cm dia internal diameter For current products: EEG-1200/9100 Packing unit: 50 pcs 
------	--	---



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
K518A		Electrode Connection Cable for Bridge Electrodes H542B/H543B, electrode side: 2 mm, input box side: 1.5 mm touch proof connector, cable length: 0.9 m For current products: EEG-1200/9100 Packing unit: 25 pcs 
K518B		Electrode Connection Cable for Bridge Electrodes H542B/H543B, electrode side: 2 mm touch proof connector, input box side: 1.5 mm touch proof connector, cable length: 0.9 m For current products: EEG-1200/9100 Packing unit: 25 pcs 
K516A	BE-403B	EEG Electrode Cable, patient side: alligator clip, input box side: 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100 Packing unit: 30 pcs 
K516B		Electrode Connection Cable for Bridge Electrodes H542B/H543A, input box side: 2 mm connector, cable length: 0.8 m For old products: EEG-4400/7300 Packing unit: 25 pcs 



Facts

Disinfection of bridge electrodes

As indicated in the user manual of our device, we recommend a thorough cleaning after each recording. This is based on hygienic recommendations pertaining to sediments and salt deposit. Chlorine differences might cause polarization and increase impedances leading to a decrease in electrolyte conductivity. This negatively affects recording results.

For cleaning bridge electrodes rinse them with clear warm water and wipe them thoroughly with a soft cotton cloth. For disinfection use 70 % ethanol solution. Never soak the sintered electrodes in the solution but wipe them with soft gauze moistened with ethanol. Afterwards dry the electrodes thoroughly with a soft cloth to avoid scratches to the electrode surface that could cause artefacts while EEG recording.

The electrode connector must be kept clean and dry to obtain reliable measurement signals. For cleaning the electrode cable and connector make sure they are not connected to the bridge electrodes any more. Then start wiping with soft cloth moistened with water or soap. Do not use ethanol for cleaning that may cause deterioration of the coating. Thoroughly dry the cable and connector after cleaning.

For a new recording, moisten the felt pads with a saline solution of maximum 0.9% to ensure a constant contact resistance for EEG recording. Please remember that salt is generally aggressive to the electrodes' surface. Therefore we recommend cleaning the bridge electrodes thoroughly after the last recording of the day and previous to a non-usage time period. Especially take care of the contact plugs which may corrode by intensive contact with saline solutions. A long lifetime of a bridge electrode depends on the maintenance and care given to the electrodes.

By providing good maintenance care to the electrodes, proper EEG recording should be able to be performed without difficulty. If there is indeed a limited functionality or persistent problem, the electrodes should be replaced completely in order to avoid recording artefacts based on different electrolyte conductivity. Please remove the felt pads for disinfection and exchange them regularly.

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

8. EEG Surface Electrodes

EEG surface electrodes are mainly used for EEG measurement because of the non-invasiveness and comfort for the patients. In contrary to minimal-invasive subdermal needle electrode, the electrodes impedances are higher. Therefore the preparation of the skin with skinPure skin preparation gel (F020) is very important to reduce the contact resistance between scalp and surface electrodes. To secure the electrodes generally some conductive cream such as Elefix (F509 or F510) and tape is used for disk electrodes. For restless patients it is however suggested to use cup electrodes with a more adhesive paste such as collodion (KOL-1) or Grass paste (NKD-EC2).

8.1. EEG Disk Electrodes


8.2. EEG Cup Electrodes

8.1. EEG Disk Electrodes

Facts

Attaching the electrodes

1. To reduce the skin impedance, clean the electrode sites with a cotton pad moistened with alcohol or rub them with skinPure skin preparation gel (F020). Thoroughly dry the skin with dry gauze.
2. Apply an optimal amount of Elefix EEG conduction paste (F509 or F510) to the cleaned sites and rub the paste around the sites in a circle not larger than 10 mm in diameter.
3. Apply Elefix EEG conduction paste about 1 mm thick to the disk (cup) of the electrode, and attach it on the pasted skin sites
4. Cover the electrode with a piece of cotton, and press the edge of the cotton. Fasten the electrode on the site with surgical tape if necessary.
5. Connect the electrode leads to the corresponding input jacks on the electrode junction box.



Supply Code	Model	Description
H503A	NE-113A	EEG Disk Electrode for EOG & EMG, Ag, size: 11 mm dia, with brown plastic covering, 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100; MEE-1000 Packing unit: 26 pcs 



H503B	NE-114A	EEG Disk Electrode for EOG & EMG, Ag, size: 11 mm dia, 1.5 mm touch proof connector, cable length: 0.8 m For current products: EEG-1200/9100; MEE-1000 Packing unit: 26 pcs 
-------	---------	---

H503	NE-103A	EEG Disk Electrode for EOG & EMG, Ag, size: 11 mm dia, 2 mm connector, cable length: 1.5 m For old products: EEG-4400/7300/8300 Packing unit: 26 pcs 
------	---------	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!





 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

8.2. EEG Cup Electrodes

Facts

Attaching the electrodes

1. Determine the electrode placement sites according to the 10/20 system.
2. Clean the area with an alcohol moistened cotton pad or skinPure skin preparation gel (F020) to remove oil from the skin. Wipe it afterwards with a dry cloth.
3. Thinly spread collodion paste (KOL-1) on electrode brim, scalp and hair roots.
4. Put Elefix conduction paste (F509 or F510) or cardioCream (F010) in a syringe and inject it into the hole of fastened electrode.
5. Cover the electrode with a small gauze pad and press the pad around the electrode. Attach the electrode with surgical tape or use the EEG net (H568) to secure the electrode position.



Supply Code	Model	Description
H526	NE-134A	EEG Cup Electrode (Collodion), Ag, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100; BSM-3000/6000/9100; PSG-1100; MEB-2300; MEE-1000 Packing unit: 12 pcs 
H527	NE-136A	EEG Cup Electrode (Collodion), Ag, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 0.7 m For current products: EEG-1200/9100; BSM-3000/6000/9100; PSG-1100; MEB-2300; MEE-1000 Packing unit: 12 pcs 
H520	NE-124A	EEG Cup Electrode (Collodion), Ag, size: 10 mm dia, 2 mm connector, cable length: 1.5 m For old products: EEG-4400/7300/8300 Packing unit: 12 pcs 
H521	NE-126A	EEG Cup Electrode (Collodion), Ag, size: 10 mm dia, 2 mm connector, cable length: 0.7 m For old products: EEG-4400/7300/8300 Packing unit: 12 pcs 






 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Facts**Attaching the electrodes**

1. Determine the electrode placement sites according the 10/20 system.
2. Clean the area with a alcohol moistened cotton pad or skinPure skin preparation gel (F020) to remove oil from the skin. Wipe it afterwards with a dry cloth.
3. Apply a small amount of Elefix EEG conduction paste (F509 or F510) or cardioCream (F010) in a 1 cm diameter circle on the cleaned skin. Avoid spreading the paste too thin or creating to large space.
4. Press the electrode gently down on the pasted skin area.
5. Cover the electrode with a small gauze pad and press the pad around the electrode. Attach the electrode with surgical tape or use the EEG net (H568) to secure the electrode position.



Supply Code	Model	Description
NKD-FE5GH		Grass Gold Cup Electrode, Au, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 10 pcs 

NKD-FE5GH-07		Grass Gold Cup Electrode, Au, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 0.7 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 10 pcs 
--------------	--	---

NKD-FE5GH-12		Grass Gold Cup Electrode, Au, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 1.2 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 10 pcs 
--------------	--	---

NKD-FE5GH-24		Grass Gold Cup Electrode, Au, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 2.4 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 10 pcs 
--------------	--	---

NKD-FE6GH60		Grass Gold Cup Electrode, Au, size: 6 mm dia, 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 10 pcs 
-------------	--	--

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

9. Subdermal Needle Electrodes

In comparison to surface electrodes, subdermal needle electrodes can be quickly and simply inserted subcutaneously in the patient's scalp without skin preparation or use of conductive paste. This rapid application possibility is especially advantageous when every second counts such as for emergency EEG on patients supposed to have serious acute brain dysfunction (ABD) because of head injuries and the threat of permanent brain damage or death. For a continuous EEG (cEEG) monitoring on ICU or on comatose patients the needle electrodes can be fixed with an adhesive paste and covered with a tape permeable to air to ensure an even better stay on place.

Supply Code	Model	Description
H537A	NE-224S	Reusable Subdermal EEG Needle Electrode, stainless steel, 1.5 mm touch proof connector, size: 0.22 x 22.5 mm, cable length: 1.5 m For current products: EEG-1200/9100; MEB-9400; MEE-1000 Packing unit: 20 pcs R



NE-110B		Disposable Subdermal Needle Electrode Prewired, single type, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 1.0 m For current products: EEG-1200/9100; MEE-1000 Packing unit: 25 pcs D
---------	--	---



NE-115B		Disposable Subdermal Needle Electrode Prewired, single type, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 1.5 m For current products: EEG-1200/9100; MEE-1000 Packing unit: 25 pcs D
---------	--	---

NE-120B		Disposable Subdermal Needle Electrode Prewired, single type, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 2.0 m For current products: EEG-1200/9100; MEE-1000 Packing unit: 25 pcs D
---------	--	---

NE-125B		Disposable Subdermal Needle Electrode Prewired, single type, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 2.5 m For current products: EEG-1200/9100; MEE-1000 Packing unit: 25 pcs D
---------	--	---

NE-215B		Disposable Subdermal Needle Electrode Prewired, twisted pair, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 1.5 m For current products: EEG-1200/9100; MEE-1000 Packing unit: 20 pcs D
---------	--	--



NE-220B		Disposable Subdermal Needle Electrode Prewired, twisted pair, stainless steel, 1.5 mm touch proof connector, size: 0.4 x 13 mm needle (27G), cable length: 2.0 m For current products: EEG-1200/9100; MEE-1000 Packing unit: 20 pcs D
---------	--	--

D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!


10. Earlobe Electrodes

During EEG measurement reference electrodes need to be placed on “neutral” locations. Besides any surface electrode, often located on the centre line, earlobes are used as references. To avoid electrical artefacts it is recommended to use the same material e.g. silver for the earlobe electrodes as the electrodes applied on the scalp.

Facts

Attaching the electrodes

1. Apply a small amount of Elefix EEG conduction paste (F509 or F510) to both earlobes. Gently clip on the earlobe electrodes.
2. Do not apply excessive pressure to the electrode because this may cause direct electrode-skin-contact (not through the paste) and generate unnecessary polarization voltage.
3. Secure the electrode cable to the neck with surgical tape to prevent the cable from movement.



Note: For EEG-1100/5500/5600/7400/9200 the adapter Y085 is required for earlobe electrodes with 2 mm connector.
For EEG-4400/7300/8300 the adapter Y085 is not required.

Supply Code	Model	Description	
H540A	NE-311A	EEG Earlobe Electrode Kit, Ag, incl. 1 x earlobe electrode clip & 1 x earlobe electrode, cable length: 1.5 m, 1.5 mm touch proof connector For current products: EEG-1200/9100 Packing unit: 1 kit R	
H540	NE-301A	EEG Earlobe Electrode Kit, Ag, incl. 1 x earlobe electrode clip & 1 x earlobe electrode, cable length: 1.5 m, 2 mm connector, requires Y085 For old products: EEG-1100/4400/5500/5600/7300/7400/8300/9200 Packing unit: 1 kit R	
H543B		Sintered EEG Earlobe Electrode Kit, Ag/AgCl, child/adult, incl. 1 x earlobe electrode clip & 1 x earlobe electrode, cable length: 0.9 m, 1.5 mm touch proof connector, requires Y085 For current products: EEG-1200/9100 Packing unit: 1 kit R	
H543A		Sintered EEG Earlobe Electrode Kit, Ag/AgCl, child/adult, incl. 1 x earlobe electrode clip & 1 x earlobe electrode, cable length: 0.9 m, 2 mm connector, requires Y085 For old products: EEG-1100/4400/5500/5600/7300/7400/8300/9200 Packing unit: 1 kit R	
H541	NE-301B	EEG Earlobe Electrode Clip for H540/H540A/H543B/H543A For current products: EEG-1200/9100 Packing unit: 10 pcs R	
Y085	6144-007171	Adapter with 2 mm Socket and 1.5 mm touch proof connector, to connect electrode cables with 2 mm connector to input boxes with 1.5 mm touch proof type sockets For current products: EEG-1200/9100; MEB-9400 Packing unit: 5 pcs R	

D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

11. ECG Measurement

During EEG measurement physiological and non-physiological caused artefacts could lead to a misinterpretation of high peaks. One of those artefacts is possible through ECG. Therefore, an ECG measurement in addition to the EEG recording is usually realized to identify those artefacts as well as heart cycle-related brain activities.

11.1. Clip Electrodes

11.2. EEG Electrodes

11.1. Clip Electrodes

Facts

The NIHON KOHDEN Fastclip ECG Electrodes allow a flexible application for routine work and emergency on every paediatric and adult patient thanks to a simple adjustability of the clip pressure (clip size). A more comfortable handling is guaranteed due to different connectors for any existing adapter type, a big contact area for reliable ECG measurements and easier cleaning after the measurement because of less asperity of the clip. In addition to this IEC coloured label sets are delivered with NIHON KOHDEN Fastclip ECG Electrodes.

By following the next steps, optimized derivations can be realized with NIHON KOHDEN Fastclip ECG Electrodes:

- 1.) Clean the according skin areas with alcohol in order to remove traces of dirt and oil on the skin.
- 2.) Apply a thin layer of cardioCream F010 on the cleaned skin areas and on the electrode part.
- 3.) Clip the electrode on the wrist or ankle. Adjust the electrode position by sliding the electrode part according to the thickness of the limb.

Adjustability

strong middle wide
weak narrow
spring

Contact pressure and position of electrodes section can be adjusted according to patients' conditions.

Easy to handle

Opens more widely*, and round form brings easy to grip.
(*compared to the previous model)

Easy cleaning

trench shape raised dots
Previous model New model

Less asperity, easy cleaning.

Three types of connectors available

Tip Clip Snap

Any of tip, clip and snap type connector is available.

Supply Code	Model	Description
H068A		Fastclip ECG Electrode, 3 mm banana plug, child/adult, adjustable size, colour: blue, coloured IEC labels: red, yellow, green, black For current products: all ECG; EEG-1200/9100 Packing unit: 4 pcs R

Fastclip



K512A	BC-112B	Connection Cable for Fastclip ECG Electrodes, 1.5 mm touch proof connector, cable length: 2.0 m For current products: EEG-1200/9100 Packing unit: 2 pcs R
-------	---------	--



K512	BC-102B	Connection Cable for Fastclip ECG Electrodes, 2 mm connector, cable length: 2.0 m For old products: EEG-1100/4400/5500/5600/7300/7400/8300/9200 Packing unit: 2 pcs R
------	---------	--



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
NKD-K10/C		<p>Clip-on Limb Electrode, AgCl, child, 4 mm banana plug, colours: red, yellow, green, black</p> <p>For current products: all ECG; EEG-1200/9100</p> <p>For old products: EEG-1100/4400/5500/5600/7300/7400/8300/9200</p> <p>Packing unit: 4 pcs </p>



NKD-KR4B2SS		<p>Connection Cable for NKD-K10/C, 2 mm connector, cable length: 1.5 m</p> <p>For current products: all ECG</p> <p>For old products: EEG-1100/4400/5500/5600/7300/7400/8300/9200</p> <p>Packing unit: 1 pcs </p>
-------------	--	--



NKD-KR4BDINS		<p>Connection Cable for NKD-K10/C, 1.5 mm touch proof connector, cable length: 1.5 m</p> <p>For current products: all ECG; EEG-1200/9100</p> <p>Packing unit: 1 pcs </p>
--------------	--	---



NKD-KR26DIN20		<p>Extension Cable, 1.5 mm DIN jack to 1.5 mm touch proof connector, cable length: 2.0 m</p> <p>For current products: EEG-1200/9100</p> <p>Packing unit: 1 pcs </p>
---------------	--	--



11.2. ECG Electrodes


NIHON KOHDEN Vitrode Disposable ECG Electrodes –

The result of over 30 years experience in disposable electrode manufacturing.

The ECG electrodes are the most important portions of the ECG/arrhythmia monitoring system. It is through these electrodes that the electrical signals are conducted from the patient and into the monitoring system. If the electrodes are left in place for an extended period of time (over 48 hours), or they are reused once they are removed, they do not adhere to the patient's skin and they cause "noisy" ECG signals because of the poor conduction. In addition, the patient's movement causes the electrodes to lose conducting capabilities intermittently, and this leads to a "wandering" ECG waveform and muscle artefact. The standard monitoring practice is to change these electrodes at least every 48 hours.

Please contact us for individual offer based on your annual volume! Use the possibility of call orders!

Supply Code	Model	Description
G203	L-150	ECG Electrode Vitrode L, Ag/AgCl, adult, size: 35 mm dia, foam tape, solid gel, general use Packing unit: 5 x 30 pcs 
		Vitrode L Vitrode
G210D	F-150M	ECG Electrode Vitrode F, Ag/AgCl, adult, size M: 25 x 45 mm, cotton tape especially for sensitive skin, adhesive gel Packing unit: 3 x 50 pcs 
		Vitrode F Vitrode
G210C	F-150S	ECG Electrode Vitrode F, Ag/AgCl, neonate/child, size S: 18 x 36 mm, cotton tape especially for sensitive skin, adhesive gel Packing unit: 3 x 50 pcs 
		Vitrode F Vitrode
NKD-15026N		ECG Snap Lead Wire Kit, 1.5 mm touch proof connector, cable length: 1.5 m, colours: red, yellow, green, black For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 4 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

12. EEG Head Caps

12.1. NIHON KHODEN Caps

12.2. ElectroCaps

12.3. BraiNet Caps

12.1. NIHON KOHDEN Caps

Facts

The NIHON KOHDEN universal EEG caps are usually used as positioning help for EEG bridge electrodes H542B and H542C. Silicone tubes are reticular arranged and reunited on both ears. The design of the ear pieces as well as the flexible material of the tubes allows an easy adjustment to different patients' scalp sizes and forms as the side and cross tubes can be simply put into the desired position. This straightforward and universal application tool is a perfect help for the classic 10/20 EEG recording as well as further EEG electrodes positioning systems.

Supply Code	Model	Description
H564A	BE-0007	Universal EEG Silicone Cap, adjustable For current products: EEG-1200/9100 Packing unit: 1 pcs 



H551A	BE-0008	Chin Strap for H564A, adjustable, with chin guard, strap length: 42 cm, incl. one-touch type adapter For current products: EEG-1200/9100 Packing unit: 1 pcs 
-------	---------	--



1144-008488C		Ear Holder for Universal EEG Silicone Cap H564A For current products: EEG-1200/9100 Packing unit: 1 pcs 
--------------	--	---



2114-073312		Tube Holder for Universal EEG Silicone Cap H564A For current products: EEG-1200/9100 Packing unit: 1 pcs 
-------------	--	--

2114-073401		Tube Stopper for Universal EEG Silicone Cap H564A For current products: EEG-1200/9100 Packing unit: 10 pcs 
-------------	--	--

BE-81005		Exchange Tube for Universal EEG Silicone Cap H564A For current products: EEG-1200/9100 Packing unit: 1 pcs 
----------	--	--



H568	2114-042676B	EEG Net, to secure EEG electrodes positions for long-term examinations For current products: EEG-1200/9100 Packing unit: 1 pcs 
------	--------------	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

12.2. ElectroCaps

Easy, fast and correct EEG electrodes application is guaranteed by the large portfolio of neonatal and adult cap sizes of ElectroCaps. Whereas for adults the 10/20 system for EEG electrodes positioning is implemented, the 10/10 system for neonates or infants up to two years is also possible. Generally adult and child caps are provided with 20 or 21 electrodes, neonatal or infant caps with 14 electrodes. Any further electrodes positions can be produced based on your request. To fix the caps special chest bands and chin straps were designed. The connection to your EEG system is facilitated by only one electrode board adapter for connecting all electrodes.

Please note that all ElectroCaps are made of spandex and need to be cleaned immediately after use. Otherwise the material elasticity will be destroyed resulting in a shorter product life time. For cleaning only use the cleaning emulsion Ivory to avoid film residual on the electrode material. It is recommended to not mix different colours together.

12.2.1. ElectroCap Systems

Facts

The ElectroCap System I, II and III are comfortable starter kits as they include different types of caps, that you can define, as well as all necessary accessories to start immediately the EEG application. Besides this a training video DVD is provided to show you how to apply the caps properly.

Note: Please specify with your order if you need a 2 mm connector for your ElectroCap System Kits!

Supply Code	Model	Description
NKD-SYS1		<p>ElectroCap System I Kit</p> <p>Kit contains: 1 medium cap at your choice, 1 electrode board adapter (NKD-E2), 1 chest band (NKD-E3M), 1 quick insert electrode (NKD-E4), 1 pair of ear electrodes (NKD-E5/9Z-26), 100 disposable foam rings (NKD-E6), 1 kit of needles and syringe (NKD-E7), 1 electrode gel cup (NKD-E9), 1 adult head measuring tape (NKD-E12), 1 bottle of cleaning emulsion (NKD-E16/700)</p> <p>For current products: EEG-1200/9100</p> <p>Packing unit: 1 kit</p>



NKD-E2



NKD-E3M



NKD-E4



NKD-E5/9Z-26



NKD-E6



NKD-E7



NKD-E9



NKD-E12



NKD-E16/700



1 medium cap at your choice

📦 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

🔄 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code **Model**
NKD-SYS2


Description

ElectroCap System II Kit

Kit contains: 1 large/1 medium cap at your choice, 1 electrode board adapter (NKD-E2), 1 chest band (NKD-E3M), 2 quick insert electrodes (NKD-E4), 2 pairs of ear electrodes (NKD-E5/9Z-26), 100 disposable foam rings (NKD-E6), 2 kits of needles and syringe (NKD-E7), 2 electrode gel cups (NKD-E9), 1 adult head measuring tape (NKD-E12), 1 bottle of cleaning emulsion (NKD-E16/700)

For current products: EEG-1200/9100

Packing unit: 1 kit



NKD-E2



NKD-E3M



NKD-E4



NKD-E5/9Z-26



NKD-E6



NKD-E7



NKD-E9



NKD-E12



NKD-E16/700



1 large/1 medium cap at your choice

ⓘ Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

Ⓜ Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code **Model**
NKD-SYS3

Description

ElectroCap System III Kit

Kit contains: 1 large/1 medium/1 small/1 extra small cap at your choice, 1 electrode board adapter (NKD-E2), 2 chest bands (NKD-E3M), 4 quick insert electrodes (NKD-E4), 2 pairs of ear electrodes (NKD-E5/9Z-26), 200 disposable foam rings (NKD-E6), 4 kits of needles and syringe (NKD-E7), 3 electrode gel cups (NKD-E9), 1 adult head measuring tape (NKD-E12), 1 bottle of cleaning emulsion (NKD-E16/700)

For current products: EEG-1200/9100

Packing unit: 1 kit



NKD-E2



NKD-E3M



NKD-E4



NKD-E5/9Z-26




NKD-E6



NKD-E7



NKD-E9



NKD-E12



NKD-E16/700



1 large/1 medium/1 small/
 1 extra small cap at your choice

ⓘ Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

Ⓜ Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

12.2.2.

ElectroCaps for Adults

Supply Code	Model	Description
NKD-E1L		ElectroCap L, adult, size: 58-62 cm, colour: blue, 20 electrodes For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-E1LM		ElectroCap L, adult, size: 56-60 cm, colour: blue & red, 20 electrodes For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-E1LMF		ElectroCap L, adult, size: 56-60 cm, colour: blue & red, 20 flat electrodes For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-E1M		ElectroCap M, adult, size: 54-58 cm, colour: red, 20 electrodes For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-E1M-F		ElectroCap M, adult, size: 54-58 cm, colour: red, 20 flat electrodes For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-E3S		Chest Band for all ElectroCaps, S, adult, colour: yellow, band length: 72 cm For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-E3M		Chest Band for all ElectroCaps, M, adult, colour: red, band length: 95 cm For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-E3L		Chest Band for all ElectroCaps, L, adult, colour: blue, band length: 1.16 m For current products: EEG-1200/9100 Packing unit: 1 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
NKD-E23/7		Chin Strap for all ElectroCaps, adults, colour: white, strap length: 17.78 cm For current products: EEG-1200/9100 Packing unit: 1 pair 



12.2.3.

ElectroCaps for Children/Adults

Supply Code	Model	Description
NKD-E1MS		ElectroCap M/S, child/adult, size: 52-56 cm, colour: red & yellow, 20 electrodes For current products: EEG-1200/9100 Packing unit: 1 pcs 

NKD-E1MS-F		ElectroCap M/S, child/adult, size: 52-56 cm, colour: red & yellow, 20 flat electrodes For current products: EEG-1200/9100 Packing unit: 1 pcs 
------------	--	---

NKD-E1S		ElectroCap S, child, size: 50-54 cm, colour: yellow, 20 electrodes For current products: EEG-1200/9100 Packing unit: 1 pcs 
---------	--	--

NKD-E1SXS		ElectroCap S, child, size: 48-52 cm, colour: yellow & green, 20 electrodes For current products: EEG-1200/9100 Packing unit: 1 pcs 
-----------	--	--



NKD-E1XS		ElectroCap XS, child, size: 46-50 cm, colour: green, 20 electrodes For current products: EEG-1200/9100 Packing unit: 1 pcs 
----------	--	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
NKD-E3S		Chest Band for all ElectroCaps, S, adult, colour: yellow, band length: 72 cm For current products: EEG-1200/9100 Packing unit: 1 pcs 



NKD-E3M		Chest Band for all ElectroCaps, M, adult, colour: red, band length: 95 cm For current products: EEG-1200/9100 Packing unit: 1 pcs 
---------	--	---

NKD-E23/5		Chin Strap for all ElectroCaps, child, colour: white, strap length: 12.7 cm For current products: EEG-1200/9100 Packing unit: 1 pair 
-----------	--	--



NKD-E23/7		Chin Strap for all ElectroCaps, adults, colour: white, strap length: 17.78 cm For current products: EEG-1200/9100 Packing unit: 1 pair 
-----------	--	--

12.2.4.

ElectroCaps for Infants

Supply Code **Model** **Description**
NKD-11/1F ElectroCap, neonate/infant, size: 42-46 cm, colour: light blue, 14 flat Sn electrodes
For current products: EEG-1200/9100
Packing unit: 1 pcs 



NKD-11/2F ElectroCap, neonate/infant, size: 38-42 cm, colour: pink, 14 flat Sn electrodes
For current products: EEG-1200/9100
Packing unit: 1 pcs 

NKD-11/3F ElectroCap, neonate/infant, size: 34-38 cm, colour: brown, 14 flat Sn electrodes
For current products: EEG-1200/9100
Packing unit: 1 pcs 

NKD-11/4F ElectroCap, neonate/infant, size: 30-34 cm, colour: yellow, 14 flat Sn electrodes
For current products: EEG-1200/9100
Packing unit: 1 pcs 

NKD-11/5F ElectroCap, neonate/infant, size: 26-30 cm, colour: blue, 14 flat Sn electrodes
For current products: EEG-1200/9100
Packing unit: 1 pcs 

NKD-11/6F ElectroCap, neonate/infant, size: 22-26 cm, colour: red, 14 flat Sn electrodes
For current products: EEG-1200/9100
Packing unit: 1 pcs 

NKD-12 Chest Band for NKD-11/xF, neonate/infant, colour: white, band length: 40 cm
For current products: EEG-1200/9100
Packing unit: 1 pcs 



NKD-123 Chin Strap for all NKD-11/xF, infant/child, colour: white, strap length: 9 cm
For current products: EEG-1200/9100
Packing unit: 1 pair 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

12.2.5.

ElectroCap Accessories

Facts

The coloured sections on the tape correlate with the colours of the different ElectroCaps and indicate therefore which size will be the appropriated one for the patient.

Supply Code	Model	Description
NKD-E12		Head Measuring Tape for Adult ElectroCaps, colour: white For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-B3		Head Measuring Tape for Infant ElectroCaps, colour: white For current products: EEG-1200/9100 Packing unit: 1 pcs 



Facts

The chest bands can be easily adapted to the patients' body size thanks to their elastic material and the Velcro fastener. Those band help to secure the cap placement and minimize eventual movement artefacts or high impedances.

Supply Code	Model	Description
NKD-I2		Chest Band for NKD-I1/xF, neonate/infant, colour: white, band length: 40 cm For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-E3S		Chest Band for all ElectroCaps S, adult, colour: yellow, band length: 72 cm For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-E3M		Chest Band for all ElectroCaps M, adult, colour: red, band length: 95 cm For current products: EEG-1200/9100 Packing unit: 1 pcs 
NKD-E3L		Chest Band for all ElectroCaps L, adult, colour: blue, band length: 1.16 m For current products: EEG-1200/9100 Packing unit: 1 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
NKD-I23		Chin Strap for all NKD-I1/xF, infant/child, colour: white, strap length: 9 cm For current products: EEG-1200/9100 Packing unit: 1 pair R



NKD-E23/5		Chin Strap for all ElectroCaps, child, colour: white, strap length: 12.7 cm For current products: EEG-1200/9100 Packing unit: 1 pair R
-----------	--	---



NKD-E23/7		Chin Strap for all ElectroCaps, adult, colour: white, strap length: 17.78 cm For current products: EEG-1200/9100 Packing unit: 1 pair R
-----------	--	--

Facts

The cap bands are applied on the ElectroCaps and chest bands to ensure that the caps do not move. Thanks to their elastic material, a flexible attachment is possible.

Supply Code	Model	Description
NKD-E13S		Cap Band for all ElectroCaps S, colour: white, length: 15-18 cm For current products: EEG-1200/9100 Packing unit: 1 pcs R



NKD-E13M		Cap Band for all ElectroCaps M, colour: white, length: 16-20 cm For current products: EEG-1200/9100 Packing unit: 1 pcs R
----------	--	--

NKD-E13L		Cap Band for all ElectroCap L, colour: white, length: 16.5-21 cm For current products: EEG-1200/9100 Packing unit: 1 pcs R
----------	--	---

Facts

This electrode is an intermediate electrode to overcome sudden signal problems of a defective electrode during an examination. It can be attached to the frame of the defective electrode and directly connected to the EEG system until the cap will be completely checked and repaired.

Supply Code	Model	Description
NKD-E4		Quick Insert Electrode, 1.5 mm touch proof connector, cable length: 1.2 m, for all ElectroCaps, Note: Please indicate if you need a 2 mm connector. For current products: EEG-1200/9100 Packing unit: 1 pcs D



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!


Facts

Disposable foam rings can be used to absorb sweat, to avoid the expansion of the conductive gel and to prevent the ElectroCaps from moving. Thanks to the adhesive backing they can be easily attached to the forehead or on any other hairless area. Their use is especially recommended for infants because of the more comfortable soft feeling of the caps on the head.

Supply Code	Model	Description	
NKD-E6		Disposable Sponge Disks for ElectroCaps, for fixing on forehead, adhesive on one side For current products: EEG-1200/9100 Packing unit: 100 pcs 	
NKD-E5/6Z-26		Ear Electrodes for Infant ElectroCaps, Sn, 1.5 mm touch proof connector, electrode size: 6 mm dia, cable length: 1.22 m For current products: EEG-1200/9100 Packing unit: 1 pair 	
NKD-E5/9Z-26		Ear Electrodes for Child/Adult ElectroCaps, Sn, 1.5 mm touch proof connector, electrode size: 9 mm dia, cable length: 1.22 m For current products: EEG-1200/9100 Packing unit: 1 pair 	
NKD-E7		Disposable Needles & Syringe Kit for ElectroCaps, incl. 2 x blunted needles & 1 x syringe For current products: EEG-1200/9100 Packing unit: 1 kit 	
NKD-E7B		Disposable Syringes for ElectroCaps For current products: EEG-1200/9100 Packing unit: 100 pcs 	
NKD-E8		Disposable Blunted Needle for NKD-E7, for all ElectroCaps, needle length: 20 mm For current products: EEG-1200/9100 Packing unit: 1 pcs 	
NKD-E8B		Disposable Blunted Needle for NKD-E7, for all ElectroCaps, needle length: 20 mm For current products: EEG-1200/9100 Packing unit: 100 pcs 	

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description	
NKD-E16/700		ElectroCap Cleaning Emulsion Ivory, 700 ml bottle, for all ElectroCaps For current products: EEG-1200/9100 Packing unit: 1 pcs 	
NKD-E9		Electrode Gel for all ElectroCaps, 480 g cup For current products: EEG-1200/9100 Packing unit: 1 pcs 	
NKD-E10		Electrode Gel for all ElectroCaps, 960 g cup For current products: EEG-1200/9100 Packing unit: 1 pcs 	
NKD-E11		Conductive Electrolyte Gel for all ElectroCaps, 3780 g cup For current products: EEG-1200/9100 Packing unit: 1 pcs 	
NKD-E25		Double-Sided Adhesive Electrode Pads for Disk Electrodes, for all ElectroCaps For current products: EEG-1200/9100 Packing unit: 100 pcs 	
NKD-E2		Electrode Board Adapter for ElectroCaps XS-L, 1.5 mm touch proof connector Note: Please indicate if you need a 2 mm connector. For current products: EEG-1200/9100 Packing unit: 1 pcs 	
NF-E2/NK-9G		ElectroCaps Adapter Cable for JE-910AG/921AG For current products: EEG-1200/9100 Packing unit: 1 pcs 	
NF-E2/NK-9K		ElectroCaps Adapter Cable for JE-910A/921A For current products: EEG-1200/9100 Packing unit: 1 pcs 	

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

12.3. BraiNet Caps

BraiNet caps are perfectly suitable for quick and fast electrodes positioning for EEG recordings in intensive care or emergency units. The latex-free elastic material avoids allergic reactions of the patients. As the caps should be only used for a single patient, no cleaning is required and cross-infections are reduced. The use of BraiNet makes a classic measuring or marking of the head redundant. They can be easily adapted to the patients' scalp with adjustable straps. The colour-coded areas correspond to the colours of the jack box overlays that can be put on the jack box inputs for simply recognition which electrode needs to be attached to which input. Thanks to the easy application a faster diagnoses is possible as the preparation time is reduced to 10 minutes.

The large portfolio of sizes ensures a proper fit for all patients:

Original BraiNet caps are intended for paediatric and adult patients in emergency and intensive care units with a reduced number of electrodes. The 10/20 BraiNet caps are designed for routine EEG, epilepsy monitoring, ambulatory EEG and continuous EEG. Baby BraiNet caps are designed for STAT and amplitude-integrated EEG in the neonatal intensive care units.

Supply Code	Model	Description
NKD-A205		Adhesive ECG Rings, size: 20 x 5 mm For current products: all ECG; EEG; BSM-3000/6000/9100 Packing unit: 500 pcs 
NKD-BN100		Original BraiNet S, child/adult, size: 46-54 cm circumference, 14 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 
NKD-BN101		Original BraiNet L, child/adult, size: 54-62 cm circumference, 14 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 
NKD-BN106		10-20 BraiNet, child/adult, size: 45.5-49 cm circumference, 21 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 
NKD-BN107		10-20 BraiNet, child/adult, size: 49.5-54 cm circumference, 21 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 
NKD-BN108		10-20 BraiNet, child/adult, size: 54.5-57 cm circumference, 21 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
NKD-BN109		10-20 BraiNet, child/adult, size: 57.5-62 cm circumference, 21 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 

NKD-BN111		Baby BraiNet 1, NICU, size: 25-29 cm circumference, 14 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 
-----------	--	---



NKD-BN112		Baby BraiNet 2, NICU, size: 29.5-33 cm circumference, 14 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 
-----------	--	---

NKD-BN113		Baby BraiNet 3, NICU, size: 33.5-37 cm circumference, 14 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 
-----------	--	---

NKD-BN114		Baby BraiNet 4, NICU, size: 37.5-41 cm circumference, 14 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 
-----------	--	---

NKD-BN115		Baby BraiNet 5, NICU, size: 41.5-45 cm circumference, 14 electrodes For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 pcs 
-----------	--	---

NKD-BN500		Baby BraiNet Head Box Overlay, adhesive coloured stickers for jack box inputs customization For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 kit 
-----------	--	---



NKD-BN505		10/20 BraiNet Head Box Overlay, adhesive coloured stickers for jack box inputs customization For current products: all EEG; BSM-3000/6000/9100 Packing unit: 1 kit 
-----------	--	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

13. Additional EEG Accessories

13.1. Pastes & Creams

13.2. Bags

13.3. Recording Paper

13.4. Recording Ink

13.5. Recording Pens

13.1. Pastes & Creams

Supply Code	Model	Description
F010	Z-101BC	cardioCream ECG Electrode Paste, 100 g tube Packing unit: 2 pcs 

cardioCream

F020	YZ-0019	skinPure Skin Preparation Gel, reduces the electrode skin impedance, 135 g tube Note: Can be used with all surface applied electrodes. Packing unit: 2 pcs 
------	---------	--

skinPure**Facts**

The NIHON KOHDEN electrode paste Elefix integrates outstanding sticking qualities with an excellent electrical conductivity. Additionally Elefix has a creamy quality so that the electrodes can be attached firmly to the patient without additional fixation, on the other hand Elefix does not harden so that the electrodes can be removed without any problems even after longer examinations.

The contact impedance between skin and electrode is greatly reduced. In addition, the signal transmission remains stable even under unfavourable external conditions.

NIHON KOHDEN Elefix can be used with all surface applied electrodes.

Supply Code	Model	Description
F509	Z-181JE	Elefix EEG Conduction Paste, 180 g tube Packing unit: 1 pcs 

Elefix

F510	Z-401CE	Elefix EEG Conduction Paste, 400 g cup Packing unit: 1 pcs 
------	---------	---

Elefix

Supply Code	Model	Description
NKD-NUPREP		Skin Preparation Gel for Neonates, 114 g tube, reduces the skin impedance Packing unit: 3 pcs 



NKD-EC2		Conductive Grass Paste EC2 for long-term EEG measurement, 100 g tube Packing unit: 10 pcs 
---------	--	--



KOL-1		Collodion Kit, 6.5 ml bottle Packing unit: 20 pcs 
-------	--	--



13.2. Bags

Supply Code	Model	Description	
NKD-POUCH-JE215		Bag for Mini Input Box JE-215/216/217AK For current products: EEG-1200 Packing unit: 1 pcs R	
NKD-POUCH-JE225		Bag for Mini Input Box JE-225/226/227AK For current products: EEG-1200 Packing unit: 1 pcs R	
NKD-POUCH-JE922		Bag for Mini Input Box JE-113/114/914/922A, incl. belt and shoulder strap For current products: EEG-1200/9100 Packing unit: 1 pcs R	
Y222		Carrying Bag for Air EEG, colour: purple & black, width: 13 cm, length: 16 cm, height: 3 cm For current products: WEE-1000 Packing unit: 1 pcs R	
Y763A		Shoulder Strap for Air EEG, width: 3.8 cm, length: 1.5 m For current products: WEE-1000 Packing unit: 1 pcs R	
Y763B		Belt for Air EEG, width: 3.8 cm, length: 0.5 m For current products: WEE-1000 Packing unit: 1 pcs R	

D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

13.3. Recording Paper

Supply Code	Model	Description
NKD-155/15		Thermal Paper, Z-fold, 155 mm x 150 m (page 150 mm), 1000 pages, with grid/queue mark, without perforation For old products: EEG-9109/7209 Packing unit: 1 pcs 
NKD-210T		Thermal Paper, Z-fold, 210 mm x 300 m (page 300 mm), 1000 pages, vertical orange lines, with queue mark/perforation For old products: EEG-5600 Packing unit: 1 pcs 
NKD-245/1	FV245-30	Recording Paper, Z-fold, 245 mm x 300 m (page 300 mm), 1000 pages, vertical orange lines, with queue mark/perforation For old products: EEG-5200/6510/7300/8300 Packing unit: 1 pcs 
A720B	FQW210-10B	Thermal Paper, Z-fold, 210 mm x 100 m (page 300 mm), 333 pages, without grid, with queue mark/perforation For old products: EEG-5600; MEB-4200; MEM-4200/4104 Packing unit: 5 pcs 
C451	FV245-10	Recording Paper, Z-fold, 245 mm x 100 m (page 300 mm), 333 pages, vertical orange lines, with queue mark/perforation For old products: EEG-4409/5109/5113/5208/5210/5214/5400/6109/6113/7213/7310/7314/7400 Packing unit: 1 pcs 
C452	FV245-30	Recording Paper, Z-fold, 245 mm x 200 m (page 300 mm), 333 pages, vertical orange lines, with queue mark/perforation For old products: EEG-4409/5109/5113/5208/5210/5214/5400/6109/6113/7213/7310/7314/7400 Packing unit: 1 pcs 
C474	FV345-30	Recording Paper, Z-fold, 345 mm x 300 m (page 300 mm), 1000 pages, vertical orange lines, with queue mark/perforation For old products: EEG-1218/1514/1518/3114/3118/4414/4418/4514/4518/6718 Packing unit: 1 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
C491	FV400-30	Recording Paper, Z-fold, 400 mm x 300 m (page 300 mm), 1000 pages, vertical orange lines, with queue mark/perforation For old products: EEG-1224/1524/3121/4412/4421/4524 Packing unit: 1 pcs 



C810	RIW114-3	Recording Paper, roll, 114 mm x 30 m, without grid, without queue mark/perforation For old products: AEPH Packing unit: 10 pcs 
------	----------	--

13.4. Recording Ink

Supply Code	Model	Description
F515	NSE-40BL	Black Ink, 400 ml bottle, with flexible pouring nozzle For old products: EEG-4400 Packing unit: 1 pcs 



F516	NSE-25BL	Black Ink, 250 ml cartridge For old products: EEG-7300/7400/8300 Packing unit: 2 pcs 
------	----------	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

13.5. Recording Pens

Supply Code	Model	Description	
M520		Recording Pen SJ-85K, ruby tip, with holder, length: 8.5 cm For old products: EEG-7100 Packing unit: 2 pcs R	
M521	YE-006A	Recording Pen SJ-85M, carbide tip, with holder, length: 8.5 cm For old products: EEG-7209/7213/8213 Packing unit: 2 pcs R	
M523		Recording Pen SJ-85S, ruby tip, with holder, length: 8.5 cm For old products: EEG-7200 Packing unit: 2 pcs R	
M524A	YE-010A	Recording Pen SJ-120AM, steel tip, with holder, length: 12 cm For old products: EEG-1214/1218/1224/1500/1714/3100/ 4400/5400/6700/7300/7400 Packing unit: 2 pcs R	
M528		Recording Pen SJ-120N, ruby tip, with holder, length: 12 cm For old products: EEG-4100/7300/7400/8300 Packing unit: 2 pcs R	
M531A	YE-017A	Recording Pen SJ-120SM, ruby tip, with holder, length: 12 cm For old products: EEG-4400/7300/7400/8300 Packing unit: 2 pcs R	
M540	6114-009395	Ink tube, for all EEG systems with ink recording system, length: 2.0 m Packing unit: 1 pcs R	
NKD-INBUS		INBUS Wrench, for mounting of the pen and for adjustment of the pen pressure, size: 1.5 mm Packing unit: 1 pcs R	

D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

III. Accessories & Consumables for amplitude-integrated Electroencephalography (aEEG)


14. aEEG Measurement

Amplitude-integrated EEG (aEEG) is a Cerebral Function Monitoring (CFM) especially used in neonatal intensive care units (NICU) giving immediately an important information about the brain activity. The 2-channel aEEG employs the electrode positions C3, P3, C4 and P4. The reference electrode is usually applied on a neutral location such as central forehead or the shoulder. The method is based on filtered and time compressed EEG reflecting the minimum and maximum amplitudes of “raw EEGs”. This enables an easy evaluation of changes and trends in electrocortical activity. Neurological brain disorders are identified with relatively simple pattern recognition between normal and abnormal brain activity such as seizures, sleep-wake cycles, low voltage, burst suppression or flat trace. Cerebral function disturbances are therefore recognized at a very early stage and allow an adequate therapy.

14.1. aEEG Starter Kits

14.2. aEEG Accessories

The aEEG Booklet



A practical quick overview

NIHON KOHDEN and Karl Florian Schettler presented a reference book about the amplitude-integrated EEG.

It fits perfectly in every doctor's pocket: "The aEEG booklet", a practical reference book in a manageable paper book format, that demonstrates the range of application and the possibilities of the amplitude-integrated EEG (aEEG) in a quick and concise manner. The book is with its unique handy size a perfect fill-in between practice and theory.

Though "The aEEG booklet" is quite more than a quick help for medics and clinical expert staff, Karl Florian Schettler also demonstrates the multi-purposed options of the appliance in the field of neonatology.

With the help of case studies, Schettler shows how valuable the method is for the monitoring of brain function of newborn babies. "My intention is to bring forward the aEEG method and to establish a neurological early-warning system", says Schettler: "For instance with the aEEG physicians can early identify neurologic problems of premature infants".

Physicians can order "The aEEG booklet" free of charge in English version at NIHON KOHDEN EUROPE (www.nihonkohden.net) or by contacting our sales representatives.

14.1. aEEG Starter Kits

Facts

“Simplify your electrode application”

In practice different aEEG application possibilities exist that are mainly dependent on the experience, preference of the user or costs: subdermal needle electrodes, cup electrodes and hydrogel electrodes. Low impedances and good signal quality combined with an easy application and long-lasting fixation are expectations to those electrodes.

NIHON KOHDEN created three dedicated aEEG Starter Kits to respond to the different preferences of the users and to simplify your electrode application. They are provided in durable boxes which contain all relevant accessories to do the aEEG recording with the chosen electrode type. The re-closable box helps you to store all aEEG material at one place, ready for use.

Supply Code	Model	Description
NKD-AEEG/K3		aEEG Kit 3 – Subdermal Needle Electrodes Kit Kit contains: 1 package of subdermal needle electrodes (NKD-NM2512), 2 aEEG positioning rulers (NF-EEG/HELP), 1 coloured head box overlay (NKD-BN505) For current products: EEG-1200/9100; BSM-6000/9100 Packing unit: 1 kit



NKD-NM2512



NF-EEG/HELP



aEEG Application Steps for Subdermal Needle Electrodes

Steps

1 Electrode positioning

Description

- Define the electrode positions according to 10/20 system by using the positioning helps (i.e. NF-EEG/Help).



2 Initial cleaning

- Clean the electrode placement areas from blood, vernix, meconium with a gauze and a skin cleanser as antiseptic. Attention: The skin of preterms is very thin. Be careful not to injure the skin.
- Please clean as for transcutaneous punctures.



3 Insertion

- Prestress the according skin area in order to reduce insertion pain for the patient.
- Do not insert needles towards each other to avoid that they will come into contact with each other under the skin and to avoid a “short circuit” = artefact.
- Quickly insert the needles with an angle of around 45° to minimize the pain for the patient.
- Insert each needle electrode subcutaneously up to the plastic hub. Avoid further insertion which could cause muscle artefacts:
 - Place the BLUE electrode (LB) on the left back part of the scalp -> P3
 - Place the VIOLET electrode (LF) on the left frontal/central part of the scalp -> C3
 - Place the GREY electrode (RF) on the right frontal/central part of the scalp -> C4
 - Place the WHITE electrode (RB) on the right back part of the scalp -> P4
 - Place the BLACK electrode (Z) as reference electrode -> Z-electrode



4 Fixation

- Once applied, fix each needle electrode using an adhesive tape.
- Plug into the data acquisition unit.



! Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

! Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code
NKD-AEEG/K1

Model

Description

aEEG Kit 1 – Hydrogel Electrodes Kit

Kit contains: 2 sizes of hydrogel electrodes (NKD-N300, NKD-N306), 3 tubes of skin preparation gel (NKD-NUPREP), 5 electrode extension cables (NKD-KR26DIN15), 2 aEEG positioning rulers (NF-EEG/HELP), 1 coloured head box overlay (NKD-BN505)

For current products: EEG-1200/9100; BSM-6000/9100

Packing unit: 1 kit



NKD-N300




NKD-N306



NKD-NUPREP



NKD-KR26DIN15



NF-EEG/HELP



NKD-BN505

Supply Code
NKD-AEEG/K2

Model

Description

aEEG Kit 2 – Cup Electrodes Kit

Kit contains: 1 package of cup electrodes (H526), 3 tubes skin preparation gel (NKD-NUPREP), 1 tube of conductive EEG paste (NKD-EC2), 2 aEEG positioning rulers (NF-EEG/HELP), 1 coloured head box overlay (NKD-BN505)

For current products: EEG-1200/9100; BSM-6000/9100

Packing unit: 1 kit




H526



NKD-NUPREP



NKD-EC2



NF-EEG/HELP



NKD-BN505

🗑 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

♻ Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

aEEG Application Steps for Hydrogel and Cup Electrodes



Steps

1 Electrode positioning

Description

- Define the electrode positions according to 10/20 system by using the positioning helps (i.e. NF-EEG/Help).



2 Initial cleaning

- Clean the scalp from blood, vernix, meconium.



3 Free the areas from hairs

- Prepare the electrode positions by exposing the scalp. Part the child's hair with e.g. using a toothbrush, a fine comb or wet cotton stick. If necessary remove hairs locally, e.g. using a nasal hair trimmer.



4 Skin preparation

- Remove remaining contaminants from scalp (see point 2). A cotton stick helps keeping your already prepared skin area.
- Apply your skin peeling cream (e.g. Nuprep).
- Spread the peeling cream continuously and strongly in the skin area.
- Do not use peeling cream with small preterm infants (23-28 week of gestational age).




5 Hydrogel electrodes fixation

- When fixing the electrodes, ensure that the cables always run upwards.
- Warm-up hydrogel electrodes before attachment for good functioning.



6 Cup electrodes fixation

- Make sure that cup electrodes will be not attached over hairs.
- Use an adhesive and conductive fixation gel (e.g. Grass paste).



14.2. aEEG Accessories

Facts

Subdermal needle electrodes for aEEG application can be quickly applied and show immediately a good signal quality. Once inserted they usually do not drop off or get out of place and are perfectly for long-term aEEG monitoring or for asphyxia detection. The additional advantage of NIHON KOHDEN aEEG needles is the user-friendly design. The unique flat cable philosophy over the bigger cable part solves the problem of entangled cables. Tensions on the needles can be easily adjusted by detaching the relevant cable from the other ones as the needle ends are not attached together. The placement of the reference needle electrode on the side allows the user to employ also another type of electrode as reference when simply detaching the needle electrode from the rest of the flat cable. The length of 1.5 m of the flat cable makes an additional connection cable for the EEG junction box inputs redundant. The standard 1.5 mm touch proof connectors can be easily attached to the EEG inputs of any EEG system. A simplified application is moreover assured by the colour-coded cables and the positioning label on the scalp at both cable ends:

Needle Electrode 1: Z=Ground (colour: black)

Needle Electrode 2: RB=Right Back=P3 (colour: white)

Needle Electrode 3: RF=Right Front=C4 (colour: grey)

Needle Electrode 4: LF=Left Front=C3 (colour: violet)

Needle Electrode 5: LB=Left Back=P3 (colour: blue)

Supplementary safety is assured by sterile pouches and five single “exchange needles” in addition to the 10 flat cables per box. Those replacement needles can be used to overcome single needles that become unsterile by accident. This solution avoids to replace the complete flat cable with five needles when only one needle needs to be replaced and is therefore a cost-saving solution. NIHON KOHDEN aEEG needles are easy in handling and genius in design.

Supply Code	Model	Description
NKD-NM2512		aEEG Needles for Preterms/Neonates with flat cable tree, size: 0.25 x 12 mm (31 G) For current products: EEG-1200/9100; BSM-6000/9100 Packing unit: 10 x 5 pcs with flat cable & 5 single exchange needles 



Facts

Hydrogel electrodes are characterized by an easy application as they are self-adhesive. To obtain low impedances, the skin preparation is very important. They are particularly used for immature preterm babies with little hair.

Supply Code	Model	Description
NKD-N300		Mini Hydrogel Electrodes, preterms > 1.5 kg, size: 25 mm dia For current products: EEG-1200/9100; BSM-6000/9100 Packing unit: 20 x 3 pcs 



NKD-N306		Micro Hydrogel Electrodes, preterms < 1.5 kg, size: 10 x 25 mm For current products: EEG-1200/9100; BSM-6000/9100 Packing unit: 20 x 3 pcs 
----------	--	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Facts

Cup electrodes are well known in conventional neuro-paediatric EEG monitoring and can be easily attached to hair-covered areas by conductive paste. As for hydrogel electrodes low impedances depend on the skin preparation before the recording.

Supply Code	Model	Description
H526	NE-134A	EEG Electrode (Collodion), Ag, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100; BSM-3000/6000/9100; PSG-1100; MEB-2300; MEE-1000 Packing unit: 12 pcs 



Application Steps for aEEG Positioning Help



Steps

Description

1 Size

- Selection of the right size of the application help.
- 2 Sizes are available at NIHON KOHDEN: SMALL and NORMAL.


2 Ear Tragus

- Identify the TRAGUS point of the ear.



3 Sagittal Suture

- Imagine the SAGITAL SUTURE of the head.



4 Start

- Take the white side of the application help.
- Apply the upper part (marked with EAR TRAGUS) to the Ear Tragus.
- Put the application help totally over the ear so that you can see the Ear Tragus.
- Apply the application help vertically over the scalp.
- Make sure that the letter at the EAR TRAGUS is the same as at the SAGITAL SUTURE.
- Once the letters are the same, mark at both sides of the flash the electrode position. The flash guarantees the same distance on both sides of the head.



5 Repeat

- Repeat the procedure from Point 4 from the other side of the head to identify the 2 other electrode positions.



6 Reference

- The reference electrode can be applied either over the shoulder blade or centrally to the front.
- Use an adhesive and conductive fixation gel (e.g. Grass paste).



Supply Code	Model	Description
NF-EEG/HELP		Positioning Ruler Kit for aEEG For current products: EEEG-1200/9100; BSM-6000/9100 Packing unit: 1 kit 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
NKD-KR26DIN15		<p>Extension Cable, 1.5 mm DIN jack to 1.5 mm touch proof connector, cable length: 1.5 m</p> <p>For current products: EEG-1200/9100; BSM-6000/9100</p> <p>Packing unit: 5 pcs R</p>
NKD-EC2		<p>Conductive Grass Paste EC2 for long-term EEG measurement, 100 g tube</p> <p>For current products: EEG-1200/9100; BSM-3000/6000/9100</p> <p>Packing unit: 10 pcs D</p>
NKD-NUPREP		<p>Skin Preparation Gel for Neonates, 114 g tube, reduces the skin impedance</p> <p>For current products: EEG-1200/9100; BSM-6000/9100</p> <p>Packing unit: 3 pcs D</p>
NKD-BN505		<p>10/20 BraiNet Head Box Overlay, adhesive coloured stickers for head box port customization</p> <p>For current products: EEG-1200/9100; BSM-3000/6000/9100</p> <p>Packing unit: 1 kit D</p>



Facts

In order to simplify the aEEG application a special Minibox for the EEG junction box JE-921A was created with only five inputs which are mandatory for aEEG. Therefore there is less confusion which inputs should be used as you only need to plug in your electrodes in the aEEG Minibox that is connected to the head box by a board adapter.

The same principle as for the subdermal aEEG needle electrodes NKD-NM2512 included in the aEEG Starter Kit 3 was kept: easy and quick indication of the position by genius labels. Beside the common neurology electrode position names the inputs are marked with easy understandable abbreviations:

Electrode 1: Z=Ground (colour: black)

Electrode 2: RB=Right Back=P3 (colour: white)

Electrode 3: RF=Right Front=C4 (colour: grey)

Electrode 4: LF=Left Front=C3 (colour: violet)

Electrode 5: LB=Left Back=P3 (colour: blue)

Supply Code	Model	Description
NKD-aEEG/ Minibox		<p>Mini Input Box for EEG Head Box JE-921A, with 5 inputs for 2-channel aEEG (C3, C4, P3, P4, Z), cable length: 4.0 m</p> <p>For current products: EEG-1200/9100</p> <p>Packing unit: 1 pcs R</p>




D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

IV. Accessories & Consumables for Polysomnography (PSG)

Polysomnography (PSG) is a multi-parametric examination of the patients' biophysiological condition during an overnight study including brain activity (EEG), eye movements (EOG), muscle activity (EMG), electrocardiogram (ECG), airflow, pulse oximetry (SpO₂) and expired carbone dioxide (ETCO₂). The overall monitoring of the patient during sleep helps to detect sleep disorders such as narcolepsy, sleep apnoea, insomnias or restless legs. Based on the recorded information sleep specialists evaluate the results and define patients' therapies.

15. Sleep Study Accessory Kits

The background of the page is a light blue gradient with several overlapping, semi-transparent, wavy bands of a slightly darker blue. These bands create a sense of movement and depth, flowing from the top left towards the bottom right. The overall effect is clean, modern, and professional.

Supply Code Model Description


NKD-SLEEPKITA

Standard Starter Kit for Sleep Diagnostics Adult for Head Box JE-912 and PSG-1100!

Kit contains: 1 thermocouple sensor adult (NKD-SLP1401), 1 kit of inductive belts with 1.2 m (NKD-SLP9003-L120), 1 kit of chest and abdomen interface cables (NKD-INDUCTKIT), 1 snore sensor (NKD-SLP1250), 1 body position sensor 3.5 mm (NKD-SLP1575), 1 box of disposable ECG electrodes adult (G203) and according snap connection cables (NKD-15026N), 1 package of EEG cup electrodes (H526), 1 tube of conduction paste (NKD-EC2), 4 extension cables (NKD-KR26DIN20), 2 tubes of skinPure skin preparation gel (F020)

For current products: EEG-1200/9100; PSG-1100

Packing unit: 1 kit



NKD-SLP1401



NKD-SLP9003-L120



NKD-SLP9101



NKD-SLP9102



NKD-SLP1250



NKD-SLP1575



G203



NKD-15026N



H526



NKD-EC2



NKD-KR26DIN20



F020

🗑 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

🔄 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code Model Description

NKD-SLEEPKITP

Standard Starter Kit for Sleep Diagnostics Paediatric for Head Box JE-912 and PSG-1100!

Kit contains: 1 thermocouple sensor child (NKD-SLP1401P), 1 kit of inductive belts with 0.6 m (NKD-SLP9003-L60), 1 kit of chest and abdomen interface cables (NKD-INDUCTKIT), 1 snore sensor (NKD-SLP1250), 1 body position sensor 3.5 mm (NKD-SLP1575), 1 box of disposable ECG electrodes child (G210C) and according snap connection cables (NKD-15026N), 1 package of EEG cup electrodes (H526), 1 tube of conduction paste (NKD-EC2), 4 extension cables (NKD-KR26DIN20), 2 tubes of skinPure skin preparation gel (F020)

For current products: EEG-1200/9100; PSG-1100

Packing unit: 1 kit



NKD-SLP1401P



NKD-SLP9003-L60



NKD-SLP9101



NKD-SLP9102



NKD-SLP1250



NKD-SLP1575



G210C



NKD-15026N



H526



NKD-EC2



NKD-KR26DIN20



F020

ⓘ Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

Ⓜ Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code Model Description

NKD-SLEEPKITA921

Standard Starter Kit for Sleep Diagnostics Adult for Head Box JE-921!

Kit contains: 1 thermocouple sensor adult (NKD-SLP1401), 1 kit of inductive belts with 1.2 m (NKD-SLP9003-L120), 1 kit of chest and abdomen interface cables (NKD-INDUCTKIT), 1 snore sensor (NKD-SLP1250), 1 body position sensor 2.5 mm (NKD-SLP1582), 1 box of disposable ECG electrodes adult (G203) and according snap connection cables (NKD-15026N), 1 package of EEG cup electrodes (H526), 1 tube of conduction paste (NKD-EC2), 4 extension cables (NKD-KR26DIN20), 2 tubes of skinPure skin preparation gel (F020)

For current products: EEG-1200/9100

Packing unit: 1 kit



NKD-SLP1401



NKD-SLP9003-L120



NKD-SLP9101



NKD-SLP9102



NKD-SLP1250



NKD-SLP1582



G203



NKD-15026N



H526



NKD-EC2



NKD-KR26DIN20



F020

🗑 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

🔄 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code Model Description

NKD-SLEEKITP921

Standard Starter Kit for Sleep Diagnostics Paediatric for Head Box JE-921!

Kit contains: 1 thermocouple sensor child (NKD-SLP1401P), 1 kit of inductive belts with 0.6 m (NKD-SLP9003-L60), 1 kit of chest and abdomen interface cables (NKD-INDUCTKIT), 1 snore sensor (NKD-SLP1250), 1 body position sensor 2.5 mm (NKD-SLP1582), 1 box of disposable ECG electrodes child (G210C) and according snap connection cables (NKD-15026N), 1 package of EEG cup electrodes (H526), 1 tube of conduction paste (NKD-EC2), 4 extension cables (NKD-KR26DIN20), 2 tubes of skinPure skin preparation gel (F020)

For current products: EEG-1200/9100

Packing unit: 1 kit



NKD-SLP1401P



NKD-SLP9003-L60



NKD-SLP9101



NKD-SLP9102



NKD-SLP1250



NKD-SLP1582



G210C



NKD-15026N



H526



NKD-EC2



NKD-KR26DIN20



F020

ⓘ Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

ⓘ Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code
NKD-CPAPKITA

Model

Description

CPAP Starter Kit for Sleep Diagnostics Adult

Kit contains: 1 DC pressure sensor kit 2.5 mm (NKD-SLP14834KIT),
1 pressure transducer (NKD-PTRANS)

For current products: EEG-1200/9100

Packing unit: 1 kit



NKD-SLP14834KIT



NKD-PTRANS

NF-SCREEN/AC2

Accessory Kit for Trackit Sleepwalker

Kit contains: 1 kit of 5 adult cannulas (NKD-SLP15805AFT-10), 1 kit of inductive belts with 1.2 m (NKD-SLP9003-L120), 2 interface cables abdomen and chest with keyhole connector (NKD-SLP9101NO, NKD-SLP9102NO), 1 kit of limb movement sensor kit with keyhole connector (NKD-SLP1771KIT), 2 Y-Splitter (NKD-SLP1726), 2 extension cables for interface cables (NKD-SLP1326)

For current products: NF-SCREENER

Packing unit: 1 kit




NKD-SLP15805AFT-10



NKD-SLP9003-L120



NKD-SLP9101NO



NKD-SLP9102NO



NKD-SLP1771KIT (NKD-SLP1771, NKD-SLP1730)



NKD-SLP1726

ⓘ Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

Ⓜ Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16. Sleep Study Accessories

- 16.1. EEG Electrodes
- 16.2. Pastes & Creams
- 16.3. ECG Electrodes
- 16.4. Snore Sensors
- 16.5. Body Position Sensors
- 16.6. Respiratory Effort Sensors
- 16.7. Flow Sensors
- 16.8. Limb Movement Sensors
- 16.9. NIHON KOHDEN SpO₂ Probes
- 16.10. NONIN SpO₂ Probes
- 16.11. ETCO₂ Measurement
- 16.12. Sensor Tester

16.1. EEG Electrodes

Supply Code	Model	Description
H526	NE-134A	EEG Cup Electrode (Collodion), Ag, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100; BSM-3000/6000/9100; PSG-1100; MEB-2300; MEE-1000 Packing unit: 12 pcs R



H527	NE-136A	EEG Cup Electrode (Collodion), Ag, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 0.7 m For current products: EEG-1200/9100; BSM-3000/6000/9100; PSG-1100; MEB-2300; MEE-1000 Packing unit: 12 pcs R
------	---------	---



NKD-FE5GH		Grass Gold Disk Electrode, Au, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 10 pcs R
-----------	--	---



NKD-FE5GH-07		Grass Gold Disk Electrode, Au, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 0.7 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 10 pcs R
--------------	--	---

NKD-FE5GH-12		Grass Gold Disk Electrode, Au, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 1.2 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 10 pcs R
--------------	--	---

NKD-FE5GH-24		Grass Gold Disk Electrode, Au, size: 10 mm dia, 1.5 mm touch proof connector, cable length: 2.4 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 10 pcs R
--------------	--	---

NKD-FE6GH60		Grass Gold Disk Electrode, Au, size: 6 mm dia, 1.5 mm touch proof connector, cable length: 1.5 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 10 pcs R
-------------	--	--

D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.2. Pastes & Creams

Supply Code	Model	Description
F010	Z-101BC	cardioCream ECG Electrode Paste, 100 g tube Packing unit: 2 pcs 

cardioCream



F020	YZ-0019	skinPure Skin Preparation Gel, reduces the electrode skin impedance, 135 g tube Note: Can be used with all surface applied electrodes. Packing unit: 2 pcs 
------	---------	--

skinPure

**Facts**

The NIHON KOHDEN electrode paste Elefix integrates outstanding sticking qualities with an excellent electrical conductivity. Additionally Elefix has a creamy quality so that the electrodes can be attached firmly to the patient without additional fixation. Elefix does not harden so that the electrodes can be removed without any problems even after longer examinations. The contact impedance between skin and electrode is greatly reduced. In addition, the signal transmission remains stable even under unfavourable external conditions. NIHON KOHDEN Elefix can be used with all surface applied electrodes.

Supply Code	Model	Description
F509	Z-181JE	Elefix EEG Conduction Paste, 180 g tube Packing unit: 1 pcs 

Elefix



F510	Z-401CE	Elefix EEG Conduction Paste, 400 g cup Packing unit: 1 pcs 
------	---------	---

Elefix



Supply Code	Model	Description
NKD-EC2		Conductive Grass Paste EC2 for long-term EEG measurement, 100 g tube Packing unit: 10 pcs 



KOL-1		Collodion Kit, 6.5 ml bottles Packing unit: 20 pcs 
-------	--	---



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.3. ECG Electrodes

Facts

NIHON KOHDEN VITRODE Disposable ECG Electrodes –

The result of over 30 years experience in disposable electrode manufacturing.

The ECG electrodes are the most important portions of the ECG/arrhythmia monitoring system. It is through these electrodes that the electrical signals are conducted from the patient and into the monitoring system. If the electrodes are left in place for an extended period of time (over 48 hours), or they are reused once they are removed, they do not adhere to the patient's skin and they cause "noisy" ECG signals because of the poor conduction. In addition, the patient's movement causes the electrodes to lose conducting capabilities intermittently, and this leads to a "wandering" ECG waveform and muscle artefact. The standard monitoring practice is to change these electrodes at least every 48 hours.

Please contact us for individual offer based on your annual volume! Use the possibility of call orders!

Supply Code	Model	Description
G203	L-150	ECG Electrode Vitrode L, Ag/AgCl, adult, size: 35 mm dia, foam tape, solid gel, general use

Packing unit: 5 x 30 pcs 

Vitrode L Vitrode



G210D	F-150M	ECG Electrode Vitrode F, Ag/AgCl, adult, size M: 25 x 45 mm, cotton tape especially for sensitive skin, adhesive gel
-------	--------	--

Packing unit: 3 x 50 pcs 

Vitrode F Vitrode




G210C	F-150S	ECG Electrode Vitrode F, Ag/AgCl, neonate/child, size S: 18 x 36 mm, cotton tape especially for sensitive skin, adhesive gel
-------	--------	--

Packing unit: 3 x 50 pcs 

Vitrode F Vitrode



NKD-15026N	ECG Snap Lead Wire Kit, 1.5 mm touch proof connector, cable length: 1.5 m, colours: red, yellow, green, black
	For current products: EEG-1200/9100; WEE-1000; PSG-1100
	Packing unit: 4 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.4. Snore Sensors

Snore sensors respond to snoring, picked up by the vibrations of the skin, and convert them to a small analogue voltage that provides a clear, reliable indication of the presence of these sounds. Piezo or dynamic sensors do not require a battery or interface. They are easily activated by plugging the output connector into the appropriate input on the head box.

Facts

The piezo snore sensor consists of a piezo-electric crystal that converts vibrations in the neck caused by surface vibrations from snoring to a small analogue voltage. To enhance the signal quality throughout movements or sweating, the piezo snore sensor has a small knob to guarantee a continuous contact with the patient's skin.

Attaching the sensor

1. The piezo snore sensor should be always fixed to the patient's on lower neck, over trachea.
2. To position the sensor, place the knob side of the piezo snore sensor against the skin and use hypoallergenic medical tape for facial skin to tape the sensor in place. Ensure that the snore sensor is well attached. Otherwise vibrations cannot be clearly transmitted to the sensor during snoring.
3. Fix the wire as a loop.
4. Before starting the measurement, check the signal transmission by asking the patient to simulate snoring.

Supply Code	Model	Description
NKD-SLP1250		Piezo Snore Sensor, size: 2.5 x 1.0 cm, cable length: 2.0 m, 1.5 mm touch proof connector For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs 



16.5. Body Position Sensors


The body position sensor is destined to report sleeping positions of patients during sleep (dorsal, ventral, lateral left, lateral right or sitting) and to diagnose sleep-disordered breathing. The sensor signal that is transmitted to the measurement system is created by gravity evoked switches when patients are changing their sleeping positions.

Facts

Sensor application

1. For an optimal positioning of the sensor, the patient should lie on his/her back.
2. Apply the body position sensor either on the effort sensor belts or the dedicated body position sensor belt. Place the position sensor at the exact centre of the body over the sternum over patients' nightclothes.
3. Tighten the body position belt with caution. Make sure that the belt is not too loose to avoid incorrect measurements during the night.
4. Secure that the body position sensor is still over the sternum.
5. Fix the cable about 12 cm from the sensor with tape to reduce tension and secure a reliable measurement.
6. Connect the output connector to the appropriate input on the measurement system.
7. Ask the patient to show all positions and check that the signals on the PSG system are correctly indicated before leaving the patient.

Supply Code	Model	Description
NKD-SLP1575KIT		DC Body Position Sensor Kit, 3.5 mm connector, incl. 1 x NKD-SLP1575 & 1 x NKD-SLP1550 For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 kit 



NKD-SLP1575



NKD-SLP1550

NKD-SLP1575	DC Body Position Sensor, 3.5 mm connector, length: 3.2 cm, width: 2.7 cm, height: 0.9 cm For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs 
--------------------	--

NKD-SLP1550	Body Position Sensor Belt, length: 1.7 m, width: 2.7 cm, adjustable, black For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs 
--------------------	--

NKD-SLP1582	DC Body Position Sensor, 2.5 mm connector, length: 3.0 cm, width: 2.7 cm, height: 1.2 cm For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 pcs 
--------------------	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.6. Respiratory Effort Sensors

During Polysomnography (PSG) the classification of respiratory events is based on the measurement of the respiratory effort to score the various types of apnoeas (obstructive with a continued respiratory effort, central with no effort, or mixed with mainly no effort followed by obstructive phases and to define an efficient treatment).

Piezo crystal effort sensors contain a piezo crystal that is floating inside the sensor box. Elastic belts around the chest and abdomen measure tensions during respiratory effort. Volume changes are then detected by generated voltage through the floating piezo sensor that is either compressed or stretched. Measurements with piezo crystal effort sensors can have artefacts when patients squeeze the belts during sleep and changing the tension on the piezo crystal sensor or when patients are lying on the sensor.

To overcome such restraints respiratory inductive plethysmography (RIP) belts have wire loops along the belt and therefore generating a magnetic field. As the inductance is precisely encircling the body (measurement area), it changes when the patients breath. The signal is generated based on the inductance changes.

16.6.1. Respiratory Piezo Crystal Effort Sensors

Facts

1. Select the proper size of the belt for the patient – 80% of the body should be covered by the belt, 20% should be done by the elastic part of the belt.
2. The thorax belt should be placed below the armpit. The abdomen belt is placed 2 – 4 cm over the belly button.
3. The sensor box should be placed at the centre of the body.
4. Connect one end of the belt to the sensor, locate the right sensor position and buckle the sensor belts around the patients' body over the nightclothes. The sensor should not touch the skin. Then connect the other belt end to the sensor and tight the band.
5. Check the correct positioning of the belts by asking the patient to cross the arms over the head.
6. Tape the cable approximately 13 cm away from the sensor.
7. Check the correct signals transmission on the EEG system before leaving the patient.

Supply Code	Model	Description
NKD-SLP1370KIT		Piezo Crystal Effort Sensor Kit, 1.5 mm touch proof connector, incl. 1 x NKD-SLP1370 & 1 x NKD-SLP1350 For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 kit 



NKD-SLP1370



NKD-SLP1350

NKD-SLP1370	Piezo Crystal Effort Sensor, Double Buckle, 1.5 mm touch proof connector, size: 13 cm from buckle to buckle, cable length: 2.3 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs 
-------------	--

NKD-SLP1350	Large Buckle Band, clip, width: 2.7 cm, length: 1.55 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs 
-------------	--

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
NKD-SLP1387KIT		<p>Piezo Crystal Effort Sensor Kit, 1.5 mm touch proof connector, incl. 1 x NKD-SLP1387, 1 x NKD-SLP1340 & 1 x NKD-SLP1341</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 kit </p>



NKD-SLP1387



NKD-SLP1341



NKD-SLP1342

NKD-SLP1387	<p>Piezo Crystal Effort Sensor, Double Loop, 1.5 mm touch proof connector, sensor size: 11 cm from loop to loop, cable length: 2.3 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>
--------------------	--

NKD-SLP1341	<p>Large Velstretch Band, width: 4 cm, length: 1.14 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>
--------------------	--

NKD-SLP1342	<p>Small Velstretch Band, width: 2.7 cm, length: 0.73 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>
--------------------	---

Supply Code	Model	Description
NKD-SLP1388KIT		<p>Piezo Crystal Effort Kit, 1.5 mm touch proof connector, incl. 1 x NKD-SLP1388, 1 x NKD-SLP1342 & 1 x NKD-SLP1343</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 kit </p>



NKD-SLP1388



NKD-SLP1342



NKD-SLP1343

NKD-SLP1388	<p>Piezo Crystal Effort Sensor, Double Velcro tabs, 1.5 mm touch proof connector, sensor size: 23 cm from end to end (Velcro tab), cable length: 2.3 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>
--------------------	--

NKD-SLP1342	<p>Small Velstretch Band, width: 2.7 cm, length: 0.73 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>
--------------------	---

NKD-SLP1343	<p>Infant Velstretch Band, width: 2.7 cm, length: 0.36 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>
--------------------	--

Supply Code	Model	Description
NKD-SLP1389KIT		<p>Piezo Film Effort Kit, 1.5 mm touch proof connector, incl. 1 x NKD-SLP1389, 1 x NKD-SLP1340 & 1 x NKD-SLP1341</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 kit </p>



NKD-SLP1389



NKD-SLP1340



NKD-SLP1341

NKD-SLP1389	<p>Piezo Film Effort Sensor, Double Loop, adult, 1.5 mm touch proof connector, sensor size: 15 mm x 37.5 mm, cable length: 2.3 m (from sensor to connectors)</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>
-------------	--

NKD-SLP1340	<p>X-Large Velstretch Band, width: 4 cm, length: 1.73 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>
-------------	---

NKD-SLP1341	<p>Large Velstretch Band, width: 4 cm, length: 1.14 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>
-------------	---

Supply Code	Model	Description
NKD-SLP1391KIT		Piezo Respiration Effort Belt Kit, keyhole connector dark grey, incl. 1 x NKD-SLP1391 & 1 x NKD-SLP1350, cable length: 0.5 m For current products: NF-SCREENER Packing unit: 1 kit 




NKD-SLP1391



NKD-SLP1350

NKD-SLP1391	Piezo Crystal Effort Sensor, Double Loop, keyhole connector, cable length: 0.5 m For current products: NF-SCREENER Packing unit: 1 pcs 
-------------	--

NKD-SLP1350	Large Buckle Band, clip, width: 2.7 cm, length: 1.55 m For current products: NF-SCREENER Packing unit: 1 pcs 
-------------	--

Supply Code	Model	Description	
NKD-SLP1390		<p>Infant Beltless Piezo Effort Sensor, 1.5 mm touch proof connector, sensor size: 2 x 2 cm, cable length: 2.0 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>	
NKD-SLP1340		<p>X-Large Velstretch Band, width: 4 cm, length: 1.73 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>	
NKD-SLP1341		<p>Large Velstretch Band, width: 4 cm, length: 1.14 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>	
NKD-SLP1342		<p>Small Velstretch Band, width: 2.7 cm, length: 0.73 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>	
NKD-SLP1343		<p>Infant Velstretch Band, width: 2.7 cm, length: 0.36 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>	
NKD-SLP1347		<p>Disposable Soft Band, roll, length: 25 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>	
NKD-SLP1350		<p>Large Buckle Band, clip, width: 2.7 cm, length: 1.55 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs </p>	

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.6.2. Respiratory Inductive Plethysmography (RIP) Belts

Facts

Respiratory Inductive Plethysmography (RIP) Belts provide an enhanced patient comfort thanks to their high-quality, soft and flexible material. An easy adjustment to different patients' size is permitted by the Velcro tabs. Sliding and shifting of the belts during night is reduced by the belt width. For easy recognition of the belt sizes, the end pieces are colour-coded. The Inductive Interface Cables are intended to connect the Inductive Plethysmography (RIP) Belts to a PSG System.

Supply Code Model Description

NKD-SLP9002S

Inductive Plethysmography Belt XS, neonate, width: 2 cm, length: 0.4 m, belt colour: white, end piece colour: black

For current products: EEG-1200/9100; WEE-1000; PSG-1100

Packing unit: 2 pcs 



NKD-SLP9003-L40

Inductive Plethysmography Belt S, child, width: 4 cm, length: 0.4 m, belt colour: blue, end piece colour: white

For current products: EEG-1200/9100; WEE-1000; PSG-1100

Packing unit: 2 pcs 



NKD-SLP9003-L60

Inductive Plethysmography Belt M, child, width: 4 cm, length: 0.6 m, belt colour: blue, end piece colour: red

For current products: EEG-1200/9100; WEE-1000; PSG-1100

Packing unit: 2 pcs 



NKD-SLP9003-L90

Inductive Plethysmography Belt L, adult, width: 4 cm, length: 0.9 m, belt colour: blue, end piece colour: dark blue

For current products: EEG-1200/9100; WEE-1000; PSG-1100

Packing unit: 2 pcs 



NKD-SLP9003-L120

Inductive Plethysmography Belt XL, adult, width: 4 cm, length: 1.2 m, belt colour: blue, end piece colour: black

For current products: EEG-1200/9100; WEE-1000; PSG-1100

Packing unit: 2 pcs 



NKD-SLP9003-L150

Inductive Plethysmography Belt XXL, adult, width: 4 cm, length: 1.5 m, belt colour: blue, end piece colour: yellow

For current products: EEG-1200/9100; WEE-1000; PSG-1100

Packing unit: 2 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
NKD-INDUCTKIT		Inductive Interface Cable Kit, incl. 1 x NKD-SLP9101 & 1x NKD-SLP9102 For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 kit R



NKD-SLP9101



NKD-SLP9102


NKD-SLP9101	Inductive Interface Cable, abdomen, 1.5 mm touch proof connector, cable length: 2.3 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs R
--------------------	--

NKD-SLP9102	Inductive Interface Cable, chest, 1.5 mm touch proof connector, cable length: 2.3 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs R
--------------------	--

Supply Code	Model	Description
NKD-SLP9101NO		Inductive Interface Cable for Trackit Sleepwalker, abdomen, keyhole connector, cable length: 0.6 m (from end to end incl. snap) For current products: NF-SCREENER Packing unit: 1 pcs R



NKD-SLP9102NO	Inductive Interface Cable for Trackit Sleepwalker, chest, keyhole connector, cable length: 0.45 m (from end to end incl. snap) For current products: NF-SCREENER Packing unit: 1 pcs R
----------------------	---



16.7. Flow Sensors

Monitoring the airflow existence or absence is an additional PSG parameter besides the respiratory effort to detect apnoeas. Different types of airflow sensors are available:

1. The thermocouple technology is based on two different metal wires that produce a voltage with increasing temperature.
2. Thermistors instead are based on dissimilar metal oxides that show a decreasing resistance when temperature increases, means material changes occur with temperature.

The temperature changes of inhaled air (cools the sensor) and exhaled air (heats sensor) are converted in a signal and present therefore airflow fluctuations. Even if the inspiratory and expiratory fluctuations are recorded by temperature from thermocouple or thermistor, these signals are proportional to the flow.

Pressure transducers, however, are the most sensitive airflow sensor type and are therefore particularly used for hypopnoea detection.

16.7.1. Reusable Thermistors

Facts

Attaching the thermistor flow sensor

1. Place the sensor under the patients' nostrils and above the patients' upper lip.
2. Curve the two nasal prongs away from body sensor. They should not enter the nostrils.
3. Direct the peak of the oral prong to the centre of the mouth/upper lip.
4. Take the wires on both sides and bind them behind the patients' ears. They can be fixed with medical tape.
5. Set the wire tight with the clasp under the chin.

Supply Code	Model	Description
NKD-SLP1468		Reusable Airflow Thermistor, adult, keyhole connector, cable length: 0.9 m For current products: NF-SCREENER Packing unit: 1 pcs 
NKD-SLP1468H		Reusable Airflow Thermistor, adult, keyhole connector, with hanger for cannula, cable length: 0.9 m For current products: NF-SCREENER Packing unit: 1 pcs 
NKD-SLP1469		Reusable Airflow Thermistor, adult, 1.5 mm touch proof connector, cable length: 2.0 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.7.2. Reusable Thermocouples

Facts**Attaching the thermocouple sensor**

1. Place the sensor under the patients' nostrils and above the patients' upper lip.
2. Curve the two nasal prongs away from body sensor. They should not enter the nostrils.
3. Direct the peak of the oral prong to the centre of the mouth/upper lip.
4. Take the wires on both sides and bind them behind the patients' head. They can be fixed with medical tape.
5. Set the wire tight with the clasp under the chin.

Facts


This thermocouple sensor was especially designed for minimal nasal respiration flow detecting of infants and preterm babies. A strong and precise signal is created by attaching the thermocouple on the forehead.

Supply Code	Model	Description
NKD-SLP1400		Reusable Thermal Flow Sensor, infant, 1.5 mm touch proof connector, sensor size: 7 mm x 5.5 cm, cable length: 2.0 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs 



**Facts**

This robust designed paediatric thermocouple sensor can be adjusted to paediatric patients' comfort by bendable prongs.

Supply Code	Model	Description
NKD-SLP1401P		Reusable Thermocouple Flow Sensor, child, 1.5 mm touch proof connector, sensor size: 2.7 x 2 cm, cable length: 2.0 m For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs 



Supply Code	Model	Description
NKD-SLP1401		<p>Reusable Thermocouple Flow Sensor, adult, 1.5 mm touch proof connector, sensor size: 3.5 cm x 4.0 cm, cable length: 2.0 m</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs R</p>
NKD-SLP1401H		<p>Reusable Thermocouple Flow Sensor, adult, 1.5 mm touch proof connector, with hanger for cannula, sensor size: 3.5 cm x 4 cm, cable length: 2.0 m</p> <p>Note: This thermocouple sensor can be used with all kind of adult nasal cannula. Thermal flow and pressure flow can be examined at the same time.</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs R</p>
NKD-SLP1413		<p>Reusable Oral Thermocouple, adult, 1 channel, 1.5 mm touch proof connector, sensor size: 8 cm x 5 mm, cable length: 2.0 m, apply to nasal pressure cannulas</p> <p>For current products: EEG-1200/9100; WEE-1000; PSG-1100</p> <p>Packing unit: 1 pcs R</p>
NKD-SLP1401NOH		<p>Reusable Thermocouple Flow Sensor, adult, keyhole connector, with hanger for cannula, sensor size: 3.5 cm x 4 cm, cable length: 0.75 m</p> <p>For current products: NF-SCREENER</p> <p>Packing unit: 1 pcs R</p>
NKD-SLP1421H		<p>Reusable Thermocouple Flow Sensor, adult, keyhole connector, with hanger for cannula, sensor size: 3.5 cm x 4 cm, cable length: 1.25 m</p> <p>For current products: NF-SCREENER</p> <p>Packing unit: 1 pcs R</p>



16.7.3.

Airflow Pressure Sensors

Supply Code	Model	Description
NKD-SLP14831KIT		AC Airflow Pressure Plus Kit, 1.5 mm touch proof connector, incl. 1 x NKD-SLP14831 & 5 x Adult Nasal Pressure Cannulas of 2.13 m length For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit

NKD-SLP14830KIT		AC Airflow Pressure Sensor Kit, 1.5 mm connector, incl. 1 x NKD-SLP14830 & 5 x Adult Nasal Pressure Cannulas of 2.13 m length For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit
-----------------	--	---



NKD-SLP14830



NKD-SLP14831



NKD-SLP15805

Facts

The airflow determination is realized by applying either an ETCO₂ or a classic pressure cannula to the sensor and the patient. The detected airflow through the nose supplies the nasal airflow signal and the snore signal. This sensor makes use of an internal piezo crystal to produce a small voltage in answer to pressure changes. Thanks to its very small internal volume it is possible that the output signal accurately tracks pressure changes even with small lumen cannulas. The changes are filtered into two bands – slow waves to the flow channel, and fast waves to the snore channel.

Supply Code	Model	Description
NKD-SLP14830		AC Airflow Pressure Sensor, 1.5 mm touch proof connector For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 pcs R



Facts

This sensor is based on the same technology as NKD-SLP14830 but provides 3 sensitivity levels for snore and flow. According to your measurement requirements, different sensitivities of snore and flow can be easily increased or decreased with the side panels of the module.

Supply Code	Model	Description
NKD-SLP14831		AC Airflow Pressure Sensor Plus, adult, 1.5 mm touch proof connector For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 pcs R



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description	
NKD-SLP14835KIT		DC Pressure Sensor Kit (Europe), incl. 1 x DC Pressure Sensor 3.5 mm & 1 x DC Power Supply, 3.5 mm male connector For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit R	
NKD-SLP14834KIT		DC Pressure Sensor Kit (Europe), incl. 1 x DC Pressure Sensor 2.5 mm & 1 x DC Power Supply, 3.5 mm male connector For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit R	
NKD-SLP14834		DC Pressure Sensor, for all patients acc. to cannula size, 2.5 mm male connector For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 pcs R	
NKD-SLP14835		DC Pressure Sensor, for all patients acc. to cannula size, 3.5 mm male connector For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 pcs R	
NKD-SLP14836EU		DC Power Supply (Europe) for DC Pressure Sensor For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 pcs R	
NKD-SLP14837		DC Oesophageal Balloon, Catheter Pressure Sensor, 3.5 mm male connector Note: Please specify with your order if you need a 1.5 mm touch proof connector! For the input box JB-281BK a 1.5 mm touch proof connector is required. For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 pcs D	

D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.7.4.

Pressure Cannulas

Supply Code	Model	Description
NKD-SLP14804AFT-5		Paediatric Nasal Pressure Cannula with Filter, cannula length: 0.6 m For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 5 pcs 



NKD-SLP15805AFT-10		Adult Nasal Pressure Cannula with Filter, cannula length: 0.6 m For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 5 pcs 
--------------------	--	--

NKD-SLP14802AFT-5		Adult Nasal & Oral Pressure Cannula with Filter, cannula length: 0.6 m For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 5 pcs 
-------------------	--	---

NKD-SLP14804BFT-5		Paediatric Nasal Pressure Cannula with Filter, cannula length: 2.13 m For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 5 pcs 
-------------------	--	---

NKD-SLP14805BFT-5		Adult Nasal Pressure Cannula with Filter, cannula length: 2.13 m For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 5 pcs 
-------------------	--	---

NKD-SLP14802BFT-5		Adult Nasal & Oral Pressure Cannula with Filter, cannula length: 2.13 m For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 5 pcs 
-------------------	--	--

NKD-SLP14809		Disposable Cannula Filter, for all NKD-SLPxxxx-Cannulas with Filter For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 5 pcs 
--------------	--	--



NKD-SLP14810		Cannula Adapter For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 5 pcs 
--------------	--	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
NKD-SLPPS036		Male 1.5 mm touch proof connector to Keyhole Adapter For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 1 pcs R



NKD-SLPPS057		Male 1.5 mm touch proof connector to SNAP Adapter For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 2 pcs R
--------------	--	---



NKD-SLP14807		Luer to O ₂ Adapter, for all listed cannulas For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 5 pcs D
--------------	--	---



NKD-SLP14808		Male to Male Luer Adapter, for all above listed cannulas For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 5 pcs D
--------------	--	--



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.7.5. Pressure Transducer

Facts

The pressure transducer is used for detection and imitation of waveforms related to respiration, snoring and further pressure signals. It is intended for use on children (more than two years) and adults.

This device does not need a calibration before use but a battery check is possibly by doing a self-test (push the self-test button). The self-test also signifies a proper function of the electronic circuit and the connection to the PSG system. A continuous green light on the top of the pressure transducer indicates no need of battery change and typical waveforms can be detected from the three outputs.

The snore output signal on the left side is detected by the snore vibrations on the nasal pressure signal. The channel has an internal 15 Hz low frequency filter to give a flat baseline between snores by taking out the airflow signal.

An unfiltered and filtered output signal can be chosen. Whereas the unfiltered nasal pressure airflow waveform is overlaid with snoring, the filtered nasal pressure airflow output signal makes use of an internal 2.5 Hz high frequency filter to eliminate this.

Once the pressure transducer is ready for use onto the patient, you can use the "Gain Set" function which is an average output level adjustment for each patient. Before putting its low level airflow outputs to about 600 μV peak to peak and the high level airflow signals to 2 volts peak to peak, some breaths of the patient are examined. During this setting the LED on top of the pressure transducer is all the time green. A manual regulation of the output levels can be done when using the up/down output buttons (indicated by the according flash directions).

The pressure transducer is a perfect tool to spot pressure related signals.

Supply Code	Model	Description
NKD-PTRANS		Pressure Transducer PTAF2, AC/DC output, incl. 1 x Sample Cannula Kit, 1 x NKD-PT1075 & 1 x NKD-PT1080, 9 V alkaline battery, width: 7 cm, length: 12.6 cm, height: 2.4 cm, weight: 130 g (incl. battery), Input pressure range: +/- 25.4 cm H ₂ O, high level output (typical): 2 V peak to peak, low level output (typical): 600 μV peak to peak For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit 



NKD-PT1075		Replacement Cable for AC Snore/AC Airflow Output, for NKD-PTRANS, 1.5 mm touch proof connector For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 pcs 
------------	--	--



NKD-PT1080		Replacement Cable for DC Airflow Output for NKD-PTRANS, 2.5 mm male connector For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 pcs 
------------	--	---



NKD-PT1327		CPAP Titration Kit for NKD-PTRANS, incl. 1 x CPAP-T-Piece & 1 x NKD-PT1109 For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit
------------	--	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
NKD-PT1107		Disposable Cannula Titration Filter Kit for NKD-PT1327, incl. 50 cannulas For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit 



NKD-PT1109		Disposable Cannula Titration Filter Kit for NKD-PT1327, incl. 10 cannulas For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit 
------------	--	---

NKD-PT1295		Disposable Adult Nasal & Oral Cannula with Filter for NKD-PTRANS, cannula length: 2.13 m, incl. 30 cannulas For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit 
------------	--	---

NKD-PT1263		Disposable Paediatric Nasal Cannula with Filter for NKD-PTRANS, cannula length: 2.13 m, incl. 60 cannulas For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit 
------------	--	--

NKD-PT1267		Disposable Paediatric Small Nasal Cannula with Filter for NKD-PTRANS, cannula length: 2.13 m, incl. 60 cannulas For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit 
------------	--	---

NKD-PT1354		Disposable Paediatric Nasal & Oral Cannula with Filter for NKD-PTRANS, cannula length: 2.13 m, incl. 30 cannulas For current products: EEG-1200/9100; WEE-1000 Packing unit: 1 kit 
------------	--	--

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.8. Limb Movement Sensors

Supply Code	Model	Description
NKD-SLP1770KIT		Limb Movement Sensor Kit (2 limbs), adult, hand/ankle/foot, 1.5 mm touch proof connector, incl. 1 x NKD-SLP1770 & 1 x NKD-SLP1730 For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 kit 

NKD-SLP1772KIT		Limb Movement Kit (2 limbs), adult, hand/ankle/foot, 1.5 mm touch proof connector, incl. 2 x NKD-SLP1770 & 1 x NKD-SLP1730 For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 kit 
----------------	--	--



NKD-SLP1770



NKD-SLP1730

NKD-SLP1771KIT		Limb Movement Sensor Kit (2 limbs), adult, hand/ankle/foot, keyhole connector, incl. 2 x NKD-SLP1771 & 1 x NKD-SLP1730 For current products: NF-SCREENER Packing unit: 1 kit 
----------------	--	--



NKD-SLP1771



NKD-SLP1730

Facts

An internal piezo film element detects every movement of the hand or leg. The created signal is transformed to a small analogue voltage providing precise and reliable information of limb movements. The sensor may be connected to any AC (low level) channel on the recorder.

Attaching the limb movement sensor

1. The piezo limb movement sensor is placed on the wrist or ankle with the limb movement cross band kit.
2. For the ankle, please apply the sensor on the front of the foot over the tendons. The cross band shall be attached over the heel. For the wrist, fix the sensor on the inside of the wrist, over the tendons. The band can be overlapped.
3. Fix the band with the Velcro side of the sensor.
4. Fix the cable with medical tape about 15 cm from the sensor.
5. Make sure that the sensor will not move during examination and that the signal is transmitted.

Supply Code	Model	Description
NKD-SLP1770		Piezo Limb Movement Sensor without Cross Band (2 limbs), adult, hand/ankle/foot, 1.5 mm touch proof connector, cable length: 2.9 m, sensor size black Velcro: 16 cm, gray silicone: 4.5 cm, band length: 29 cm For current products: EEG-1200/9100; WEE-1000; PSG-1100 Packing unit: 1 pcs 



NKD-SLP1771		Limb Movement Sensor without Cross Band, adult, hand/ankle/foot, keyhole connector, cable length: 2.0 m, sensor size black Velcro: 16 cm, gray silicone: 4.5 cm, band length: 29 cm, cable colour: white For current products: NF-SCREENER Packing unit: 1 pcs 
-------------	--	--



Supply Code	Model	Description
NKD-SLP1730		Limb Movement Cross Band Kit for NKD-SLP1770KIT/-SLP1771KIT/-SLP1772KIT, band length: 29 cm For current products: EEG-1200/9100; WEE-1000; PSG-1100; NF-SCREENER Packing unit: 1 kit 



NKD-SLP1726		Y-Splitter Cable for NKD-SLP1771, connects 2 leg sensors to 1 system input For current products: NF-SCREENER Packing unit: 1 pcs 
-------------	--	--



NKD-SLP1326		Extension Cable for Inductive Interface Cable with Keyhole Connector, cable length: 0.5 m For current products: NF-SCREENER Packing unit: 1 pcs 
-------------	--	---

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.9. NIHON KOHDEN SpO₂ Measurement



A Pulse Oximeter is a non-invasive device used to percutaneously measure oxygen saturation in the blood by use of a clip-on finger probe. Pulse oximeter works with the integration of a probe placed onto patient's fingertip or toe, obtains signals from the probe and the main device calculates oxygen saturation in the body.

Oxygen saturation measured by a Pulse Oximeter is SpO₂. "P" stands for "pulse", and the measurement principle is based on beat-by-beat changes in blood volume. In general, arterial oxygen saturation measured by spectroscopic analysis is called "SaO₂", and oxygen saturation measured by pulse oximeter is called "SpO₂".

The method for determination of arterial oxygen saturation, based on the photoelectric plethysmography, was conceptualised and implemented in 1972 by the Japanese bio-engineer Takuo Aoyagi (NIHON KOHDEN Corporation). This idea came about while using an ear densitometer for recording dilution curves for dyes (pigments). The first commercial oximeter "OLV-5100" was presented by NIHON KOHDEN back in 1975.

Dr. Aoyagi's innovative technology, pulse oximetry, has rapidly spread worldwide and become essential in critical care. In 2002, Dr. Aoyagi was awarded the Medal of Honor in Japan, in recognition of his great contribution to clinical medicine. The NIHON KOHDEN SpO₂ Sensors of the BluPRO® Series are at the cutting edge of this progressive technology.



16.9.1. NIHON KOHDEN SpO₂ Connection Cables

Supply Code	Model	Description
Y095A	JL-550T2	SpO ₂ Connection Cable, cable length: 2.5 m, round connector For current products: TEC-5500/7700; JE-912A/921A/921AG Packing unit: 1 pcs 



K931	JL-900P	SpO ₂ Connection Cable, cable length: 2.5 m, rectangular connector For current products: all BSM; PVM; TEC-8300; PSG-1100 Packing unit: 1 pcs 
------	---------	--



Note

SpO₂ measurement starts when the SpO₂ connection cable and SpO₂ probe are connected to the device, and the SpO₂ probe is attached to the patient.

When not measuring SpO₂, disconnect the SpO₂ connection cable from the device. Otherwise, noise from the probe sensor may interfere and incorrect data is displayed on the screen.

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.9.2. NIHON KOHDEN SpO₂ Probes

BluPRO® SpO₂ Probes – A Unique Sensing Technology

Reliability and high quality of our probes provide a continuous accurate and safe monitoring with highest patient comfort. The com-fortable, sure fitting and lightweight characteristic reduces the effect of hand movement for highly reliable SpO₂ measurement.

This ensures patient safety, reduces the costs of care and increases therefore clinical efficiency and cost effectiveness. NIHON KOHDEN BluPRO® SpO₂ Probes are comfortable, durable, accurate and needs-oriented.

Note that accurate and safe measurement is assured only when you use probes, connection cords and parts manufactured by NIHON KOHDEN. Some probes in the market may advertise as „compatible with NIHON KOHDEN probes.“ We strongly recommend you to avoid these probes as a correct measurement cannot be guaranteed on NIHON KOHDEN devices. NIHON KOHDEN is not responsible for problems caused by using other manufacturer’s probes.

Please contact us for an individual offer based on your annual volume! Use the possibility of call orders!

	Supply Code (Model)	Cable length (m)	Packing unit	Weight							
				Low birth weight	Neonate	Infant	Child	Adult			
				1000 g	2 kg	3 kg	10 kg	20 kg	30 kg	40 kg	50 kg
Disposable	 P203A (TL-271T)	0.8	24								
	P203E (TL-271T3)	1.6									
	 P204A (TL-271T)	0.8	5								
	P204E (TL-271T3)	1.6									
	 P203B (TL-272T)	0.8	24								
	P203F (TL-272T3)	1.6									
	 P204B (TL-272T)	0.8	5								
	P204F (TL-272T3)	1.6									
	 P203C (TL-273T)	0.8	24								
	P203G (TL-273T3)	1.6									
	 P204C (TL-273T)	0.8	5								
	P204G (TL-273T3)	1.6									
	 P203D (TL-274T)	0.8	24								
	P203H (TL-274T3)	1.6									
	 P204D (TL-274T)	0.8	5								
	P204H (TL-274T3)	1.6									
 P228A (TL-051S)	0.8	5									
P228B (TL-052S)	1.6										
 P229A (TL-061S)	0.8	5									
P229B (TL-062S)	1.6										
 P205A (TL2601)	1.6	5			tape S	tape L	tape				
 P260A (tape S)		24			tape S	tape L	tape				
P260B (tape L)											
Reusable	 P225H (TL-201T1)	0.6									
	 P225G (TL-220T)	1.6	1								
	 P311C (TL-631T3)	1.6									
	P310C (TL-630T3)	1.6									

16.9.2.1. Reusable SpO₂ Probes **BluPRO®**

Reusable, but always clean – NIHON KOHDEN washable probes.

The reusable probes of the BluPRO® series are durable and washable, so you can clean them under running water or even immerse them in a soapy alkaline or disinfectant solution (glutaraldehyde solution 2.0 % or hydrochloric alkyl diaminoethylene glycine 0.5 %)!



P225F/P225H



P225G



P310C/P311C

Durable, washable, demand-oriented!

Facts

NIHON KOHDEN reusable SpO₂ finger-clips probes have parallel mounted upper and lower sensor shells, thus, providing an exact alignment of the optical sensor components and markedly enhanced comfort during long-term attachment. In addition, flexible lateral flaps shield the measuring signals from the external light, thus resulting to minimize the interference by additional artificial light sources! Bright LEDs provide more accurate measurement even for thicker fingers or dark skin tones. The finger-clip reusable probe is the most popular type of SpO₂ probes. Instead of the typical hinge near the finger tip, P225F/P225H have a unique sure fitting and comfortable mechanism using two sliding pieces for the front and back of the finger with two rubber flaps along both sides of the finger. When you press the button, the two halves of the probe slide directly open, parallel to each other. No additional fixation tape is needed. For a stable and comfortable measurement P225F/P225H, always keep the light emitters and detector aligned. NIHON KOHDEN provides unique reusable SpO₂ probes.

Supply Code	Model	Description
P225F	TL-201T	SpO ₂ BluPRO® Finger-Clip Probe, child/adult (> 20 kg) finger, cable length: 1.6 m, reusable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 1 pcs R



P225H	TL-201T1	SpO ₂ BluPRO® Finger-Clip Probe, child/adult (> 20 kg) finger, cable length: 0.6 m, reusable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 1 pcs R
-------	----------	---

Facts

The multi-site probe P225G can be used for neonatal feet as well as adult fingers as its curved contact surface is designed to fit on different attachment sites. The light emitter and detector are separate which allows especially an accurate alignment of the LED on the infant/neonatal instep. A unique sliding adjuster gives the best fit on the attachment site secured by taping.

Supply Code	Model	Description
P225G	TL-220T	SpO ₂ BluPRO® Multi-Site Probe, neonate (< 3 kg) instep, child/adult (> 3 kg) finger/toe, cable length: 1.6 m, reusable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 1 pcs R



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Facts

Our special designed finger-tip probes, P310C and P311C, fit immediately on the finger and are perfect for long-term measurement. A less covered finger area and cushioned materials lead to decrease pressure to the attachment site. The 3D molding attachment ensures a comfortable fit.

In combination with the soft and lightweight design the patient stress or discomfort is reduced and the probes can be remained up to 8 hours to one attachment site. A stable measurement is assured by gentle fixation tapes P263 or P267 that avoid the falling off of the probes and reduce artifact from cable movement.

Reusable but with the comfort of a disposable probe, P310C and P311C are also applicable for sleep studies, for patients that are often bending their fingers or obese patients.

Supply Code	Model	Description
P310C	TL-630T3	SpO ₂ BluPRO® Finger-Tip Probe L, child/adult (> 50kg) finger, cable length: 1.6 m, long-term use, reusable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 1 pcs R
P311C	TL-631T3	SpO ₂ BluPRO® Finger-Tip Probe M, child/adult (> 20kg) finger, cable length: 1.6 m, long-term use, reusable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 1 pcs R
P263		Disposable Probe Cable Fixing Tape for P310C/P311C, adhesive, size: 15 x 130 mm For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 3 x 30 pcs D
P267	YS-093P2	Reusable Probe Fastener Tape for P310C/P311C, adhesive, size: 15 x 140 mm For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 30 pcs R



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.9.2.2. Disposable SpO₂ Probes **BluPRO**[®]

NIHON KOHDEN disposable probes, from neonate to adult sizes, place special emphasis on the exact positioning of the LEDs and detector. For this reason the two optical components are designed separately and can thus be brought into precise alignment with each other by the user. All probes have breathable material and can be remained up to 8 hours on one attachment site. Paediatric probes have a cute cartoon character printed on the tape to reassure children.

Self-adhesive, soft and thin tape provides a secure fit. The special design allows an easy application to the attachment site. To get the most comfortable fit, please bend the probe before attaching, do not strongly pull the attachment tape when attaching and do not push the tip of the finger with the attachment tape.

16.9.2.2.1. Disposable SpO₂ Probes 24 pcs/box – cable length: 0.8 m

Supply Code	Model	Description
P203A	TL-271T	SpO ₂ BluPRO [®] Probe, adult (> 30 kg) finger/toe, cable length: 0.8 m, disposable For current products: all BSM/PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 24 pcs 
P203B	TL-272T	SpO ₂ BluPRO [®] Probe, child (10 – 50 kg) finger/toe, cable length: 0.8 m, disposable For current products: all BSM/PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 24 pcs 
P203C	TL-273T	SpO ₂ BluPRO [®] Probe, neonate (< 3 kg) instep, adult (> 40 kg) finger/toe, cable length: 0.8 m, disposable For current products: all BSM/PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 24 pcs 
P203D	TL-274T	SpO ₂ BluPRO [®] Probe, infant (3 – 20 kg) finger/toe, cable length: 0.8 m, disposable For current products: all BSM/PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 24 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.9.2.2.2. Disposable SpO₂ Probes 24 pcs/box – cable length: 1.6 m

Supply Code	Model	Description
P203E	TL-271T3	SpO ₂ BluPRO® Probe, adult (> 30 kg) finger/toe, cable length: 1.6 m, disposable For current products: all BSM/PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 24 pcs 



P203F	TL-272T3	SpO ₂ BluPRO® Probe, child (10 – 50 kg), finger/toe, cable length: 1.6 m, disposable For current products: all BSM/PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 24 pcs 
-------	----------	---



P203G	TL-273T3	SpO ₂ BluPRO® Probe, neonate (< 3 kg) instep, adult (> 40 kg) finger/toe, cable length: 1.6 m, disposable For current products: all BSM/PVM/TEC/EEG/PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 24 pcs 
-------	----------	---



P203H	TL-274T3	SpO ₂ BluPRO® Probe, infant (3 – 20 kg) finger/toe, cable length: 1.6 m, disposable For current products: all BSM/PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 24 pcs 
-------	----------	--

16.9.2.2.3. Disposable SpO₂ Probes 5 pcs/box – cable length: 0.8 m

Supply Code	Model	Description
P204A	TL-271T	SpO ₂ BluPRO® Probe, adult (> 30 kg) finger/toe, cable length: 0.8 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 



P204B	TL-272T	SpO ₂ BluPRO® Probe, child (10 – 50 kg) finger/toe, cable length: 0.8 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 
-------	---------	--



P204C	TL-273T	SpO ₂ BluPRO® Probe, neonate (< 3 kg) instep, adult (> 40 kg) finger/toe, cable length: 0.8 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 
-------	---------	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description
P204D	TL-274T	SpO ₂ BluPRO® Probe, infant (3 – 20 kg) finger/toe, cable length: 0.8 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 



Facts

Absolutely cost-effective, the series P228 and P229 are intended for long-term measurement on a single patient. The probes do not need to be replaced with a new one before 96 hours by changing the attachment site every 8 hours. So your running costs can be reduced!

The light emitter and detector are covered by sponge and the entire attachment site has a flat surface. This reduces the risk of skin burn and irritation.

Since the probe itself is not self-adhesive, you can easily remove and attach the probe to change the attached finger. For fixation, wrap the probe around the finger with attachment tape P259.

P228A	TL-051S	SpO ₂ BluPRO® Probe, neonate (< 3 kg) instep, adult (> 50 kg), finger/toe, cable length: 0.8 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 
-------	---------	--



P229A	TL-061S	SpO ₂ BluPRO® Probe, infant/child (3 – 15 kg) toe, child (15 – 50 kg) finger, cable length: 0.8 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 
-------	---------	---



P259		COTTONY Tape for fixing disposable SpO ₂ Probes, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 20 pcs 
------	--	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.9.2.2.4. Disposable SpO₂ Probes 5 pcs/box – cable length: 1.6 m

Supply Code	Model	Description
P204E	TL-271T3	SpO ₂ BluPRO® Probe, adult (> 30 kg) finger/toe, cable length: 1.6 m, disposable For current products: all BSM; PVM; TEC; EEG PSG-1100; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 



P204F	TL-272T3	SpO ₂ BluPRO® Probe, child (10 – 50 kg), finger/toe, cable length: 1.6 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 
-------	----------	--



P204G	TL-273T3	SpO ₂ BluPRO® Probe, neonate (< 3 kg) instep, adult (> 40 kg) finger/toe, cable length: 1.6 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 
-------	----------	---



P204H	TL-274T3	SpO ₂ BluPRO® Probe, infant (3 – 20 kg) finger/toe, cable length: 1.6 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 
-------	----------	---



Facts

Absolutely cost-effective, the series P228 and P229 are intended for long-term measurement on a single patient. The probes do not need to be replaced with a new one before 96 hours by changing the attachment site every 8 hours. So your running costs can be reduced! The light emitter and detector are covered by sponge and the entire attachment site has a flat surface. This reduces the risk of skin burn and irritation. Since the probe itself is not self-adhesive, you can easily remove and attach the probe to change the attached finger. For fixation, wrap the probe around the finger with attachment tape P259.

P228B	TL-052S	SpO ₂ BluPRO® Probe, neonate (< 3 kg) instep, adult (> 50 kg) finger, cable length: 1.6 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 
-------	---------	---



P229B	TL-062S	SpO ₂ BluPRO® Probe, infant/child (3 – 15 kg) toe, child (15 - 50 kg) finger, cable length: 1.6 m, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 5 pcs 
-------	---------	---



P259		COTTONY Tape for fixing disposable SpO ₂ Probes, disposable For current products: all BSM; PVM; TEC; EEG; PSG-1100; ZS-930/530/940PG; ZM-930/940PG Packing unit: 20 pcs 
------	--	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.9.2.2.5. Disposable Multi-Site SpO₂ Probes 5 pcs/box – cable length: 2.0 m

Facts

The probe P205A is specially designed for extremely sensitive skin. It uses soft and breathable sponge for the attachment tape. The attachment tape and sensor parts are separate. You can select from two sizes of attachment tape: P260A or P260B. When using the probes without those tapes for prematures, SpO₂ measurement may be incorrect and the skin at the attachment site may be irritated.

The attachment tape consists of two parts: adhesive and non-adhesive part. To get best fit, place the tape at the required attachment site. Wrap the sponge attachment tape (non-adhesive part) around the attachment site. Confirm that the light emitters and detector face to each other before removing backing paper of the adhesive part. Then continue wrapping to fasten the probe. Replace the sponge attachment tapes every 24 hours.

Supply Code	Model	Description
P205A	TL-260T	SpO ₂ BluPRO Multi-Site Probe, cable length: 2.0 m, neonate (< 3 kg) instep, child/adult (> 3kg) finger/toe, fixed with P260A/B or for P256, single-patient use For current products: all BSM; PVM; TEC; EEG Packing unit: 5 pcs 



P260A	6143-011854	Sponge Attachment Tape Type S for P205A, premature (< 1 kg) instep, neonate/infant (> 3 kg) finger/toe, disposable For current products: all BSM; PVM; TEC; EEG; ZS-930/530/940PG; ZM-930/940PG Packing unit: 24 pcs 
-------	-------------	--



P260B	6143-011863	Sponge Attachment Tape Type L, neonate (> 3 kg) instep, disposable For current products: all BSM; PVM; TEC; EEG; ZS-930/530/940PG; ZM-930/940PG Packing unit: 24 pcs 
-------	-------------	--



Handling of the probe affects measurement accuracy and patient safety



For the proper placement of the probe on the finger, it is advisable that the light emitter and detector are placed just opposite to each other. When using a finger-clip probe, it is necessary to ensure that the probe covers the finger properly, otherwise some portion of the light never passes through the body tissues and may reach the light-receiving part directly, leading to incorrect measurement results.

The more distant the light-emitter and detector gets, the weaker the transmitted light intensity becomes. Additionally, a thinner body part tends to have small amount of blood, which makes it difficult to extract the arterial blood signal component. The most suitable thickness of the probe site is about 10 mm.

Too loose probes may interfere with alignment of LEDs and detector, come off easily by body movement, cause larger motion artifacts or be effected by venous pulse.

Too tight probes may cause venous pulsation, cause skin injuries such as burn and necrosis or unstable measurement because of a lower pulsatile signal.

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.10. NONIN SpO₂ Probes

Supply Code
NF-N/8000AA

Model

Description

NONIN reusable SpO₂ Finger-Clip Sensor, child/adult (> 30 kg), cable length: 1.0 m

For current products: NF-SCREENER; WEE-1000

Packing unit: 1 pcs 



NF-N/8000J

NONIN reusable SpO₂ Flex Sensor, child/adult (> 20 kg), cable length: 1.0 m

For current products: NF-SCREENER; WEE-1000

Packing unit: 1 pcs 



NF-N/8008J

NONIN reusable SpO₂ Flex Sensor, infant/child (2 - 20 kg), cable length: 1.0 m

For current products: NF-SCREENER; WEE-1000

Packing unit: 1 pcs 



JL-101A

SpO₂ Adapter for NONIN SpO₂ Sensors, 3.5 mm connector, cable length: 2.0 m

For current products: WEE-1000

Packing unit: 1 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.11. ETCO₂ Measurement (Capnography)




Reasons to measure ETCO₂ (the level of carbon dioxide in the air exhaled from the body – normal value is about 38 mm Hg):
Capnography helps to detect life-threatening problems in:

- Airway
- Breathing
- Circulation and
- Metabolism

Capnography is a fast and reliable help to determinate respiratory difficulties, wrong connections of airway adapters, leaks in the ventilator circuit due to its placement in the respiratory circuit. An early detection by capnography might help to avoid inadequate oxygen supply of the body.

The American Academy of Sleep Medicine (AASM) recommends monitoring hypoventilation as standard parameter. For this purpose the continuous measurement of arterial PCO₂ by transcutaneous PCO₂ or end-tidal PCO₂ for detecting hypoventilation is advised during sleep studies. Both, transcutaneous and end-tidal PCO₂ are reflecting correlations to the arterial PCO₂.


ETCO₂ shows clinical and technical problems




Breathing circuit abnormality

Respiratory abnormality


Trend




5 phases of the capnograph




Conditions reflected by ETCO₂



Early detection of trouble in breathing or ventilation equipment



ETCO₂ and PaCO₂ show ventilation and pulmonary blood flow condition



CO₂ Measurement is useful in many hospital sites

Nihon Kohden Unique Sensing Technology



Both intubated and non-intubated cap-ONE CO₂ Sensor

16.11.1. Intubated Expired and Inspired ETCO₂ “Opening a New Era in Capnography”

cap-ONE
ORAL NASAL EXPIRATION

Facts

NIHON KOHDEN CO₂ Sensor Sensation – P909 World’s smallest and lightest mainstream ETCO₂ sensor (quantitative method) with only 4 g (±2 g). This avoids tube pressures and increases patient comfort. Especially for this reason as well as the dead space of approximately 1 ml it is perfectly suited to the most sensible intubated patients – neonates in critical care (application on neonates only with patient monitors). Besides this P909 can be also applied to children and adults. Easy and simple handling is guaranteed by its waterproofness (IPX7 standard) and our Smart Cable™ Technology so that no additional module is needed. Just plug in and start a continuous measurement after only 10 seconds warm-up time. Those features make the P909 an unbeatable help for daily work routine in OR, ER, ICU or Ward.

Supply Code	Model	Description
P909	TG-970P	cap-ONE CO ₂ Sensor Kit for Intubated Patients, quantitative measurement, incl. P924/K987, height: 8 mm, width: 13 mm, length: 35 mm, cable length: 3.5 m For Airway Adapters: R805/R806 For current products: all BSM; PVM-2703; OLG-2800; TEC-8300; PSG-1100 Packing unit: 1 pcs 



R805	YG-211T	Airway Adapter for P909, intubated, quantitative method, adult (> 7 kg), dead space: 4 ml For current products: all BSM; PVM-2703; OLG-2800; TEC-8300; PSG-1100 Packing unit: 30 pcs 
------	---------	--



R806	YG-213T	Airway Adapter for P909, intubated, quantitative method, neonate/child (2-7 kg), dead space: 1 ml For current products: all BSM; PVM-2703; OLG-2800; PSG-1100 Packing unit: 30 pcs 
------	---------	--



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

cap-ONE
ORAL NASAL EXPIRATION

16.11.2. **Intubated and Non-Intubated Expired ETCO₂ CapOne**
“Opening a New Era in Capnography”

Facts

NIHON KOHDEN cap-ONE – when every breath counts!

World’s first wearable mainstream ETCO₂ sensor P907 (semi-quantitative method) for intubated and non-intubated patients which is applicable for children and adults.

Unique oral-nasal airway adapters catch oral as well as nasal expired CO₂ – this increases accuracy. The cap-ONE is not only small but has also a mini-weight of only 4 g. This avoids tube pressures and increases patient comfort. Easy and simple handling is guaranteed by our Smart Cable™ Technology so that no additional module is needed, just plug in and start measuring after only 5 seconds. A continuous, stable and accurate measurement is ensured due to no lag time and no need to change tubes because of occlusions. Its waterproofness and outstanding shock resistance are an unbeatable advantage in daily work routine.

Those features make the cap-ONE P907 an excellent measurement tool in Sleep Labs, OR, ICU, Recovery Room, PACU, Transport, ER, Cardiac Catheter Laboratory, GI Endoscopy, Pain Management, Abdominal Operation, Obstetric/Orthopedic/Eye Surgery, Plastic or Dental Surgery.



Supply Code	Model	Description
P907	TG-920P	cap-ONE CO ₂ Sensor Kit for Intubated/Non-Intubated Patients, semi-quantitative measurement, incl. P923/K984, size: 11 x 11 x 40 mm, cable length: 3.5 m For Airway Adapters: Intubated patients: R804 Non-intubated patients: V928/V929 with flow measurement; V921/V922/V923 For current products: all BSM; PVM-2703; OLG-2800; TEC-8300; PSG-1100 Packing unit: 1 pcs 
P908	TG-921T3	cap-ONE CO ₂ Sensor Kit for Intubated Patients, quantitative measurement, incl. P924/K987, cable length: 3.5 m For Airway Adapters: Intubated patients: R804 Non-intubated patients: V928/929 with flow measurement; V921/V922/V923 For current products: TEC-5500/7700; EEG-1200/9100/9200 with JE-921A/921AG Packing unit: 1 pcs 
R804	YG-111T	Airway Adapter for P907/P908/P923, intubated, child/adult (> 7 kg), dead space: 4 ml for P907: all BSM; PVM-2703; TEC-8300; OLG-2800; PSG-1100 for P908: TEC-5500/7700 For current products: all BSM; PVM-2703; OLG-2800; TEC-5500/7700/8300; PSG-1100 Packing unit: 30 pcs 



 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Supply Code	Model	Description	
V921	YG-120T	CO ₂ Nasal Adapter for Nasal Breathing, non-intubated, child/adult (> 10 kg) for P907: all BSM; PVM-2703; TEC-8300; OLG-2800; PSG-1100 for P908: TEC-5500/7700; EEG-1200/9100/9200 with JE-921A/921AG For current products: all BSM; PVM-2703; OLG-2800; TEC-8300; EEG-1200/9100/9200 with JE-921A/921AG; PSG-1100 Packing unit: 30 pcs 	
V922	YG-121T	CO ₂ Nasal Adapter for Nasal/Oral Breathing, non-intubated, child/adult (> 10 kg) for P907: all BSM; PVM-2703; TEC-8300; OLG-2800; PSG-1100 for P908: TEC-5500/7700; EEG-1200/9100/9200 with JE-921A/921AG For current products: all BSM; PVM-2703; OLG-2800; TEC-8300; EEG-1200/9100/9200 with JE-921A/921AG; PSG-1100 Packing unit: 30 pcs 	
V923	YG-122T	CO ₂ Nasal Adapter for Nasal/Oral Breathing with Oxygen Cannula Holder, non-intubated, child/adult (> 10 kg) for P907: all BSM; PVM-2703; TEC-8300; OLG-2800; PSG-1100 for P908: TEC-5500/7700; EEG-1200/9100/9200 with JE-921A/921AG For current products: all BSM; PVM-2703; OLG-2800; TEC-8300; EEG-1200/9100/9200 with JE-921A/921AG; PSG-1100 Packing unit: 30 pcs 	
V928	YG-125T	Adult Nasal/Oral CO ₂ Adapter for Flow Measurement, child/adult (> 10 kg), dead space: 1.2 ml for P907: PSG-1100 for P908: EEG-1200/9100/9200 with JE-921A/921AG Note: This adapter is to be used for sleep studies for adults and includes an integrated pressure transducer cannula. This adapter needs to be connected to the pressure transducer NKD-PTRANS. For current products: EEG-1200/9100/9200 with JE-921A/921AG; PSG-1100 Packing unit: 10 pcs 	
V929	YG-135T	Paediatric Nasal/Oral CO ₂ Adapter for Flow Measurement, child (> 7 kg), dead space: 1.2 ml for P907: PSG-1100 for P908: EEG-1200/9100/9200 with JE-921A/921AG Note: This adapter is to be used for sleep studies and has a smaller mouth shape for children. It includes an integrated pressure transducer cannula and needs to be connected to the pressure transducer NKD-PTRANS. For current products: EEG-1200/9100/9200 with JE-921A/921AG; PSG-1100 Packing unit: 10 pcs 	

 Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

 Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

16.12. Sensor Tester

Facts

The continuity impedance tester for electrodes, sensor cables, thermocouples etc. is recommended for each sleep laboratory. No calibration is necessary. The tester needs a 9 Volt battery. For testing disconnect sensors and equipment from patient before testing. Never connect it to any electrical power.

Technical features:

1. Hi/Lo/Off switch

- Hi: is for testing thermocouples through the 1.5 mm touch proof connectors.
- Lo: is for all other testing.
- Off: turns power to the unit off and preserves battery life.

2. Sensor test

The sensor test can be done with sensor having 1.5 mm touch proof connectors. Two inputs are available on the sensor tester. The sensor tester can be used for respiratory effort sensors, thermocouples, piezo snore sensors or piezo limb movement sensors.

3. Continuity test

For this test a cable with 1.5 mm touch proof connector and alligator clip is required for testing single wire electrodes to find cable breaks. In the case the electrode is not broken (positive continuity) the sensor tester is giving a sound. If no tone is audible, a cable break is possible.

4. Square wave

The sensor tester can be used to check cable connections and head boxes by sending the square wave output signal into the PSG system. Please note that the required output cables are not included because of various inputs of the different data acquisition units. Please choose the lowest low frequency filter of the PSG system to examine the square wave at its best resolution. Usually the amplitude should be +/- 50 μ V whereas 0 – 1 volt output will be visible with high level DC output cables.

5. Display

The LED display shows a bar chart. Any stimulation of the attached sensors will display a change of the bar. When extending an effort sensor, the light will swing accordingly or it will move quickly when a snore sensor is slightly knocked off.

Supply Code	Model	Description
NKD-PTST3		Continuity Tester Kit for Electrodes & Wires, incl. 1 x Sensor Tester & 1 x NKD-PT1315, 9 V alkaline battery, height 2.4 cm, width: 7 cm, length: 12.6 cm, weight: 130 g, square wave signal output: 100 μ V peak-to-peak, polarity Packing unit: 1 kit R



Sensor Tester



NKD-PT1315

NKD-PT1315	Wire Kit for NKD-PTST3 Packing unit: 1 kit R
------------	---



D Disposable accessories are intended for single use only, on one patient, during a single examination! Do not reuse disposable accessories!

R Reusable accessories are intended for repeated patient use on different patients during different examinations. Please make sure to clean reusable accessories after each patient or between different applications to prevent or reduce the risk of cross infections!

Index as per Supply Code

Supply Code	page	Supply Code	page	Supply Code	page	Supply Code	page
1144-008488C	63	H542B	48	H662A	33	K516B	49
2114-073312	63	H542C	48	H663	34	K518A	49
2114-073401	63	H543A	57	H664	34	K518B	49
2114-080955	26	H543B	57	H664A	34	K611	13, 14
A720B	40, 80	H551A	63	H665	17	K612	15
A725	40	H561	48	H666	17	K616	17
A726	40	H564A	63	H667	17	K617	20, 41
BE-81005	63	H568	63	H668	17	K617A	20, 41
C451	80	H569	48	H669	17	K617B	20, 41
C452	80	H630	13	H690	29	K617C	20, 41
C474	80	H631	13	H691	29	K619	14
C491	81	H632	15	H692	29	K620	17
C810	81	H633	15	H693	33	K621A	43
F010	39, 78, 101	H634	31	H694A	29	K622A	43
F020	39, 78, 101	H635	26	H694B	29	K623	13, 14
F509	39, 78, 101	H636	26	H851	19	K624	15
F510	39, 78, 101	H638	31	H851A	19	K625A	30
F515	81	H639	27	H852	19	K625B	30
F516	81	H646	26	H852A	19	K625C	30
G203	61, 102	H647	26	H854	19	K625D	30
G210C	61, 102	H650	13	HEAD BAND	36	K626	20, 43
G210D	61, 102	H650A	14	HEAD-BAND/C	36	K627	43
H068A	59	H651	13	JL-101A	133	K629A	30
H503	51	H652	13	K101	37	K634	43
H503A	21, 51	H653	28	K101B	37	K635A	37
H503B	21, 51	H654	34	K297A	42	K635B	37
H520	52	H655	34	K297B	42	K636	30
H521	52	H656	34	K297C	42	K637	31, 42
H526	21, 52, 90, 100	H657	34	K297D	42	K638A	31
H527	21, 52, 100	H658	33	K297E	42	K638B	31
H528	19	H658A	33	K297F	42	K931	124
H537A	17, 55	H659	28	K297G	42	KOL-1	78, 101
H540	57	H660	14	K512	59	LS-101J	38
H540A	57	H661	14	K512A	59	LS-102J	38
H541	57	H662	33	K516A	49	M520	82

Index as per Supply Code

Supply Code	page	Supply Code	page	Supply Code	page	Supply Code	page
M521	82	NKD-BN107	75	NKD-E3S	67, 69, 71	NKD-N300	89
M523	82	NKD-BN108	75	NKD-E4	72	NKD-N306	89
M524A	82	NKD-BN109	76	NKD-E5/6Z-26	73	NKD-NCSD1	30
M528	82	NKD-BN111	76	NKD-E6	73	NKD-NCSD2	30
M531A	82	NKD-BN112	76	NKD-E7	73	NKD-NM2512	89
M540	82	NKD-BN113	76	NKD-E7B	73	NKD-NM3525	12
NE-110B	16, 55	NKD-BN114	76	NKD-E8	73	NKD-NM4525	12
NE-115B	16, 55	NKD-BN115	76	NKD-E8B	73	NKD-NM4525P	11
NE-120B	16, 55	NKD-BN500	76	NKD-E9	74	NKD-NM4537	12
NE-125B	16, 55	NKD-BN505	76, 91	NKD-EC2	78, 91, 101	NKD-NM4537P	11
NE-215B	16, 55	NKD-CPAPKITA	98	NKD-EPD1	16	NKD-NM4540	12
NE-220B	16, 55	NKD-E1-LMF	67	NKD-EPD2	16	NKD-NM4545	12
NF-E2/NK-9G	74	NKD-E10	74	NKD-EPDK	16	NKD-NM4550	12
NF-E2/NK-9K	74	NKD-E11	74	NKD-FE5GH	53, 100	NKD-NM4550P	11
NF-EEG/HELP	90	NKD-E12	71	NKD-FE5GH-07	53, 100	NKD-NM5565	12
NF-N/8000AA	133	NKD-E13L	72	NKD-FE5GH-12	53, 100	NKD-NUPREP	78, 91
NF-N/8000J	133	NKD-E13M	72	NKD-FE5GH-24	53, 100	NKD-O-RING	48
NF-N/8008J	133	NKD-E13S	72	NKD-FE6GH60	53, 100	NKD-P242B	38
NF-SCREEN/AC2	98	NKD-E16/700	74	NKD-FILZ	48	NKD-POUCH-JE215	79
NKD-110-150	40	NKD-E1L	67	NKD-I1/1F	70	NKD-POUCH-JE225	79
NKD-15026N	61, 102	NKD-E1LM	67	NKD-I1/2F	70	NKD-POUCH-JE922	79
NKD-155/15	80	NKD-E1M	67	NKD-I1/3F	70	NKD-PT1075	120
NKD-210-10/295	40	NKD-E1M-F	67	NKD-I1/4F	70	NKD-PT1080	120
NKD-210T	80	NKD-E1-LMF	67	NKD-I1/5F	70	NKD-PT1107	121
NKD-245/1	80	NKD-E1MS	68	NKD-I1/6F	70	NKD-PT1109	121
NKD-A205	75	NKD-E1MS-F	68	NKD-I2	70, 71	NKD-PT1263	121
NKD-aEEG/ Minibox	91	NKD-E1S	68	NKD-I23	70, 72	NKD-PT1267	121
NKD-AEEG/K1	87	NKD-E1SXS	68	NKD-I3	71	NKD-PT1295	121
NKD-AEEG/K2	87	NKD-E1XS	68	NKD-INBUS	82	NKD-PT1315	138
NKD-AEEG/K3	86	NKD-E2	74	NKD-INDUCTKIT	112	NKD-PT1327	120
NKD-BE-413AG	46	NKD-E23/5	69, 72	NKD-K10/C	60	NKD-PT1354	121
NKD-BM113A	11, 12	NKD-E23/7	68, 69, 72	NKD-KR26DIN15	91	NKD-PTRANS	120
NKD-BN100	75	NKD-E25	74	NKD-KR26DIN20	60	NKD-PTST3	139
NKD-BN101	75	NKD-E3L	67, 71	NKD-KR4B2SS	60	NKD-SLEEPKITA	94
NKD-BN106	75	NKD-E3M	67, 69, 71	NKD-KR4BDINS	60	NKD-SLEEPKITA921	96

Index as per Supply Code

Supply Code	page	Supply Code	page	Supply Code	page	Supply Code	page
NKD-SLEEPKITP	95	NKD-SLP14805BFT-5	118	NKD-SLP9102	112	P228B	131
NKD-SLEEPKITP921	97	NKD-SLP14807	119	NKD-SLP9102NO	112	P229A	130
NKD-SLP1250	103	NKD-SLP14808	119	NKD-SLPPS036	119	P229B	131
NKD-SLP1326	123	NKD-SLP14809	118	NKD-SLPP057	119	P259	130, 131
NKD-SLP1340	108, 110	NKD-SLP14810	118	NKD-STL	27	P260A	132
NKD-SLP1341	106, 108, 110	NKD-SLP14830	116	NKD-STL-BAND	27	P260B	132
NKD-SLP1342	106, 107, 110	NKD-SLP14830KIT	116	NKD-STLF	27	P263	127
NKD-SLP1343	107, 110	NKD-SLP14831	116	NKD-SYS1	64	P267	127
NKD-SLP1347	110	NKD-SLP14831KIT	116	NKD-SYS2	65	P310C	127
NKD-SLP1350	105, 109	NKD-SLP14834	117	NKD-SYS3	66	P311C	127
NKD-SLP1350	110	NKD-SLP14834KIT	117	NKD-TEMPHEARTS4009	38	P907	137
NKD-SLP1370	105	NKD-SLP14835	117	NKD-TEMPHEARTS4012	38	P908	137
NKD-SLP1370KIT	105	NKD-SLP14835KIT	117	NKD-TR4537	12	P909	136
NKD-SLP1387	106	NKD-SLP14836EU	117	NKD-TR4537P	11	R804	137
NKD-SLP1387KIT	106	NKD-SLP14837	117	P203A	128	R805	136
NKD-SLP1388	107	NKD-SLP1550	104	P203B	128	R806	136
NKD-SLP1388KIT	107	NKD-SLP1575	104	P203C	128	V921	138
NKD-SLP1389	108	NKD-SLP1575KIT	104	P203D	128	V922	138
NKD-SLP1389KIT	108	NKD-SLP15805AFT-10	118	P203E	129	V923	138
NKD-SLP1390	110	NKD-SLP1582	104	P203F	129	V928	138
NKD-SLP1391	109	NKD-SLP1726	123	P203G	129	V929	138
NKD-SLP1391KIT	109	NKD-SLP1730	123	P203H	129	Y085	43, 57
NKD-SLP1400	114	NKD-SLP1770	123	P204A	129	Y095A	124
NKD-SLP1401	115	NKD-SLP1770KIT	122	P204B	129	Y222	79
NKD-SLP1401H	115	NKD-SLP1771	123	P204C	129	Y763A	79
NKD-SLP1401NOH	115	NKD-SLP1771KIT	122	P204D	130	Y763B	79
NKD-SLP1401P	114	NKD-SLP1772KIT	122	P204E	131	Y848A	36
NKD-SLP1413	115	NKD-SLP9002S	111	P204F	131	Y849A	36
NKD-SLP1421H	115	NKD-SLP9003-L120	111	P204G	131	Y849B	36
NKD-SLP1468	113	NKD-SLP9003-L150	111	P204H	131	Y849C	36
NKD-SLP1468H	113	NKD-SLP9003-L40	111	P205A	132	Y849D	36
NKD-SLP1469	113	NKD-SLP9003-L60	111	P225F	126	Y855A	37
NKD-SLP14802AFT-5	118	NKD-SLP9003-L90	111	P225G	126	Y855B	37
NKD-SLP14804AFT-5	118	NKD-SLP9101	112	P225H	126	Y855C	37
NKD-SLP14804BFT-5	118	NKD-SLP9101NO	112	P228A	130	Y855D	38

Supply Code	page	Model Code	page	Model Code	page	Model Code	page
Y855E	38	BM-710S	17	NE-124A	52	NM-451B	28
YE-102J	37	BM-800S	17	NE-126A	52	NM-501B	33
YE-103J	36	BM-840S	14	NE-131B	19	NM-502B	34
YZ-0172	28, 43	BM-900B	31	NE-132B	19	NM-503B	34
		BM-901B	31	NE-134A	21, 52, 90, 100	NM-505B	34
Index as per Model Code		BM-910B	20, 43	NE-136A	21, 52, 100	NM-506B	33
Model Code	page	BM-911B	43	NE-136B	19	NM-511S	34
2114-042676B	63	BM-990B	30	NE-142B	19	NM-522S	34
2114-064171	20, 41	DR-531B-10B	36	NE-224S	17, 55	NM-531S	34
6114-009395	82	DR-531B-14	36	NE-301A	57	NM-540B	34
6143-009402	26	DR-531B-15	36	NE-301B	57	NM-550B	33
6143-009411	26	EX212BK20	42	NE-311A	57	NM-551B	33
6143-011854	132	EX212BW20	42	NE-420B	20, 41	NM-640S	14
6143-011863	132	EX212GR20	42	NE-430B	20, 41	NM-640T	14
6144-007171	43, 57	EX212PR20	42	NE-432B	20, 41	NM-710T	17
BC-102B	59	EX212RE20	42	NM-121T	13	NM-715S	17
BC-112B	59	EX212WH20	42	NM-130S	14	NM-730S	17
BE-0007	63	EX212YE20	42	NM-131T	13	NM-830S	17
BE-0008	63	F-150M	61, 102	NM-151T	13	NM-850S	17
BE-403B	49	F-150S	61, 102	NM-220T	15	NSE-25BL	81
BM-001B	30	FQW110-2-150	40	NM-250T	15	NSE-40BL	81
BM-002B	30	FQW110-3-150	40	NM-310Y	33	RC-86014-H5P	37
BM-003B	30	FQW210-10B	40, 80	NM-312S	31	RC-86014-R4C	37
BM-004B	30	FV245-10	80	NM-314YL	29	RIW114-3	81
BM-005B	30	FV245-30	80	NM-314YS	29	RP-86004	36
BM-115S	13, 14	FV345-30	80	NM-315S	31	TG-920P	137
BM-121B	43	FV400-30	81	NM-316Y	29	TG-921T3	137
BM-121S	13, 14	JL-550T2	124	NM-317Y3	29	TG-970P	136
BM-211S	15	JL-900P	124	NM-319Y	29	TL-051S	130
BM-215S	15	L-150	61, 102	NM-320T	13	TL-052S	131
BM-312B	31, 42	NE-103A	51	NM-330T	13	TL-061S	130
BM-503B	43	NE-113A	21, 51	NM-410S	26	TL-062S	131
BM-504B	43	NE-114A	21, 51	NM-420S	26	TL-201T	126
BM-701B	37	NE-121B	19	NM-422B	27	TL-201T1	126
BM-702B	37	NE-121J	19	NM-450S	28	TL-220T	126

Index as per Model Code

Model Code	page
TL-260T	132
TL-271T	128, 129
TL-271T3	129, 131
TL-272T	128, 129
TL-272T3	129, 131
TL-273T	128, 129
TL-273T3	129, 131
TL-274T	128, 130
TL-274T3	129, 131
TL-630T3	127
TL-631T3	127
YE-006A	82
YE-010A	82
YE-017A	82
YG-111T	137
YG-120T	138
YG-121T	138
YG-122T	138
YG-125T	138
YG-135T	138
YG-211T	136
YG-213T	136
YS-093P2	127
YSI-409B	38
YZ-0019	39, 78, 101
YZ-002A3	37
YZ-002A4	37
YZ-002A5	37
YZ-002A6	38
YZ-002A7	38
Z-101BC	39, 78, 101
Z-181JE	39, 78, 101
Z-401CE	39, 78, 101



For Germany, Austria, Switzerland and Benelux
NIHON KOHDEN EUROPE GmbH, Vertrieb Deutschland
Albert-Einstein-Ring 9, 14532 Kleinmachnow, Deutschland
Fax: +49 (0) 33203/573-19, E-Mail: bestellung@nke.de

For Spain and Portugal
NIHON KOHDEN IBERICA S.L.
C/Ulises, 75A, 28043 Madrid, Spain
Fax: +34 (0) 91 300 46 76, E-Mail: info@nkib.es

For United Kingdom
NIHON KOHDEN UNITED KINGDOM Ltd.
Tolworth Tower – 19th Floor, Ewell Road, Surbiton,
Surrey KT6 7EL, United Kingdom
Fax: +44 (0) 20 83 90 46 75, E-Mail: mail@nihonkohden.co.uk

For France
NIHON KOHDEN FRANCE SARL
8, rue François Delage, 94230 Cachan, France
Fax: +33 (0) 1 49 08 93 32, E-Mail: info@nkfrance.fr

For Italy
NIHON KOHDEN ITALIA S.r.l.
Via Fratelli Bronzetti, 28, 24124 Bergamo, Italy
Fax: +39 (0) 35-232546, E-Mail: info@nkitaly.com

For all other European Countries, Russia, Turkey and Israel
NIHON KOHDEN EUROPE GmbH
Raiffeisenstraße 10, 61191 Rosbach v. d. H., Germany
Fax: +49 (0) 6003/827-599, E-Mail: info@nke.de